Surgical Holdings
Parkside Centre, Temple Farm Industrial Estate, Southend-on-Sea, Essex S52 5SJ

T. +44 (0)1702 602050
E. info@surgicalholdings.co.uk
www.surgicalholdings.co.uk

en

F-114 (Issue 3, July 2023)

Surgical Holdings – Non-Sterile Implants for Osteosynthesis: Kirschner / Arthrodesis Wires & Steinmann Pins

Intended Use
Single use implants used for the temporary fixation, correction or

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Til kirurgn

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stabilization of bones in the hands and feet to support normal healing processes.

Indications

- Use in hands and feet for:
 - Fixation of bone fractures and fragments.
 - Corrective osteotomies.
 - Joint arthrodesis.

Contraindications

There are no specific contraindications for these devices, though please pay attention to the warnings, precautions and possible adverse effects, below.

Product Instructions


The following instructions should be followed in order for your Class IIb Medical Device to remain safe.

If in a sterile state, in order to remain in this state, it must be kept in a sterile environment and handled with sterile rubber gloves and other necessary clothing. It must be used in a sterile operating theatre. All necessary precautions should be taken to avoid contamination, and for the device to remain in a sterile state.

The device should be handled with care to avoid injury as it may have sharp edges.

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These products should be decontaminated using machines that comply with EN ISO 15883-1 and ISO 15883-2.

Sterilisation shall be via a steriliser which conforms to EN 285 on a cycle of 134 -137°C for a minimum holding time of 3 minutes.

Attention Operating Surgeon

• These products are only to be implanted by a qualified medical practitioner, fully trained in the surgical technique, using compatible CE Marked instruments. Operative techniques (F-146) may be requested by emailing info@surgicalholdings.co.uk

- Skeletal pins and wires are temporary implants. These devices are specified for use for 30 – 120 days, after which they should be removed. Devices must not be reused. The skeletal wires should never be misused in any way.

Possible Adverse Effects

Infection, subclinical nerve damage, arterial puncture, tendon risk, metal sensitivity, loosening, migration, bending, or fracture of the device and/or bone.

This product contains Nickel. It is recommended that hospitals perform nickel allergy tests if there is any doubt regarding patient sensitivity.

Warnings and Precautions

- Do not use implant if the shelf life has expired.
- Do not use if the packaging is damaged, or there is a flaw in the sterile barrier.
- Implants must not be re-used or re-sterilized.
- Surgical Holdings skeletal pins and wires are non-magnetic. However, they have not been evaluated for safety and compatibility in Magnetic Resonance (MR) environment and have not been tested for heating or migration in the MR environment, unless otherwise specified on the product labels.

The surgeon should:

- Discuss all physical and psychological limitations inherent in the use of these devices with the patient.
- Select the correct size and type of implant.
- Warn the patient against premature weight-bearing activity levels, and the necessity of periodic medical follow up and cleaning of the wound around the wires.
- Select the patient to report any unusual changes around the operated site to his physician, e.g. sign of infection.

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Disposal
Products should be disposed in clinical waste as per hospital protocol.

REF **Product Range**
Kirschner and Arthrodesis wires are prefixed with codes KWH5xxxx-S and SHxxxx-S where xxxx = catalogue number.

Kirschner wires typically have one trocar point/one plain end and Arthrodesis wires two trocar points.

Steinmann pins are prefixed with codes SHxxxx-S where xxxx = catalogue number. Steinmann Pins typically have one trocar point and either a plain, square or triangular end.

Any serious incident which occurs involving this medical device must be reported to the manufacturer, Surgical Holdings, and the Competent Authority in the Country it has been used.

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