 Surgical Holdings
Parkside Centre, Temple Farm
Industrial Estate, Southend-on-Sea,
Essex SS2 5SJ

T. +44 (0)1702 602050
E. info@surgicalholdings.co.uk
www.surgicalholdings.co.uk

EC REF

Advena Ltd. Tower Business Centre,
2nd Flr., Tower Street, Swatar,
BKR 4013 Malta

CE 1639

F-115 (Issue 1, October 2021)

en

Surgical Holdings – Sterile Implants for Osteosynthesis: Kirschner / Arthrodesis Wires & Steinmann Pins

Product Instructions

The following instructions should be followed in order for your Class IIb Medical Device to remain safe.

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Hygienic Safety

In order to maintain hygiene, the device should be cleaned and sterilised as described. It should be handled with care to avoid injury as it may have sharp edges.

Product supplied in a sterile state, in order to remain in this state, it must be kept in a sterile environment and handled with sterile rubber gloves and other necessary clothing. It must be used in a sterile operating theatre.

All necessary precautions should be taken to avoid contamination, and for the device to remain in a sterile state.

Cleaning and Sterilisation

These products are supplied sterile. They have been sterilised using gamma radiation.

Attention Operating Surgeon

- These products are only to be implanted by a qualified medical practitioner using compatible CE Marked instruments.
- Skeletal pins and wires should only be used for fusion, skeletal traction, and fixation (to fixate bones while they are healing) in the areas of the body for which they were intended, namely hands and feet.
- Skeletal pins and wires are non-permanent implants: they should

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be removed after 120 days, and SHOULD NEVER BE RE-USED. The skeletal wires should never be misused in any way.

Possible Adverse Effects

This product contains Nickel. It is recommended that hospitals perform nickel allergy tests if there is any doubt regarding patient sensitivity.

Infection, subclinical nerve damage, arterial puncture, tendon risk, metal sensitivity, loosening, migration, bending, or fracture of the device and/ or bone.

Warnings and Precautions

The surgeon should:

- Discuss all physical and psychological limitations inherent in the use of these devices with the patient.
- Select the correct size and type of implant.
- Warn the patient against premature weight-bearing activity levels, and the necessity of periodic medical follow up and cleaning of the wound around the wires.
- The patient should be instructed to report any unusual changes around the operated site to his physician, e.g. sign of infection.
- Surgical Holdings skeletal pins and wires are non-magnetic.

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However, they have not been evaluated for safety and compatibility in Magnetic Resonance (MR) environment and have not been tested for heating or migration in the MR environment, unless otherwise specified on the product labels.

Disposal

Products should be disposed in clinical waste as per hospital protocol.

REF Product Range

Kirschner and Arthrodesis wires are prefixed with codes KWH5xxxx-S and SHxxxx-S where xxxx = catalogue number. Kirschner wires typically have one trocar point/one plain end and Arthrodesis wires two trocar points.

Steinmann pins are prefixed with codes SHxxxx-S where xxxx = catalogue number. Steinmann Pins typically have one trocar point and either a plain, square or triangular end.

Any serious incident which occurs involving this medical device must be reported to the manufacturer, Surgical Holdings, and the Competent Authority in the Country it has been used.

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