en
**Surgical Holdings – Sterile
Implants for Osteosynthesis:
Kirschner / Arthrodesis Wires &
Steinmann Pins**

Product Instructions
The following instructions should be
followed in order for your Class IIb
Medical Device to remain safe.

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Hygienic Safety
In order to maintain hygiene, the device
should be cleaned and sterilised as
described. It should be handled with
care to avoid injury as it may have sharp
edges.

Product supplied in a sterile state, in
order to remain in this state, it must
be kept in a sterile environment and
handled with sterile rubber gloves and
other necessary clothing. It must be
used in a sterile operating theatre.

All necessary precautions should be
taken to avoid contamination, and for
the device to remain in a sterile state.

Cleaning and Sterilisation
These products are supplied sterile.
They have been sterilised using gamma
radiation.

- Attention Operating Surgeon**
- These products are only to be
implanted by a qualified medical
practitioner using compatible CE
Marked instruments.
 - Skeletal pins and wires should only be
used for fusion, skeletal traction, and
fixation (to fixate bones while they
are healing) in the areas of the body
for which they were intended, namely
hands and feet.
 - Skeletal pins and wires are non-
permanent implants: they should be
removed after 120 days, and SHOULD
NEVER BE RE-USED. The skeletal
wires should never be misused in any
way.

Possible Adverse Effects
This product contains Nickel. It is
recommended that hospitals perform
nickel allergy tests if there is any doubt

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regarding patient sensitivity.
Infection, subclinical nerve damage,
arterial puncture, tendon risk, metal
sensitivity, loosening, migration,
bending, or fracture of the device and/
or bone.

- Warnings and Precautions**
The surgeon should:
- Discuss all physical and psychological
limitations inherent in the use of these
devices with the patient.
 - Select the correct size and type of
implant.
 - Warn the patient against premature
weight-bearing activity levels, and the
necessity of periodic medical follow
up and cleaning of the wound around
the wires.
 - The patient should be instructed to
report any unusual changes around the
operated site to his physician, e.g.
sign of infection.
 - Surgical Holdings skeletal pins and
wires are non-magnetic. However,
they have not been evaluated for
safety and compatibility in Magnetic
Resonance (MR) environment and
have not been tested for heating or
migration in the MR environment,
unless otherwise specified on the
product labels.

Disposal
Products should be disposed in clinical
waste as per hospital protocol.

REF Product Range
Kirschner and Arthrodesis
wires are prefixed with codes
KWHsxxxx-S and SHxxxx-S
where xxxx = catalogue number.

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Kirschner wires typically have
one trocar point/one plain end
and Arthrodesis wires two trocar
points.

Steinmann pins are prefixed with
codes SHxxxx-S where xxxx =
catalogue number. Steinmann
Pins typically have one trocar
point and either a plain, square or
triangular end.

Any serious incident which occurs
involving this medical device must be
reported to the manufacturer, Surgical
Holdings, and the Competent Authority
in the Country it has been used.

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