


urology

retractors	PAGE 301
forceps	PAGE 301
urethral bougies	PAGE 302
urethral sounds	PAGE 302
miscellaneous	PAGE 304

RETRACTORS / FORCEPS


Walton bladder retractor

SC-UR/001 parallel action,
with convex lateral blades - 25mm x 70mm
and curved central blade - 70mm x 215mm



Thompson Walker bladder retractor

SC-UR/002 one posterior blade
95mm x 30mm, one anterior wall blade
175mm x 25mm tapering to 20mm,
two pairs of lateral swivel blades,
shallow: 55mm x 22mm
deep: 75mm x 25mm



Millin retropubic retractor



self retaining

SC-UR/003

two short blades 25 x 55mm,
two long lateral blades 25 x 80mm,
one centre blade 115 x 445mm

Millin retractor blades

- SC-UR/004a** short 25 x 55mm
- SC-UR/004b** long 25 x 80mm
- SC-UR/004c** centre 115 x 445mm



Millin vulsellum forceps

screw joint for capsule, 225mm long

SC-UR/018 8-9 teeth, jaws 10mm wide




Millin dissecting spring forceps

230mm long

SC-UR/014 3-4 teeth


URETHAL BOUGIES AND SOUNDS



Winsbury White urethral bougie

French gauge, child size, 175mm long

- SC-UR/063** 1-4 to 3-6
- SC-UR/064** 5-8 to 7-10
- SC-UR/065** 9-12 to 11-14



Canny Ryall metal urethral bougie

English gauge, double-ended, tapering
160mm long

- | | |
|----------------------------------|----------------------------------|
| SC-UR/075 1-4 and 3-6 | SC-UR/079 17-20 and 19-22 |
| SC-UR/076 5-8 and 7-10 | SC-UR/080 21-23 and 22-24 |
| SC-UR/077 9-12 and 11-14 | SC-UR/081 26-28 and 29-31 |
| SC-UR/078 13-16 and 15-18 | |



Clutton urethral sounds

French gauge
tapering, 265mm long

- SC-UR/007** set of 12 contains
SC-UR/008 to SC-UR/008k

- | | |
|---------------------------|---------------------------|
| SC-UR/008 6 - 10 | SC-UR/008f 18 - 22 |
| SC-UR/008a 8 - 12 | SC-UR/008g 20 - 24 |
| SC-UR/008b 10 - 14 | SC-UR/008h 22 - 26 |
| SC-UR/008c 12 - 16 | SC-UR/008i 24 - 28 |
| SC-UR/008d 14 - 18 | SC-UR/008j 26 - 30 |
| SC-UR/008e 16 - 20 | SC-UR/008k 28 - 32 |



Lister urethral sounds

English gauge
with olive ends, straight,
280mm long

- | | |
|---------------------------|--------------------------|
| SC-UR/010 0 - 1 | SC-UR/010f 4 - 7 |
| SC-UR/010a 0 - 2 | SC-UR/010g 5 - 8 |
| SC-UR/010b 0.5 - 3 | SC-UR/010h 6 - 9 |
| SC-UR/010c 1 - 4 | SC-UR/010i 7 - 10 |
| SC-UR/010d 2 - 5 | SC-UR/010j 8 - 11 |
| SC-UR/010e 3 - 6 | SC-UR/010k 9 - 12 |

URETHAL BOUGIES

Buckston Browne urethral bougie

English gauge, 295mm long

SC-UR/020	4-6	SC-UR/025	9-11	SC-UR/030	14-16
SC-UR/021	5-7	SC-UR/026	10-12	SC-UR/031	15-17
SC-UR/022	6-8	SC-UR/027	11-13	SC-UR/032	16-18
SC-UR/023	7-9	SC-UR/028	12-14		
SC-UR/024	8-10	SC-UR/029	13-15		



Clutton urethral bougie

French gauge, 260mm long

SC-UR/033	6 - 10	SC-UR/037	14 - 18	SC-UR/041	22 - 26
SC-UR/034	8 - 12	SC-UR/038	16 - 20	SC-UR/042	24 - 28
SC-UR/035	10 - 14	SC-UR/039	18 - 22	SC-UR/043	26 - 30
SC-UR/036	12 - 16	SC-UR/040	20 - 24	SC-UR/044	28 - 32



Guyon urethral bougie

French gauge, 275mm long

SC-UR/046	6-24
------------------	------



Lister urethral bougie

English gauge, 270mm long

SC-UR/047	4-7	SC-UR/051	8-11	SC-UR/055	12-15
SC-UR/048	5-8	SC-UR/052	9-12	SC-UR/056	13-16
SC-UR/049	6-9	SC-UR/053	10-13	SC-UR/057	14-17
SC-UR/050	7-10	SC-UR/054	11-14	SC-UR/058	15-18




Wyndham Powell urethral bougie

French gauge, 170mm long

SC-UR/059	12	SC-UR/062a	19	SC-UR/065	26
SC-UR/060	13	SC-UR/059b	20	SC-UR/066a	27
SC-UR/061	14	SC-UR/060b	21	SC-UR/067a	28
SC-UR/062	15	SC-UR/061b	22	SC-UR/068a	29
SC-UR/059a	16	SC-UR/062b	23	SC-UR/069a	30
SC-UR/060a	17	SC-UR/063	24	SC-UR/070a	31
SC-UR/061a	18	SC-UR/064	25	SC-UR/071a	32




MISCELLANEOUS



Millin bladder neck spreader

self-retaining, angled jaw

SC-UR/015 240mm long



Millin boomerang needle holder

with two needles, flat shaft, slot eye


SC-UR/016 230mm long

Millin needle

SC-UR/017a small

SC-UR/017b medium

SC-UR/017c large



Sieve

SC-GI/5071 203mm x 50mm
