

general instruments

artery forceps	PAGE 3
clamps	PAGE 8
spring forceps	PAGE 9
tissue forceps	PAGE 12
sponge holding and sterilising forceps	PAGE 14
other forceps	PAGE 15
scissors	PAGE 16
needle holders	PAGE 24
combined scissor and needle holders	PAGE 28
bag and towel clips	PAGE 29
retractors	PAGE 30
probes and dissectors	PAGE 39
miscellaneous	PAGE 40
hollow ware	PAGE 41
silicone matting	PAGE 42
wire products / sterilising baskets	PAGE 43
optical fibre cables	PAGE 45

ARTERY FORCEPS

Adson artery forceps

box joint, 180mm long

serrated jaws


SC-GI/034 straight

SC-GI/034a curved

1-2 teeth

SC-GI/034b straight

SC-GI/034c curved



Cairns artery forceps


box joint, half serrated jaws

145mm long

SC-GI/415 straight

SC-GI/020 curved

SC-GI/021 curved to side



Crile artery forceps

box joint, serrated jaws

straight


SC-GI/014 140mm long

SC-GI/014a 160mm long

curved

SC-GI/015 140mm long

SC-GI/015a 160mm long



Cushing artery forceps


box joint, fine, half serrated jaws, 140mm long

straight

SC-GI/001

curved

SC-GI/002



Dunhill artery forceps

box joint, fine serrations

straight


SC-GI/065 135mm long

SC-GI/065a 190mm long


curved

SC-GI/043 135mm long

SC-GI/468 180mm long



ARTERY FORCEPS





serrated jaw
straight
SC-GI/007 125mm long
curved
SC-GI/008 125mm long
fully curved
SC-GI/009 125mm long

Halstead mosquito forceps
standard pattern, haemostatic, box joint


1-2 teeth
straight
SC-GI/010 125mm long
SC-GI/012 190mm long
curved
SC-GI/011 125mm long
SC-GI/013 190mm long




micro, 100mm long
SC-GI/003 straight, serrated
SC-GI/004 curved, serrated
SC-GI/005 straight, 1-2 teeth
SC-GI/006 curved, 1-2 teeth

Heiss artery forceps
box joint, 205mm long
serrated jaws
SC-GI/016 straight
SC-GI/017 curved




1-2 teeth
SC-GI/016a straight
SC-GI/017b curved


Kelly artery forceps box joint, half serrated jaws
straight
SC-GI/022 140mm long
SC-GI/023 160mm long



curved
SC-GI/024 140mm long
SC-GI/025 160mm long



Fraser Kelly artery forceps curved
SC-GI/468a 180mm





Kocher artery forceps
box joint, 1-2 teeth

straight
SC-GI/052 130mm long
SC-GI/053 150mm long
SC-GI/054 180mm long
SC-GI/055 200mm long
SC-GI/056 230mm long



curved
SC-GI/057 150mm long
SC-GI/058 180mm long
SC-GI/059 200mm long
SC-GI/060 230mm long

Lawson Tait artery forceps
box joint, straight jaw
140mm long
SC-GI/SH5924


ARTERY FORCEPS

Moynihan artery forceps

box joint,
half-serrated jaws

straight
SC-GI/018 145mm long
SC-GI/068 180mm long



curved
SC-GI/019 145mm long
SC-GI/069 180mm long



Roberts artery forceps

screw joint, curved on flat, serrated

SC-GI/070 230mm long




Rochester-Pean artery forceps

box joint

straight
SC-GI/044 165mm long
SC-GI/045 190mm long
SC-GI/046 205mm long
SC-GI/047 240mm long

curved
SC-GI/048 165mm long
SC-GI/049 180mm long
SC-GI/050 205mm long
SC-GI/051 240mm long




Spencer Wells artery forceps

box joint,
serrated jaws

straight
SC-GI/035 127mm long
SC-GI/036 152mm long
SC-GI/037 177mm long
SC-GI/038 203mm long

curved
SC-GI/039 127mm long
SC-GI/040 152mm long
SC-GI/041 177mm long
SC-GI/042 203mm long


OTHER SIZES AVAILABLE



Dandy artery forceps

box joint, curved in plane


SC-GI/076 140mm long




ARTERY FORCEPS



non-traumatic jaw, straight




non-traumatic jaw, curved



Carmalt angiotribe forceps
box joint, non-traumatic,
light pattern

straight
SC-GI/062 165mm long
SC-GI/063 200mm long
SC-GI/064 230mm long


curved
SC-GI/061 165mm long



Cooley artery forceps
box joint, non-traumatic,
curved

SC-GI/071 160mm long
SC-GI/072 200mm long


SC-GI/073 150mm long
large, curved 60°
SC-GI/074 140mm long
large, strongly curved 90°
SC-GI/075 160mm long
curved 30°



De Bakey Bainbridge forceps

box joint, non-traumatic, haemostatic,
180mm long


SC-GI/078 straight
SC-GI/077 curved



De Bakey Bainbridge forceps

box joint, non-traumatic,
curved


SC-GI/079 165mm long



De Bakey forceps

box joint, non-traumatic,
curved

SC-GI/080 165mm long
SC-GI/081 190mm long
SC-GI/082 203mm long
SC-GI/083 215mm long
SC-GI/084 220mm long



De Bakey forceps

box joint, non-traumatic,
curved


SC-GI/086 230mm long

ARTERY FORCEPS

De Bakey forceps

box joint, non-traumatic,
curved


- SC-GI/088** 230mm long
- SC-GI/089** 240mm long
- SC-GI/090** 245mm long



De Bakey forceps

box joint, non-traumatic,
curved


- SC-GI/091** 250mm long
- SC-GI/092** 270mm long
- SC-GI/093** 280mm long



De Bakey forceps

box joint, non-traumatic,
curved


- SC-GI/094** 150mm long



De Bakey forceps

box joint, non-traumatic,
curved


- SC-GI/096** 180mm long
- SC-GI/097** 203mm long




De Bakey forceps

box joint, non-traumatic,
straight

- SC-GI/100** 160mm long
- SC-GI/104** 165mm long
- SC-GI/101** 203mm long
- SC-GI/098** 215mm long
- SC-GI/102** 220mm long
- SC-GI/103** 230mm long
- SC-GI/099** 245mm long



ARTERY FORCEPS AND CLAMPS



De Bakey Bainbridge forceps

box joint, non-traumatic,
haemostatic
155mm long

SC-GI/105 curved
SC-GI/106 straight



Dieffenbach bulldog clamp
straight, serrated

SC-GI/112
1.2mm fine jaw, 40mm long
SC-GI/113
1.2mm fine jaw, 45mm long
SC-GI/108
1.8mm jaw, 30mm long
SC-GI/109
1.8mm jaw, 35mm long
SC-GI/110
2.5mm jaw, 40mm long
SC-GI/111
2.5mm jaw, 45mm long
SC-GI/107
6mm heavy jaw, 60mm long
SC-GI/455
20mm jaw, 75mm long
SC-GI/455a
30mm jaw, 85mm long
SC-GI/455b
45mm jaw, 105mm long
SC-GI/455c
65mm jaw, 120mm long

Dieffenbach bulldog clamp
straight, plain

SC-GI/129
1.2mm fine jaw, 40mm long
SC-GI/131
1.2mm fine jaw, 45mm long
SC-GI/121
1.8mm jaw, 30mm long
SC-GI/123
1.8mm jaw, 35mm long
SC-GI/125
2.5mm jaw, 40mm long
SC-GI/127
2.5mm jaw, 45mm long



Dieffenbach bulldog clamp
curved, serrated

SC-GI/119
1.2mm fine jaw, 40mm long
SC-GI/120
1.2mm fine jaw, 45mm long
SC-GI/115
1.8mm jaw, 30mm long
SC-GI/116
1.8mm jaw, 35mm long
SC-GI/117
2.5mm jaw, 40mm long
SC-GI/118
2.5mm jaw, 45mm long
SC-GI/114
6mm heavy jaw, 60mm long
SC-GI/456
20mm jaw, 75mm long
SC-GI/456a
30mm jaw, 85mm long
SC-GI/456b
45mm jaw, 105mm long
SC-GI/456c
65mm jaw, 120mm long

Dieffenbach bulldog clamp
curved, plain

SC-GI/130
1.2mm fine jaw, 40mm long
SC-GI/132
1.2mm fine jaw, 45mm long
SC-GI/122
1.8mm jaw, 30mm long
SC-GI/124
1.8mm jaw, 35mm long
SC-GI/126
2.5mm jaw, 40mm long
SC-GI/128
2.5mm jaw, 45mm long

Blalock clamp

serrated, with hole in handle

straight
SC-GI/480
11mm jaw, 38mm long
SC-GI/482
14mm jaw, 48mm long
SC-GI/493
17mm jaw, 60mm long
SC-GI/481
27mm jaw, 70mm long



curved
SC-GI/478
11mm jaw, 38mm long
SC-GI/479
14mm jaw, 48mm long
SC-GI/491
17mm jaw, 60mm long
SC-GI/492
27mm jaw, 70mm long

Blalock clamp

plain, with hole in handle

straight
SC-GI/486
11mm jaw, 38mm long
SC-GI/487
14mm jaw, 48mm long
SC-GI/488
17mm jaw, 60mm long
SC-GI/489
27mm jaw, 70mm long

curved
SC-GI/483
11mm jaw, 38mm long
SC-GI/484
14mm jaw, 48mm long
SC-GI/485
17mm jaw, 60mm long
SC-GI/490
27mm jaw, 70mm long

SPRING FORCEPS

Adsons dissecting forceps


127mm long

SC-GI/257 1-2 teeth

SC-GI/258 lightweight, micro

SC-GI/259 serrated

SC-GI/260 serrated, lightweight, micro




Block end dissecting forceps

serrated

SC-GI/255 127mm long

SC-GI/256 152mm long

SC-GI/256a 203mm long




Bonney dissecting forceps

178mm long

SC-GI/267 serrated, 1-2 teeth

SC-GI/268 1-2 teeth



Debakey atraumatic dissecting forceps

straight, non traumatic jaws

1.5mm

SC-GI/238a 150mm long

SC-GI/239a 196mm long

SC-GI/240a 240mm long

SC-GI/241a 180mm long

2.5mm

SC-GI/238b 150mm long

SC-GI/239b 196mm long

SC-GI/240b 240mm long

SC-GI/241b 180mm long

2mm

SC-GI/238 150mm long

SC-GI/239 196mm long

SC-GI/240 240mm long

SC-GI/241 180mm long


2.8mm

SC-GI/238c 150mm long

SC-GI/239c 196mm long

SC-GI/240c 240mm long


SC-GI/241c 180mm long



Charnley-McIndoe dissecting forceps

150mm long

SC-OR/N045a 2-3 teeth



SPRING FORCEPS



Dissecting forceps

turnover end, serrated jaws

SC-GI/242 127mm long

SC-GI/243 152mm long

SC-GI/244 178mm long

SC-GI/245 203mm long



Dunhill dissecting forceps

light, fine serrated jaws

SC-GI/457 145mm long



Fredrick Coole dissecting forceps

extra fine serrated jaws

125mm long

SC-GI/066 straight

SC-GI/067 curved



Gillies dissecting forceps

1-2 teeth

SC-PL/052 152mm long



Jackson Burrows dissecting forceps

140mm long

SC-OR/N037 tongue and groove, 1-2 teeth

SC-OR/N037a tongue and groove only



SPRING FORCEPS


Jeans dissecting forceps

serrated jaws , 140mm long

SC-GI/269 fine

SC-GI/270 1-2 teeth

SC-GI/271 2-3 teeth 178mm long



Lanes dissecting forceps

1-2 teeth on serrated jaws


SC-GI/246 127mm long

SC-GI/247 152mm long

SC-GI/248 178mm long

SC-GI/249 2-3 teeth 178mm long

SC-GI/250 2-3 teeth 203mm long




McIndoe dissecting forceps

152mm long

SC-GI/476 serrated


SC-GI/477 1-2 teeth



McIndoe modified dissecting forceps

152mm long, 1mm tip, modified for fine powerful dissection with strong action tip

SH/5961 serrated



Treves dissecting forceps


1-2 teeth

SC-GI/251 127mm long

SC-GI/252 152mm long

SC-GI/253 2-3 teeth 127mm long

SC-GI/254 2-3 teeth 152mm long



Waugh dissecting forceps

fine, serrated

SC-GI/261 152mm long

SC-GI/262 178mm long

SC-GI/263 203mm long


Waugh dissecting forceps

fine, 1-2 teeth


SC-GI/264 152mm long

SC-GI/265 178mm long

SC-GI/266 203mm long



TISSUE FORCEPS




Allis tissue forceps
box joint, 150mm long

SC-GI/272 3-4 teeth

SC-GI/273 4-5 teeth


SC-GI/414 5-6 teeth



Littlewood tissue forceps

box joint, 2-3 teeth

SC-GI/275 190mm long




Babcock tissue forceps

box joint, serrated jaws

SC-GI/276 160mm long

SC-GI/277 200mm long

SC-GI/278 230mm long




Duval tissue forceps

box joint

SC-GI/279

13mm wide, 203mm long



TISSUE FORCEPS


Lane tissue forceps

box joint, 1-2 teeth,
fenestrated

SC-GI/280 127mm long

SC-GI/281 152mm long


SC-GI/282 197mm long



Rutherford Morrison tissue forceps

box joint


SC-GI/283 4-5 teeth, 152mm long



Shaw tissue forceps

box joint

SC-GI/458 1-2 teeth, 178mm long




Guys Gallabine suture-holding forceps

screw joint, angled flat blades

SC-GI/450 125mm long

SC-GI/4402 150mm long



Russian dissecting forceps


with fenestrated end

SC-GI/SH-5503b 150mm long


SC-GI/SH-5503 200mm long

SC-GI/403 230mm long

SC-GI/SH-5503a 250mm long




SPONGE HOLDING AND STERILISING FORCEPS



Rampley sponge holding forceps
box joint, straight, serrated jaws


SC-GI/353 178mm long

SC-GI/354 241mm long



Foersters sponge holding forceps
box joint, straight, serrated jaws

SC-GI/354a 178mm long




Lister sinus forceps
box joint, serrated tip

SC-GI/367 127mm long


SC-GI/368 152mm long

SC-GI/369 178mm long



Cheattle sterilizer forceps
screw joint

SC-GI/364 267mm long



Harrison bowl sterilizer forceps
screw joint, angled,
double jaw

SC-GI/365 356mm long

SC-GI/366 457mm long

OTHER FORCEPS

Hunter splinter forceps

serrated grooved pointed jaws


SC-GI/360 115mm long



Michel skin holding forceps

with magazine for Michel clips


SC-GI/358 1-2 teeth, 127mm long



Ward bone cutting forceps

slightly angled on flat



SC-GI/459 180mm long



Childe approximating forceps

with magazine for Michel clips


SC-GI/469 180mm long



Michel forceps

box joint, double-ended, for inserting and removing clips stainless steel

SC-GI/357 130mm long



Michel clip extracting forceps

box joint


SC-GI/359 125mm long

Robins anchoring forceps


140mm long

SC-GI/453 5mm, with holder for diathermy leads

SC-GI/453a 13mm, with holder for suction tubing



SCISSORS premium quality stainless steel



Dressing scissors

blunt blunt

straight

SC-GI/133 130mm long

SC-GI/136 146mm long

SC-GI/139 165mm long



curved

SC-GI/142 130mm long

SC-GI/143 146mm long

SC-GI/144 165mm long

SC-GI/145 185mm long

SC-GI/404 147mm long
with one serrated blade

blunt sharp

straight

SC-GI/134 130mm long

SC-GI/137 146mm long

SC-GI/140 165mm long

sharp sharp

straight


SC-GI/135 130mm long

SC-GI/138 146mm long

SC-GI/141 165mm long

curved

SC-GI/146 130mm long



Dressing scissors, left-handed

with one black bow, 145mm long

blunt/blunt

SC-GI/192 straight

SC-GI/193 curved

blunt/sharp


SC-GI/194 straight

SC-GI/195 curved

sharp/sharp

SC-GI/196 straight


SC-GI/197 curved



Littler dissecting scissors

fine pattern,
with eye for suture

SC-PL/001f 115mm long



Kilner dissecting scissors

fine pattern

SC-PL/001a 115mm long


Metzenbaum scissors

straight

- SC-GI/405** 180mm long
- SC-GI/406** 147mm long
- SC-GI/406a** 200mm long
- SC-GI/406b** 230mm long
- SC-GI/406c** 250mm long
- SC-GI/406d** 280mm long

curved on flat

- SC-GI/407** 185mm long
- SC-GI/408** 147mm long
- SC-GI/408a** 200mm long
- SC-GI/408b** 230mm long
- SC-GI/408c** 250mm long
- SC-GI/408d** 280mm long




Mayo scissors

straight

- SC-GI/409** 140mm long
- SC-GI/409a** 160mm long
- SC-GI/411** 190mm long
- SC-GI/411a** 230mm long

curved on flat


- SC-GI/410** 140mm long
- SC-GI/410a** 160mm long
- SC-GI/412** 190mm long
- SC-GI/412a** 230mm long



Fergusson's scissors


angled on flat

- SC-GI/148** 180mm long



Paynes suture scissors


- SC-GI/152** 100mm long




Ligature stitch scissors

fine pattern

- SC-GI/153** 115mm long



SCISSORS



Lister bandage scissors

SC-GI/154 147mm long

SC-GI/155 178mm long




Lister bandage scissors

angled, with probe end

SC-GI/361 147mm long

SC-GI/362 178mm long




McIndoe scissors

178mm long

SC-GI/149 straight

SC-GI/150 curved on flat

SC-GI/150sc supercut, curved




Thomas Strabismus scissors

114mm long


SC-OP/030 straight

SC-OP/029 curved




Heath stitch scissors

SC-GI/151 152mm long



"Careless" scissors

SC-GI/474 165mm long



Tubing scissors

SC-GI/147 130mm long



Utility scissors

SC-GI/466 190mm long

with two gold bows for identification **SCISSORS TUNGSTEN CARBIDE HARD INSERT**

Dressing scissors


tungsten carbide hard insert,
blunt/blunt

straight

- SC-GI/430** 130mm long
- SC-GI/432** 145mm long
- SC-GI/434** 165mm long
- SC-GI/436** 185mm long

curved

- SC-GI/431** 130mm long
- SC-GI/433** 145mm long
- SC-GI/435** 165mm long



Dressing scissors


tungsten carbide hard insert,
blunt/sharp

straight

- SC-GI/423** 130mm long
- SC-GI/425** 145mm long
- SC-GI/427** 165mm long
- SC-GI/429** 185mm long

curved

- SC-GI/424** 130mm long
- SC-GI/426** 145mm long
- SC-GI/428** 165mm long



Dressing scissors


tungsten carbide hard insert
sharp/sharp

straight

- SC-GI/437** 130mm long
- SC-GI/439** 145mm long
- SC-GI/441** 165mm long
- SC-GI/443** 185mm long

curved

- SC-GI/438** 130mm long
- SC-GI/440** 145mm long
- SC-GI/442** 165mm long



Mayo dressing scissors


tungsten carbide hard insert

straight


- SC-GI/419** 140mm long
- SC-GI/419a** 160mm long
- SC-GI/421** 190mm long
- SC-GI/421a** 230mm long

curved

- SC-GI/420** 140mm long
- SC-GI/420a** 160mm long
- SC-GI/422** 190mm long
- SC-GI/422a** 230mm long



SCISSORS TUNGSTEN CARBIDE HARD INSERT with two gold bows for identification



Mayo-Harrington dressing scissors

tungsten carbide hard insert

230mm long


SC-GI/417 straight

SC-GI/418 curved

200mm long

SC-GI/162a straight

SC-GI/413b curved



Metzenbaum scissors

tungsten carbide hard insert

straight

SC-GI/445 145mm long

SC-GI/448 185mm long

SC-GI/471 200mm long

SC-GI/471a 230mm long

SC-GI/471b 250mm long

SC-GI/471c 280mm long

curved

SC-GI/446 145mm long


SC-GI/449 185mm long

SC-GI/449a 200mm long

SC-GI/449b 230mm long

SC-GI/449c 250mm long

SC-GI/449d 280mm long




Baby Metzenbaum scissors

tungsten carbide,
hard insert, 145mm long

SC-GI/445a straight

SC-GI/446a curved




Ligature scissors

tungsten carbide hard
insert, curved


SC-GI/444 145mm long

SC-GI/447 185mm long


with one black and one gold bow for identification **SUPERCUT TUNGSTEN CARBIDE SCISSORS**
Mayo Harrington scissors
 supercut, tungsten carbide

 230mm long
 straight **SC-GI/162**
 curved **SC-GI/413**

 200mm long
 straight **SC-GI/162b**
 curved **SC-GI/413b**
Mayo scissors
 supercut, tungsten carbide

 straight **SC-GI/163** 140mm long
SC-GI/163a 160mm long
SC-GI/165 190mm long
SC-GI/165a 230mm long

 curved **SC-GI/164** 140mm long
SC-GI/164a 160mm long
SC-GI/166 190mm long
SC-GI/166a 230mm long

Metzenbaum dissecting scissors
 supercut, tungsten carbide

 straight **SC-GI/190tc** 147mm long
SC-GI/188tc 185mm long
SC-GI/188tca 200mm long
SC-GI/188tcb 230mm long
SC-GI/188tcc 250mm long
SC-GI/188tcd 280mm long

 curved **SC-GI/191tc** 147mm
SC-GI/189tc 185mm
SC-GI/186tca 200mm
SC-GI/188tcb 230mm
SC-GI/188tcc 250mm
SC-GI/188tcd 280mm

Dressing scissors
 supercut, tungsten carbide


 blunt/blunt
 straight **SC-GI/174** 130mm long
SC-GI/176 145mm long
SC-GI/178 165mm long
SC-GI/180 185mm long

 blunt/sharp
 straight **SC-GI/167** 130mm long
SC-GI/169 145mm long
SC-GI/171 165mm long
SC-GI/173 185mm long


 sharp/sharp
 straight **SC-GI/181** 130mm long
SC-GI/183 145mm long
SC-GI/185 165mm long
SC-GI/187 185mm long

 curved **SC-GI/175** 130mm long
SC-GI/177 145mm long
SC-GI/179 165mm long
SC-GI/180a 185mm long

 curved **SC-GI/168** 130mm long
SC-GI/170 145mm long
SC-GI/172 165mm long
SC-GI/173a 185mm long

 curved **SC-GI/182** 130mm long
SC-GI/184 145mm long
SC-GI/186 165mm long
SC-GI/187a 185mm long


SUPER CUT OPERATING SCISSORS with two black bows for identification




Dressing scissors

super cut, blunt/blunt
145mm long

SC-GI/160 straight

SC-GI/161 curved



Dressing scissors

super cut, blunt/sharp
145mm long

SC-GI/156 straight

SC-GI/157 curved




Dressing scissors

super cut, sharp/sharp
145mm long

SC-GI/158 straight

SC-GI/159 curved



Mayo scissors

super cut

straight

SC-GI/166a 140mm long

SC-GI/166e 160mm long


SC-GI/166b 170mm long

curved

SC-GI/166c 140mm long

SC-GI/166f 160mm long

SC-GI/166d 170mm long



Boyd scissors

super cut, blunt/blunt

SC-THO/047 180mm long

with two black bows for identification **SUPER CUT OPERATING SCISSORS**

Metzenbaum dissecting scissors


super cut

straight (length)

SC-GI/190 147mm **SC-GI/188b** 230mm
SC-GI/188 185mm **SC-GI/188c** 250mm
SC-GI/188a 200mm **SC-GI/188d** 280mm

curved (length)

SC-GI/191 147mm **SC-GI/186b** 230mm
SC-GI/189 185mm **SC-GI/186c** 250mm
SC-GI/186a 200mm **SC-GI/186d** 280mm



Stevens tenotomy scissors

super cut

curved/blunt


SC-OP/424 114mm long
SC-OP/424a 125mm long
SC-OP/424gb 135mm long

curved/sharp

SC-OP/425 114mm long
SC-OP/425a 140mm long
SC-OP/425b 160mm long

SC-OP/426 straight/blunt, 114mm long
SC-OP/427 straight/sharp, 114mm long
SC-OP/427c straight/sharp, 180mm long
with cutout, 110mm long


SC-OP/427a straight/blunt
SC-OP/427b curved/blunt



Strabismus scissors

super cut
114mm long


SC-GI/475 straight
SC-GI/475a curved




Iris scissors

super cut
114mm long

SC-OP/476 straight
SC-OP/476a curved




NEEDLE HOLDERS



Bruce Clarke's needle holder

screw joint


SC-GI/205 127mm long



Crile-Wood needle holder

straight , baby, fine

SC-GI/199 150mm long




Kilner needle holder

box joint, cross-serrated jaw

SC-GI/200 165mm long

SC-GI/201 133mm long




Lawrence needle holder

fine jaw, 150mm long

SC-PL/012a

SC-PL/012b tungsten carbide jaw



Needle holder

Surgical Holdings pattern, box joint,
lightweight, with curved plain jaws
and one offset bow for thumb

SC-GI/209 152mm long

NEEDLE HOLDERS


MacPhail's needle holder

screw joint, with one copper jaw

SC-GI/206 159mm long

SC-GI/207 178mm long

SC-GI/208 209mm long




Mathieu needle holder

screw joint

SC-GI/210 172mm long

SC-GI/211 209mm long

SC-GI/212 243mm long



Mayo needle holder

box joint, straight



wide cross serrated jaws


SC-GI/224 159mm long

SC-GI/398 203mm long



narrow cross serrated jaws

SC-GI/225 159mm long



Naughton Morgan needle holder

box joints, cross serrated jaws


SC-GI/399 185mm long

SC-GI/400 254mm long

Tungsten carbide

SC-GI/401 185mm long


SC-GI/402 254mm long




Rochester-Mayo needle holder

screw joint

SC-GI/213 178mm long




TUNGSTEN CARBIDE NEEDLE HOLDERS with two gold bows for identification



Berry needle holder
tungsten carbide jaws

SC-GI/467 200mm long




Crile-Wood needle holder

tungsten carbide insert, box joint, satin finish, straight

SC-GI/217 159mm long

SC-GI/218 178mm long


SC-GI/198 150mm long



Derf needle holder

tungsten carbide insert

SC-GI/BUPHF002 115mm long




Finochietto needle holder

tungsten carbide jaws

SC-GI/222 215mm long


SC-GI/223 265mm long



Halsey needle holder

tungsten carbide insert,
box joint, satin finish, straight

SC-GI/220 127mm long



Kilner needle holder

tungsten carbide insert


SC-GI/221 140mm long

with two gold bows for identification **TUNGSTEN CARBIDE NEEDLE HOLDERS**

Masing needle holder

tungsten carbide insert,
one angled bow

SC-GI/009R 130mm long



Mayo-Hegar needle holder

tungsten carbide insert,
box joint, satin finish, straight


SC-GI/472 130mm long

SC-GI/214 159mm long

SC-GI/215 178mm long

SC-GI/216 209mm long


LONGER SIZE AVAILABLE



Nievert needle holder

tungsten carbide insert,
box joint, satin finish, straight

SC-GI/219 127mm long



Ryder needle holder

tungsten carbide insert, box joint, satin finish, straight

0.2mm jaw pitch

SC-GI/SH-5197a 135mm long

SC-GI/SH-5197b 155mm long

SC-GI/SH-5197c 175mm long

SC-GI/SH-5197d 210mm long


0.4mm jaw pitch

SC-GI/SH-5198a 135mm long

SC-GI/SH-5198b 155mm long

SC-GI/SH-5198c 175mm long


SC-GI/SH-5198d 210mm long




Sarot needle holder

diamond jaw, tungsten carbide, gold plated bows

SC-GI/451 181mm long



COMBINED SCISSOR AND NEEDLE HOLDERS




Olsen-Hegar scissor/needle holder

screw joint

SC-GI/226 140mm long

SC-GI/227 165mm long

SC-GI/228 190mm long




Olsen-Hegar scissor/needle holder

with tungsten

carbide insert

SC-GI/229 147mm long


SC-GI/230 165mm long



Gillies scissor/needle holder

right hand, with tungsten carbide insert

SC-GI/235 158mm long




Gillies scissor/needle holder

158mm long

SC-GI/231 right hand

SC-GI/232 left hand




Gillies scissor/needle holder

heavy pattern

158mm long

SC-GI/233 right hand

SC-GI/234 left hand



Foster Gillies scissor/needle holder

fine pattern


SC-GI/236 127mm long

BAG AND TOWEL CLIPS

Spring clip

serrated pattern, stainless steel,
for fixing slings to leg splints,
also for towel clipping

SC-GI/344 63mm wide



Spring bag clip

serrated


SC-GI/345 26mm wide



Jones towel clip

superior quality
spring forcep type


SC-GI/346 89mm long



Shardles towel clip

spring curved, cross action

SC-GI/347 77mm long




Backhaus towel clip

superior quality, box joint

SC-GI/348 89mm long

SC-GI/348a 102mm long


SC-GI/349 130mm long



Moynihan towel clip

screw joint, 2-2 teeth, for sides
of incision

SC-GI/470 191mm long




Ball and socket towel forceps

box joint


SC-GI/350 89mm long

SC-GI/351 115mm long

SC-GI/352 140mm long




RETRACTORS



Goligher retractor
SC-GI/305

complete with 2 centre blades:
1 x 90mm wide x 50mm deep
1 x 90mm wide x 90mm deep
and 2 pairs of lateral blades:
1 pair 25mm wide x 45mm deep
1 pair 50mm wide x 90mm deep



Miles retractor
SC-GI/376 340mm long

complete with 2 centre blades:
1 x 90mm wide x 50mm deep
1 x 90mm wide x 75mm deep
and 2 pairs of lateral blades:
1 pair 80mm wide x 35mm deep
1 pair 80mm wide x 75mm deep
other blade sizes available


RETRACTORS

Balfour abdominal retractor

bivalve, baby

SC-GI/372

fully open 89mm



Balfour abdominal retractor

trivalve, frame 203mm

SC-GI/373


child, 102mm opening, 51mm deep lateral blades, 38mm deep centre blade

SC-GI/374

adult, 178mm opening, 64mm deep lateral blades, 44mm deep centre blade

SC-GI/375


extra large, 178mm opening, 102mm deep lateral blades, 108mm deep centre blade




Berkeley abdominal retractor

parallel action

SC-GI/461 blades 35mm x 60mm deep, 215mm long



RETRACTORS




Collin retractor

110mm fully open, 235mm long

SC-GI/322 complete with centre blade
55 x 40mm, and one pair of blades 38mm x 60mm

SC-GI/322a without centre blade




Gelpi retractor

self-retaining

SC-GY/013a 140mm long


SC-GY/013b 178mm long



Gossetts abdominal retractor

SC-GI/303 large, fully open 185mm long

SC-GI/304 small, fully open 95mm long




Joll thyroid retractor

self retaining, with lock-in nut


SC-GI/363 152mm long

RETRACTORS

Pozzi abdominal retractor



SC-GI/370
swivel blades: 38 x 60mm,
bivalve, 254mm long




SC-GI/371
swivel blades: 38 x 59mm,
trivalve, 203mm long

Bowlby dilator

curved


SC-GI/495 145mm long



Norfolk and Norwich retractor

self retaining, with internal rack,
4 x 5 teeth, 51mm blades,
216mm long


SC-GI/325a 4 x 5 teeth,
51mm blades




Irwin retractor

self retaining, with cam rack, 4 x 4 teeth

SC-GI/324 178mm long




RETRACTORS



Travers retractor

self retaining with internal rack, 210mm long

SC-GI/323 4 x 5 teeth, 95mm fully open




West retractor

self retaining, 140mm long

SC-GI/326 3 x 4 teeth, 60mm fully open

SC-GI/326cr 3 x 4 teeth, 60mm fully open, with cam rack



Weitlander retractor

3-4 teeth

135mm long

SC-GI/496 blunt

SC-GI/496a sharp

SC-GI/496b sharp, 2-3 teeth

SC-GI/496c blunt, 2-3 teeth

165mm long

SC-GI/497 blunt

SC-GI/497a sharp



Kocher retractor

230mm long

SC-GI/378 36 x 20mm blade


250mm long

SC-GI/378a 67 x 25mm blade

SC-GI/378b 80 x 42mm blade

SC-GI/378c 80 x 50mm blade

SC-GI/378d 80 x 60mm blade



Kocher Langenbeck retractor

215mm long

SC-GI/336sp 25 x 6mm blade

SC-GI/379a 35 x 8mm blade

SC-GI/379b 35 x 11mm blade

SC-GI/379c 35 x 15mm blade

SC-GI/357sp 41 x 11mm blade

SC-THO/cutsp015 55 x 11mm blade

SC-GI/379d 70 x 15mm blade

SC-GI/379e 80 x 12mm blade


SC-GI/379f 80 x 16mm blade

RETRACTORS

Langenbeck retractor


single ended

- SC-GI/336** blade 22mm x 6mm
- SC-GI/337** blade 45mm x 13mm
- SC-GI/337a** blade 65mm x 13mm
- SC-GI/338a** blade 75mm x 25mm
- SC-GI/338a1** blade 65mm x 25mm



London retractor


SC-GI/494 230mm



Canny Ryall navy retractor

double-ended, with one single hook blade 5mm wide, and one fenestrated blade


- SC-GI/343** fenestrated blade 16mm wide, 190mm long
- SC-GI/343a** fenestrated blade 19mm wide, 190mm long
- SC-GI/343b** fenestrated blade 25mm wide, 220mm long
- SC-GI/343c** fenestrated blade 30mm wide, 220mm long



Canny Ryall army retractor

double prong, 13mm wide,
165mm long


- SC-GI/342** large, 19mm wide
- SC-GI/342a** small, 11mm wide




Czerny retractor

double ended

SC-GI/341 171mm long



RETRACTORS




Lane retractor

double ended blade
229mm long

SC-GI/339 13mm x 19mm wide

SC-GI/340 25mm x 32mm wide




Morris retractor

single blade, 229mm long

SC-GI/330 38 x 51mm

SC-GI/331 51 x 51mm

SC-GI/332 64 x 51mm




Morris retractor

double ended, 216mm long

SC-GI/327 32 x 51mm and 38 x 51mm

SC-GI/328 44 x 51mm and 51 x 51mm

SC-GI/329 64 x 51mm and 38 x 51mm



Moynihan retractor

SC-GI/462 single ended,
blade 13 x 22mm ,
215mm long



Deaver abdominal retractors

SC-GI/306 set of 5, 315mm overall length

SC-GI/307 25mm wide

SC-GI/308 25mm wide

SC-GI/309 38mm wide

SC-GI/310 51mm wide

SC-GI/311 76mm wide

RETRACTORS

Deaver abdominal retractors

with round handles

SC-GI/312 set of 5, 305mm overall length


SC-GI/313 25mm wide

SC-GI/314 25mm wide

SC-GI/315 38mm wide

SC-GI/316 51mm wide

SC-GI/317 76mm wide




Dyball and Lock retractor

for liver

SC-GI/335 64 x 144mm, 241mm long


OTHER SIZES AVAILABLE - PRICE ON APPLICATION



Elliot (modified Deaver) retractor

with flat handle,

SC-GI/308M




Kelly retractor blade

267mm long

SC-GI/333 38 x 191mm

SC-GI/334 57 x 191mm




RETRACTORS



Durham retractor
with fenestrated blade
216mm long


SC-GI/300 10mm wide
SC-GI/301 14mm wide
SC-GI/302 20mm wide



Ballance (Volkman) retractor
203mm long

sharp
SC-GI/463 1 prong
SC-GI/463a 2 prong
SC-GI/463b 3 prong
SC-GI/463c 4 prong
SC-GI/463d 6 prong

blunt
SC-GI/463e 1 prong
SC-GI/463f 2 prong
SC-GI/463g 3 prong
SC-GI/463h 4 prong
SC-GI/463i 6 prong



Volkman retractor

SC-OR/400 sharp, 220mm long
SC-OR/400a blunt, 220mm long
SC-OR/400b sharp, 245mm long



Kocher retractor

SC-OR/N036c sharp, 205mm long
SC-OR/N036 sharp, 220mm long
SC-OR/N036b blunt, 220mm long



Senn-Miller retractor
170mm long

SC-SH/5881



Single hook retractor
sharp

SC-GI/SH5671 152mm long



Gillies skin hook
solid handle, 190mm long

SC-PL/065a standard, fine
SC-PL/065b half size, fine

DISSECTORS AND PROBES

Macdonald dissector double ended

SC-GI/285 191mm long



Brodie director probe ended, silver plated

SC-GI/286 165mm long



Watson Cheyne dissector

fine, double-ended, with probe end

SC-GI/287 127mm long

SC-GI/288 191mm long

SC-GI/288a 230mm long

SC-GI/288b 255mm long



Probe, malleable silver plated

SC-GI/289a 127mm long

SC-GI/289 152mm long

SC-GI/290 178mm long

SC-GI/291 203mm long

SC-GI/292 255mm long

SC-GI/292a 305mm long



Probe, sinus silver

SC-GI/464 152mm long




Probe, fish tail silver plated

SC-GI/293 127mm long

SC-GI/294 152mm long

SC-GI/295 203mm long

SC-GI/296 254mm long



Probe, Jobson Horne dressing (aural)

with one flat tapered end, serrated edge and one small ring end
stainless steel

SC-GI/297 140mm long

SC-GI/298 178mm long




MISCELLANEOUS



Aneurysm needle


small, for infusions

SC-GI/472a 146mm long



Syme aneurysm needle

SC-GI/473 165mm long




Ring cutter

SC-GI/452 circular saw
blade

165mm long


SC-GI/452a spare blade



Mayo instrument pin

SC-GI/355 102mm long

SC-GI/356 152 mm long



Trocar and cannula

for hormone implant, 5mm bore

SC-GY/072 115mm long



Baron scalpel handle

130mm long

SC-PM/003 No 3

SC-PM/004 No 4



Scalpel handle

127mm long

SC-PM/007 No 3

SC-PM/008 No 4


SC-PM/008a No 4 L

SC-PM/001a No 7

HOLLOW WARE


Kidney dish
stainless steel

- SC-GI/380** 152mm
- SC-GI/381** 203mm
- SC-GI/382** 254mm
- SC-GI/383** 305mm




Lotion bowl
without lid, stainless steel

- SC-GI/384** 102mm
- SC-GI/385** 152mm
- SC-GI/386** 203mm
- SC-GI/387** 254mm
- SC-GI/388** 305mm




Instrument tray
stainless steel


- SC-GI/389** 254 x 203 x 50mm
- SC-GI/390** Lid only
- SC-GI/391** 305 x 254 x 50mm
- SC-GI/392** Lid only
- SC-GI/393** 356 x 305 x 50mm
- SC-GI/394** Lid only



HOLLOW WARE/SILICONE MATTING




Cheatles jar
SC-GI/397
stainless steel
203mm x 50mm



Measuring jug
SC-GI/395
1 litre capacity




Petri dish
SC-SH/002149
100mm, stainless steel,
with lid




Gallipot
SC-GI/396
stainless steel
63mm

Cushion, protect, and organise delicate surgical instruments during storage, transportation, and sterilisation. Perforations allow for ventilation during the sterilisation process. The silicone mat can be cut to fit any size container.


Silicone matting Medical grade for storage, sterilisation and transportation of instruments



BES/1004M
502mm long x 242mm wide
7mm holes, 10mm between holes



BES/1001M
515mm long x 240mm wide
3mm holes, 6mm between holes




BES/1002M
570mm long x 550 wide
4mm holes, 13mm between holes

silicone matting available for all baskets - please ask for price **WIRE PRODUCTS - STERILISING BASKETS**

Din basket

wire mesh base with perforated sides and two carrying handles




SC-GI/CO006 230 x 130 x 60mm
SC-GI/CO007 240 x 130 x 50mm
SC-GI/CO008 240 x 130 x 60mm
SC-GI/CO009 240 x 250 x 50mm
SC-GI/CO010 240 x 250 x 60mm
SC-GI/CO014 300 x 300 x 50mm

SC-GI/CO015 350 x 350 x 60mm
SC-GI/CO011 480 x 250 x 50mm
SC-GI/CO012 480 x 250 x 60mm
SC-GI/CO013 480 x 250 x 70mm
SC-GI/CO016 480 x 400 x 60mm
SC-GI/CO017 480 x 480 x 50mm

SC-GI/CO018 480 x 480 x 60mm
SC-GI/CO019 500 x 300 x 50mm
SC-GI/CO020 550 x 430 x 50mm
SC-GI/CO021 600 x 300 x 60mm

Din basket

wire mesh basket with two carrying handles



SC-GI/CO022 240 x 250 x 50mm
SC-GI/CO023 240 x 250 x 60mm
SC-GI/CO024 240 x 250 x 70mm
SC-GI/CO031 300 x 300 x 60mm


SC-GI/CO025 480 x 250 x 50mm
SC-GI/CO026 480 x 250 x 60mm
SC-GI/CO027 480 x 250 x 70mm
SC-GI/CO028 480 x 250 x 100mm

SC-GI/CO029 480 x 480 x 50mm
SC-GI/CO030 480 x 480 x 70mm


Din basket lids

for perforated and wire mesh baskets

SC-GI/CO032 240 x 250mm
SC-GI/CO033 300 x 300mm
SC-GI/CO034 480 x 250mm
SC-GI/CO035 480 x 480mm



WIRE PRODUCTS - STERILISING BASKETS silicone matting available for all baskets - please ask for price



Arthroscope basket

wire mesh base with two fixing clamps

SC-GI/CO001 290 x 80 x 52mm

Endoscope basket

with small lens basket


SC-GI/CO002 6400 x 150 x 77mm

wire mesh base with two fixing clamps

SC-GI/CO003 S - 460 x 80 x 52mm

SC-GI/CO004 D - 480 x 80 x 52mm

SC-GI/CO005 D - 680 x 80 x 52mm




ENT basket

wire mesh base with perforated sides and two carrying handles, supplied with lid

SC-GI/CO036 480 x 250 x 160mm

SC-GI/CO037 fixing clamps



Fine mesh basket

with hinged take off / lockable lid

SC-GI/CO041 120 x 155 x 40mm

SC-GI/CO044 120 x 155 x 50mm

SC-GI/CO042 235 x 155 x 40mm

SC-GI/CO045 235 x 155 x 50mm


SC-GI/CO043 375 x 260 x 40mm

with separate lid

SC-GI/CO038 120 x 155 x 40mm

SC-GI/CO039 235 x 155 x 40mm

SC-GI/CO040 375 x 260 x 40mm






Sieve

SC-GI/5071 203mm x 50mmø

OPTICAL FIBRE CABLES

3.5mm ø		4.8mm ø	
SC-LL/001SM	1.8 metres long	SC-LL/002SM	1.8 metres long
SC-LL/001	2.3 metres long	SC-LL/002	2.3 metres long
SC-LL/001XL	3 metres	SC-LL/002XL	3 metres long

light source end



Basic raw material	Silicate
Core to clad ratio	80:20
Numerical aperture	0.54 (65°) 0.66 (82°)
Wavelength range for max transmission	380 - 1500nm
Index	step index
Length	up to 30 feet
Max temperature	250° and upwards

TO ORDER

1 SELECT THE DIAMETER AND LENGTH REQUIRED,
2 ADD THE LETTER FOR THE LIGHT SOURCE END,
FOLLOWED BY THE LETTER FOR THE INSTRUMENT END.

eg to order a 4.8mm x 3metre Olympus to Olympus
quote code SC-LL/002XL/d/o

OPTICAL FIBRE CABLES

instrumental end




j
AFI, ASPEN, CABOT, CODMAN, CUDA, EDER, STRYKER, ZIMMER



k
FEMALE



L
DESIGNS FOR VISION HEADLIGHT



m
LUXTEC HEADLIGHT




n
PILLING INSTRUMENT



o
KARL STORZ OLYMPUS



p
WOLF




q
BFW-WEHMERLITE


angled end



r
NO ANGLE



s
45° ANGLE



t
90° ANGLE

TO ORDER

1 SELECT THE DIAMETER AND LENGTH REQUIRED,
2 ADD THE LETTER FOR THE LIGHT SOURCE END,
FOLLOWED BY THE LETTER FOR THE INSTRUMENT END.

eg to order a 4.8mm x 3metre Olympus to Olympus
quote code SC-LL/002XL/d/o