




































# 13. Protective clothing

APPEARANCE	FEATURES	INTENDED USE & BENEFIT
 <p>09368</p>	<ol style="list-style-type: none"> <li><b>Heat protective glove</b> <ul style="list-style-type: none"> <li>Made of silicone</li> <li>Washable in washer-disinfector</li> <li>Sterilizable (heat resistant up to 240°C)</li> <li>REF: 09368</li> </ul> </li> </ol>	<ol style="list-style-type: none"> <li><b>Heat protective glove</b> <ul style="list-style-type: none"> <li>For secure unloading of steam sterilizers or washer-disinfectors</li> </ul> </li> </ol>
 <p>09455</p>	<ol style="list-style-type: none"> <li><b>Disposable nitrile examination &amp; reprocessing gloves</b> <ul style="list-style-type: none"> <li>powder-free, latex-free</li> <li>high puncture- and cut-resistance (EN388) / high chemical protection (EN374)</li> <li>Total length: 40cm</li> <li>Various sizes (S-XL)</li> <li>REF: 09455 et seqq.</li> </ul> </li> </ol>	<ol style="list-style-type: none"> <li><b>Disposable nitrile examination &amp; reprocessing gloves</b> <ul style="list-style-type: none"> <li>Extra long 40cm gloves provide more security during examinations and the cleaning of surgical instruments</li> <li>Good grip on wet products due to rough fingertips - excellent tactility</li> </ul> </li> </ol>
 <p>09401 / 09402</p>	<ol style="list-style-type: none"> <li><b>Disposable protective gowns</b> <ul style="list-style-type: none"> <li>Made of reinforced (front and arms) non-woven material</li> <li>Water resistant / breathable</li> <li>size: universal</li> <li>length: approx. 130-136 cm</li> <li>colour: bright yellow / orange</li> <li>material: 30 g/m<sup>2</sup> (yellow) / 15 g/m<sup>2</sup> (orange)</li> <li>REF: 09401 / 09402</li> </ul> </li> </ol>	<ol style="list-style-type: none"> <li><b>Disposable protective gowns</b> <ul style="list-style-type: none"> <li>The Interlock gowns are designed to protect the user against contaminations</li> <li>The gowns cover the entire front and back of the upper part of the body</li> </ul> </li> </ol>
 <p>09081 / 09082</p>	<ol style="list-style-type: none"> <li><b>Disposable face protection</b> <ul style="list-style-type: none"> <li>Various options (Goggles / Face Shields / Mouth-Eye protection)</li> <li>Sizes: universal</li> <li>REF: 09081 / 09082 / 09417 / 09440 / 09441</li> </ul> </li> </ol>	<ol style="list-style-type: none"> <li><b>Disposable face protection</b> <ul style="list-style-type: none"> <li>Goggles are designed to be used during examinations or cleaning of instruments</li> <li>Face Shields are designed to be used when a protection of the complete face is needed (e.g. for cleaning heavily soiled instruments, working with high risk patients etc.)</li> <li>Mouth-Eye protection is used in when a protection of the complete face is needed and furthermore the environment shall be protected against spittle and bacteria</li> </ul> </li> </ol>
 <p>09417</p>		

## INTERLOCK USP

- Silicone heat protective gloves offer a long durability at low costs
- Wide variety of protective clothing for CSSDs and examinations available
- Economical prices