


Symmetry Surgical® Vein Strippers



Symmetry Surgical's single-use **Vein Strippers** are designed to provide proven, cost effective removal of saphenous varicose veins. Our vein strippers bring a 30-year heritage of quality and eliminate the need to clean and process them after use.

Each box of Symmetry Surgical Vein Strippers includes:

- 12 sterile packs per box
- Each individual sterile pack contains:
 - o Nylon, x-ray detectable components
 - o 2 x 91cm (36") cables
 - o 1 x handle
 - o 1 x 9.5mm ϕ Olive
 - o 1 x 12.7mm ϕ Olive
 - o 1 x 15.9mm ϕ Olive



Sterile, Each Package Includes Three Olives (S, M, L), One Handle and Two Stripper Cables

CODE	DESCRIPTION	LENGTH
63-4031	Disposable Vein Stripper	3/8 in 10 mm

Distributed in the UK by


surgicalholdings

7-9 Parkside Centre, Potters Way, Southend-on-Sea
Essex SS2 5SJ T 01702 602050
E info@surgicalholdings.co.uk



symmetry surgical®

Instruments | Sterilization Solutions | Retractors | Energy

©2013 Symmetry Surgical, Inc. All rights reserved. Cat No. SM0043