






13. Protective clothing

APPEARANCE	FEATURES	INTENDED USE & BENEFIT
 <p>09368</p>	<ol style="list-style-type: none"> Heat protective glove <ul style="list-style-type: none"> Made of silicone Washable in washer-disinfector Sterilizable (heat resistant up to 240°C) REF: 09368 	<ol style="list-style-type: none"> Heat protective glove <ul style="list-style-type: none"> For secure unloading of steam sterilizers or washer-disinfectors
 <p>09455</p>	<ol style="list-style-type: none"> Disposable nitrile examination & reprocessing gloves <ul style="list-style-type: none"> powder-free, latex-free high puncture- and cut-resistance (EN388) / high chemical protection (EN374) Total length: 40cm Various sizes (S-XL) REF: 09455 et seqq. 	<ol style="list-style-type: none"> Disposable nitrile examination & reprocessing gloves <ul style="list-style-type: none"> Extra long 40cm gloves provide more security during examinations and the cleaning of surgical instruments Good grip on wet products due to rough fingertips - excellent tactility
 <p>09401 / 09402</p>	<ol style="list-style-type: none"> Disposable protective gowns <ul style="list-style-type: none"> Made of reinforced (front and arms) non-woven material Water resistant / breathable size: universal length: approx. 130-136 cm colour: bright yellow / orange material: 30 g/m² (yellow) / 15 g/m² (orange) REF: 09401 / 09402 	<ol style="list-style-type: none"> Disposable protective gowns <ul style="list-style-type: none"> The Interlock gowns are designed to protect the user against contaminations The gowns cover the entire front and back of the upper part of the body
 <p>09081 / 09082</p>	<ol style="list-style-type: none"> Disposable face protection <ul style="list-style-type: none"> Various options (Goggles / Face Shields / Mouth-Eye protection) Sizes: universal REF: 09081 / 09082 / 09417 / 09440 / 09441 	<ol style="list-style-type: none"> Disposable face protection <ul style="list-style-type: none"> Goggles are designed to be used during examinations or cleaning of instruments Face Shields are designed to be used when a protection of the complete face is needed (e.g. for cleaning heavily soiled instruments, working with high risk patients etc.) Mouth-Eye protection is used in when a protection of the complete face is needed and furthermore the environment shall be protected against spittle and bacteria
 <p>09417</p>		

INTERLOCK USP

- Silicone heat protective gloves offer a long durability at low costs
- Wide variety of protective clothing for CSSDs and examinations available
- Economical prices