

What is new in our range of Laparoscopic Instruments and Devices...

① New Jaw Parts

» H00-533-631



» H00-533-721



» H00-533-732



» H00-533-822



» H00-533-672



» H00-533-708



» H00-533-722



» H00-533-954



Some jaw parts have been optimized in jaw length and others have been developed completely new – for better and safer application.

(see Laparoscopy catalogue/pages 14-31)

② Extended information – e.g. approx. value of jaw length

③ Extended standards – Instrument length 450 mm

» now all scissors, forceps, needle holders and electrodes are available with length 450mm

(only applicable on $\varnothing 5$ / $\varnothing 10$ mm)

④ Needle holders – extended range



- » handles – stainless steel or aluminium anodized
 - straight or curved type
- » flushing port Luer Lock
- » different types of jaw parts

(see Laparoscopy catalogue/page 34)

5 Trocar sleeves – new models



» optimized slim version of stainless steel trocar sleeve



» housing made of PPSU
» silicone valve

These new version are both reduced in weight and improved in working length – for better handling.

Shafts are with color code for easy and fast assembling.

(see Laparoscopy catalogue/pages 38-39)

6 Cases for scopes and instruments



- » storage / sterilization trays for endoscopes and light cable – plastic or stainless steel available
- » Instrument holding kit – with special frame for laparoscopic instruments

(see Laparoscopy catalogue/page 59)