

Bovie Medical Corporation's electrosurgery
business is now a part of



Suite of **DERMATOLOGY PRACTICE PRODUCTS** *designed for your portfolio of procedures*

Precision, versatility, and safety
working seamlessly together

Expand patient care through more
in-office procedures that provide
diagnostic and treatable solutions

Full suite of products offered in
one convenient place



Bovie

Quality Surgical Products
procedure-focused design

The Next Generation of High Frequency Desiccators and Electrosurgical Generators

HIGH FREQUENCY DESICCATORS & ELECTROSURGICAL GENERATORS



Unsurpassed 4-year warranty on all Bovie Electrosurgery Generators!



A942 – Bovie Derm 942

Feature: User-friendly and fewer accessories needed for additional functionality.

Benefit: Saves time and money during procedures and less accessories to purchase.

Feature: Adjust power from 0.1 watt – 40 watts.

Benefit: Higher resolution in bi-polar mode.

Feature: 3-Button Pencil works with or without optional footswitch.

Benefit: Cost-effective and user has complete control including up/down power from the handpiece as well as activation.

Feature: Digital Error Detection

Benefit: Continuously monitors every aspect of the unit's output. At the sign of any problem, the machine instantly disables the output and displays the appropriate error code on the display. No grounding pad needed!



A952 – Bovie Bantam | PRO

Feature: Designed and engineered to be the most reliable and durable generator available today.

Benefit: Unit's life is prolonged by incorporating the latest in chip technology and reducing the number of internal parts.

Feature: 50 watts of power

Benefit: No grounding pad needed when in desiccation mode.

Feature: High frequency desiccator with cutting ability.

Benefit: Provides a 2-in-one generator system offering 50 watts of cutting power and coag capabilities.

Feature: 3-Button Handpiece (included)

Benefit: Allows the surgeon to adjust the power output up and down as well as activate the unit, even when in foot control mode.

Feature: Split Pad Bovie NEM Technology

Benefit: Complete patient safety to prevent burns during procedures.

Perfect option for MOHS procedures



A1250S – Specialist | PRO

Feature: The A1250S Specialist|PRO is a very affordable solution for most in-office electrosurgery procedures.

Benefit: The surgeon receives a much faster return on their investment.

Feature: Cut and coag activation is on the hand piece.

Benefit: The surgeon can move seamlessly from cut to coag without leaving the sterile field.

Feature: The push-buttons controls allow easier use of the generator.

Benefit: The push-button front saves time and motion in adjusting the generator. The smoothly embossed push-button front allows the staff to clean the generator more quickly since there are no knobs.

Feature: Split Pad Bovie NEM Technology

Benefit: Complete patient safety to prevent burns during procedures.

Perfect option for MOHS procedures

COMPARISON CHART

FEATURES FOR:	A942	A952	A1250S
Cut (# of modes)	N/A	✓ 50 W	✓ 120 W (1)
Blend (# of steps)	N/A	✓ 50 W	✓ 90 W (1)
Coag	✓ 40 W	✓ 50 W	✓ 80 W (1)
Fulguration/Spray Coag	✓ 40 W	✓ 40 W	✓ 40 W
Standard Bipolar	✓ 40 W	✓ 40 W	✓ 30 W

FEATURES FOR:	A942	A952	A1250S
Pad Alarm	N/A	✓	✓
Pad Sensing (tissue impedance)	N/A	✓ NEM	✓
ACCESSORIES/OTHER FOR:	A942	A952	A1250S
Foot Control	✓ A803	✓ A803	✓
Pencils/Handpiece	✓ A902	✓ A902	✓ ESPR2

Illuminating the Path to Better Healthcare

MEDICAL LIGHTING

Bovie offers energy efficient LED Medical Lights to provide 50,000 hours of cool, clear visibility and are ideal for patient exams or small surgical procedure rooms. Control your lighting options with one hand. Bovie Medical Lighting provides a multi-dimming and illumination field that is controlled with a sterilizable handle keeping you in a sterile field.



MI-550



MI-750



MI-1000

Recommended Complementary Products

SMOKE EVACUATION SYSTEM



Smoke Shark™ II

The Bovie Smoke Shark II is designed to be easy to use and is effective for controlling surgical smoke generated from the use of any energy devices including electrosurgical, laser, and ultrasonic products.



SERS2 – EZ link remote activation switch for Smoke Shark II

DIAGNOSE AND TREAT



DermaScout™ II

The DermaScout II is an efficient and dependable instrument for visualization of dermal lesions.



High-Temperature Cauteries

The world leader in battery-operated cauteries.

CRYOSURGERY



IceOut™

Bovie IceOut

IceOut cryosurgery system is patient friendly and convenient to use while providing efficient and safe patient outcomes in 40 seconds or less. Perfect option to treat various skin lesions including warts and skin tags and provides 60+ treatments per system.



Bovie Freezpoint™

Cryosurgery is a popular treatment for simple skin lesions that can be performed by physicians, in-office. The Bovie Freezpoint is a safer, easier, and colder device, that delivers an improved, low cost per procedure, cryosurgical option.

Recommended Complementary Products

PENS AND SHEATHS



A902 – Three-button power control & autoclavable handpiece



ESP1 – Sterile disposable electro-surgical push button pencil



A910ST – Disposable handpiece and cord sheath, sterile



A910 – Disposable handpiece sheath, non-sterile

BIPOLAR ACCESSORIES



A824 – Reusable bipolar forceps



A827V – Bipolar cable



A803 – Foot switch

SPECIALTY ELECTRODES



Tungsten wire loop electrodes



ES61 – Supercut tungsten needle electrode



ES20 – Ball electrode



ES02 – Needle electrode

DERMAL TIP ELECTRODES



A804 (non-sterile) & A805 (sterile) – Sharp dermal tip



A806 (non-sterile) & A807 (sterile) – Blunt dermal tip



A806DE (non-sterile) & A807DE (sterile) – Derm-Elite premium blunt dermal tip



H10012 (non-sterile) & H10008 (sterile) – H-type sharp tip disposable electrode



H10112 (non-sterile) & H10108 (sterile) – H-type blunt tip disposable electrode