


# Bovie®

Operating with precision™

## INTERNATIONAL PRODUCT CATALOGUE

  
**surgical**holdings

UK Distributor  
T 01702 602050  
E [info@surgicalholdings.co.uk](mailto:info@surgicalholdings.co.uk)



# BOVIE – THE INVENTOR OF ELECTROSURGERY

The first use of the Bovie® electrosurgical generator in an operating room was on October 1, 1926, at Peter Bent Brigham Hospital in Boston, Massachusetts. A few days earlier (September 28, 1926) Dr. Cushing had been unsuccessful in removing a mass from a patient's head due to its highly vascular nature. In Dr. Cushing's notes from October 1st, "...with Dr. Bovie's help I proceeded to take off most satisfactorily the remaining portion of tumor with practically none of the bleeding which was occasioned in the preceding operation."

Over seventy-five years later, we continue to improve the work of Dr. Bovie. The following pages represent our product line, including four Bovie electrosurgical generators and three high frequency desiccators. The work continues today, but our sincere thanks go out to Dr. Cushing and William T. Bovie for their incredible invention.

Bovie Medical Corporation is a public company on the NYSE Amex Market and can be found under the symbol BVX.

**Bovie Medical Corporation**  
Clearwater, FL 33760-4004, USA  
Int'l Phone +1 727 384 2323 • Fax +1 727 347 9144  
[www.boviemedical.com](http://www.boviemedical.com) • [sales@boviemed.com](mailto:sales@boviemed.com)

**Investor Relations Contact:**  
Bovie Medical Corporation  
[investor.relations@boviemed.com](mailto:investor.relations@boviemed.com)

**EU Authorized Representative:**  
Emergo Europe  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands



RoHS Compliant



## ELECTROSURGERY

IDS-310 and OR   PRO 300 .....	4
IDS-210 and Surgi-Center   PRO .....	5
Specialist   PRO .....	6
Specialist   PRO-G .....	7
Bantam   PRO™ .....	8
Bantam   PRO-G .....	9
DERM 942 .....	10
DERM 101 & DERM 102 .....	11
Electrosurgery Pencils .....	12
Suction Coagulators .....	13
Disposable Electrodes, Blades, Balls, Needles, Dermal Tips .....	14
Resistick II™ Coated Electrodes .....	15
Laparoscopic Electrodes .....	16
Disposable Loops .....	17
Tungsten SuperCut™ Needles & Loop, Arthroscopic Electrodes .....	18
Reusable Electrodes, Blades, Balls, Needles .....	19
Reusable Loops .....	20
Reusable Bipolar Forceps .....	21
Mobile Stands, Footswitches, Wall Mount Kits .....	22
Return Electrodes, Active Cable, Cords, Adapters .....	23

## SMOKE

Smoke Shark™ II .....	24
Smoke Evacuation Accessories .....	25

## COLPOSCOPES, ACUDOP™ II & CLEARSPEC®

Colpo-Master™ I .....	26
Colpo-Master II.....	27
Colpo-Master III.....	28
AcuDop II .....	29
ClearSpec .....	30

## CRYOSURGERY & DERMASCOUT™ II

Freezpoint™ .....	31
IceOut™ .....	31
DermaScout II .....	32

## CAUTERIES

High-Temperature Cauteries.....	33
Low-Temperature Cauteries.....	34
Change-A-Tip™ Cauteries .....	34
Replacement Tips.....	35

## SPECIALTY LIGHTS

Flexible Lights, Stylets .....	36
Specialty Lights, Penlights, UV Woods .....	37

## OTHER PRODUCTS

Burr Handles, Burr Tips, Eye Bubble .....	38
Neuro-Pulse™ .....	38



## IDS-310 AND BOVIE® OR I PRO 300 ELECTROSURGICAL GENERATOR

The IDS-310 and Bovie OR I PRO 300 systems were designed in conjunction with OR clinicians to provide safe, user-friendly generators to meet all operating room procedure performance demands. These innovative energy-based devices offers 13 modes of surgical energy delivery. Special modalities include the advanced Bovie Bipolar, Auto Bipolar options and special Laparoscopic modes for surgeons to tailor all electrosurgery needs for their variety of surgical applications and user preferences.



### Style # IDS-310, A3350

#### ▼ Products Included

IDS-310 or Bovie OR I PRO 300 unit  
Hospital-grade power cord  
(120 VAC, 240 VAC)\*  
User's guide (MC-55-230-001)\*

#### Dimensions

Height: 16.5cm (6.5 in.)  
Width: 37.5cm (14.75 in.)  
Depth: 46cm (18.1 in.)  
Weight: <9.07 kg (<20 lbs)

#### Output Power

Cut I & II: 300 watts  
Blend (1-4): 200 watts  
Pinpoint Coagulation: 120 watts  
Spray Coagulation: 120 watts  
Gentle Coagulation: 120 watts  
Macro Bipolar: 150 watts  
Micro Bipolar: 80 watts  
Standard Bipolar: 50 watts  
Bovie Bipolar: 225  
Auto Bipolar: 80

#### Output Frequency

490 kHz ± 5

#### Line Voltage

100 – 240V~ ± 10%

#### Line Frequency

50 – 60 Hz

#### Warranty: **FOUR YEARS**

\* Not Pictured  
† Sold Separately

### ★ RECOMMENDED ACCESSORIES †

To avoid incompatibility and unsafe operation, we recommend using the following accessories with the IDS-310 and Bovie OR I PRO 300:

 BV-1253B	 BV-1254B	 A1255A	 A827V	 BV-IDS-CS
 ESRE	 A821	 ESP1, ESP6, or ESPR2		

## IDS-210 AND BOVIE® SURGI-CENTER I PRO ELECTROSURGICAL GENERATOR

The IDS-210 and Bovie Surgi-Center I PRO Electrosurgical Generators are intended for cutting, coagulation, ablation of tissue in general, gynecologic, orthopedic, ENT and urological procedures performed in an operating suite and procedure room. These units were designed to provide electrosurgery solutions capable of most operating room demands in affordable 200 watt generators. These units offer 12 options of energy delivery for a variety of surgical application; 2 cut, 4 blend, 3 coagulation and 3 bipolar modalities.



### Style # IDS-210, A2350

#### ▼ Products Included

IDS-210 or Bovie Surgi-Center I PRO unit  
Hospital-grade power cord  
(120 VAC, 240 VAC)\*  
User's guide (MC-55-232-001)\*

#### Dimensions

Height: 16.5cm (6.5 in.)  
Width: 37.5cm (14.75 in.)  
Depth: 46cm (18.1 in.)  
Weight: <9.07 kg (<20 lbs)

#### Output Power

Cut I & II: 200 watts  
Blend (1-4): 200 watts  
Pinpoint Coagulation: 120 watts  
Spray Coagulation: 120 watts  
Gentle Coagulation: 120 watts  
Macro Bipolar: 80 watts  
Micro Bipolar: 80 watts  
Standard Bipolar: 50 watts

#### Output Frequency

490 kHz  $\pm$  5

#### Line Voltage

100 – 240V  $\sim$   $\pm$  10%

#### Line Frequency

50 – 60 Hz

#### Warranty: **FOUR YEARS**

\* Not Pictured  
† Sold Separately


### ★ RECOMMENDED ACCESSORIES †

To avoid incompatibility and unsafe operation, we recommend using the following accessories with the IDS-210 and Bovie Surgi-Center I PRO:

 BV-1253B	 BV-1254B	 A1255A	 A827V	 BV-IDS-CS
 ESRE	 A821	 ESP1, ESP6, or ESPR2		

## SPECIALIST | PRO HIGH FREQUENCY ELECTROSURGICAL GENERATOR

The Specialist | PRO is the affordable solution for the busy physician office or specialty facility. Capable of most monopolar and bipolar surgical procedures and packed with dependable safety features, the Specialist | PRO satisfies the physician's surgical demands with safety, reliability and convenience. Paired with the Smoke Shark™ II Smoke Evacuator, Bovie® has developed a total system solution for procedures and facilities that require smoke evacuation.



### Style # A1250S

#### Dimensions

Height: 15.2cm (6 in.)  
Width: 26cm (10.25 in.)  
Depth: 30.5cm (12 in.)  
Weight: < 6.5 kg (< 14 lbs.)

#### Output Power

Cut: 120 watts  
Blend: 90 watts  
Coagulation: 80 watts  
Fulguration: 40 watts  
Bipolar: 30 watts

#### Output Frequency

350 – 600 kHz

#### Line Voltage

A1250U: 100-240 VAC

#### Line Frequency

50 – 60 Hz

**Warranty:**  
**FOUR YEARS**

### ★ RECOMMENDED ACCESSORIES†



BV-1253B



A1255A

### ▼ PRODUCTS INCLUDED



ESP1



ES01, ES02, ES20



ESRE



A1252C

### ▼ Products Included

**A1250S** (1) electrosurgical generator  
**ESP1** (1) sterile disposable pencil  
**ES01** (1) blade sterile disposable electrode  
**ES02** (1) needle sterile disposable electrode  
**ES20** (1) ball sterile disposable electrode  
**ESRE** (5) disposable split return electrodes  
**A1252C** (1) reusable grounding cord  
110 VAC Hospital-grade power cord\*  
220 VAC Hospital-grade power cord\*  
User's guide (**MC-55-237-001**)\*

\* Not Pictured

† Sold Separately

## SPECIALIST I PRO-G OB/GYN TOTAL SYSTEM SOLUTION

The A1250S-G combines the Specialist I PRO package with a smoke evacuation system and the latest mobile stand technology to produce a compact, mobile system for the office or specialty facility. As an added benefit, the A1250-G includes a loop electrode kit (5 of our most popular loops), 5 extended ball electrodes and a smoke evacuator hose reducer fitting for attaching to a vaginal speculum.

### ▼ PRODUCTS INCLUDED

 <p>A1250S</p>	 <p>ESP1</p>		
	 <p>ES01, ES02, ES20</p>	 <p>ESRE</p>	 <p>A1252C</p>
 <p>ESMS2</p>	 <p>ESLK</p>	 <p>SE02</p>	
	 <p>ES07</p>	 <p>786T</p>	 <p>SERF</p>

### Style # A1250S-G

#### ▼ Products Included

##### • A1250S electrocautery generator

- + **ESP1** (1) disposable pencil
- + **ES01** (1) blade sterile disposable electrode
- + **ES02** (1) needle sterile disposable electrode
- + **ES20** (1) ball sterile disposable electrode
- + **ESRE** (5) disposable split return electrodes
- + **A1252C** (1) reusable grounding cord
- + 110 VAC Hospital-grade power cord\*
- + 220 VAC Hospital-grade power cord\*
- + User's guide (**MC-55-237-001**)\*

##### • Smoke Shark™ II system

- + **SE02** (1) Smoke Shark with pneumatic footswitch
- + **SF35** (1) 35-hour filter
- + **786T** (1) 7/8" (22.225mm) tube 6' (182.9cm) long
- + **SERF** (1) reducer fitting (attaches to vaginal speculum)

##### • ESMS2 electrocautery mobile stand

##### • ESLK electrocautery loop kit

- + **ES09** (1) 10mm x 10mm loop
- + **ES11** (1) 15mm x 10mm loop
- + **ES12** (1) 20mm x 8mm loop
- + **ES13** (1) 20mm x 15mm loop
- + **ES16** (1) 10mm x 8mm square

##### • ES07 (5) extended ball electrodes

\* Not Pictured

† Sold Separately



**BOVIE® BANTAM | PRO™**

The Bovie Bantam | PRO, A952 is a High Frequency Desiccator with cutting ability. You might call it a split personality. One half of the unit is a high frequency desiccator, an upgrade to our very popular Bovie A940. The other half is a Bovie Generator, giving you 50 watts of cutting power and coag capabilities – plenty of power for most cutting procedures performed in the doctor's office.

**▼ Products Included**

- A952** (1) electrosurgical generator
  - A902** (1) handpiece
  - A804** (5) non-sterile sharp dermal tips
  - A805** (2) sterile sharp dermal tips
  - A806** (5) non-sterile blunt dermal tips
  - A807** (2) sterile blunt dermal tips
  - A806DE** (2) non-sterile DERM-Elite™ Premium blunt dermal tips
  - A807DE** (2) sterile DERM-Elite Premium blunt dermal tips
  - ES01** (3) blade sterile disposable electrodes
  - ES02** (1) needle sterile disposable electrodes
  - ES20** (1) ball sterile disposable electrodes
  - A910** (2) non-sterile disposable handpiece sheaths
  - A910ST** (2) sterile disposable handpiece sheaths
  - A1252C** (1) reusable grounding cord\*
  - ESRE** (5) disposable split grounding pads\*
  - A837** (1) wall mount kit
- 110 VAC Hospital-grade power cord\*  
User's guide (**MC-55-238-001**)\*

\*Not Pictured

**Style # A952****Dimensions**

Height: 188mm (7.40 in.)  
Width: 228mm (8.98 in.)  
Depth: 105mm (4.13 in.)  
Weight: < 2.26 kg (< 5 lbs.)

**Output Power**

Cut: 50 watts  
Blend: 50 watts  
Coagulation: 50 watts  
Fulguration: 40 watts  
Bipolar: 40 watts  
Micro Bipolar: 40 watts

**Output Frequency:** 304-451kHz

**Line Frequency:** 50Hz

**Line Voltage:** 100-240 VAC

**Warranty:** FOUR YEARS

**▼ PRODUCTS INCLUDED**

					
A902	A804, A805, A806, A807, A806DE, A807DE	ES01, ES02, ES20	A837	A910	A910ST

## BOVIE® BANTAM I PRO™ A952-G

The A952-G combines the Bovie Bantam I PRO A952 with the smoke evacuation system and a complete mobile stand to produce a compact, mobile system for the office. As an added benefit, the A952-G includes a loop electrode kit ( 5 of our most popular loops), 5 extended ball electrodes and a smoke evacuation hose reducer fitting for attaching to a vaginal speculum.

### ▼ PRODUCTS INCLUDED

				
	A902	ES01, ES02, ES20	A804, A805, A806, A807, A806DE, A807DE	
				
	A1252C	ESRE	A910ST	A910
				
	ESLK	A837	SE02	
				
	ES07	786T		SERF

### Style # A952-G

#### ▼ Products Included

##### • A952 electrocautery generator

- + **A902** (1) handpiece
- + **A804** (5) non-sterile sharp dermal tips
- + **A805** (2) sterile sharp dermal tips
- + **A806** (5) non-sterile blunt dermal tips
- + **A807** (2) sterile blunt dermal tips
- + **A806DE** (2) non-sterile DERM-Elite™ Premium blunt dermal tips

- + **A807DE** (2) sterile DERM-Elite Premium blunt dermal tips
- + **ES01** (3) blade sterile disposable electrodes
- + **ES02** (1) needle sterile disposable electrode
- + **ES20** (1) ball sterile disposable electrode
- + **A910** (2) non-sterile disposable handpiece sheaths\*
- + **A910ST** (2) sterile disposable handpiece sheaths
- + **A1252C** (1) reusable grounding cord
- + **ESRE** (5) disposable split grounding pads
- + **A837** (1) wall mount kit
- + 110 VAC Hospital-grade power cord\*
- + User's guide (**MC-55-238-001**)\*

##### • Smoke Shark™ II system

- + **SE02** (1) Smoke Shark with pneumatic footswitch
- + **SF35** (1) 35-hour filter
- + **786T** (1) 7/8" (22.23mm) tube 6' (182.9cm) long
- + **SERF** (1) reducer fitting\* (attaches to vaginal speculum)

##### • A812-C electrocautery mobile stand

- + **A812** (1) mobile stand
- + **A812-E** (1) height extender
- + **A808-T** (1) instrument tray
- + **A812-BT** (1) bottom tray

##### • ESLK electrocautery loop kit

- + **ES09** (1) 10mm x 10mm loop
- + **ES11** (1) 15mm x 10mm loop
- + **ES12** (1) 20mm x 8mm loop
- + **ES13** (1) 20mm x 15mm loop
- + **ES16** (1) 10mm x 8mm square

##### • ES07 (5) extended ball electrodes

\* Not Pictured

**BOVIE® DERM 942**

The DERM 942 High Frequency Desiccator and its patent pending handpiece are designed and engineered to be the most reliable and durable desiccator available today. The unit's life is prolonged by incorporating the latest in chip technology and reducing the number of internal parts, thereby generating less internal heat.

**▼ Products Included**

- A942** (1) high frequency desiccator
- A902** (1) 3-Button autoclavable handpiece
- A804** (5) non-sterile sharp dermal tips
- A805** (2) sterile sharp dermal tips
- A806** (5) non-sterile blunt dermal tips
- A807** (2) sterile blunt dermal tips
- A806DE** (2) non-sterile DERM-Elite™ Premium blunt dermal tips
- A807DE** (2) sterile DERM-Elite Premium blunt dermal tips
- A910ST** (2) sterile disposable handpiece sheath
- A910** (2) non-sterile disposable handpiece sheath
- A837** (1) wall mount kit  
110 VAC Hospital-grade power cord\*  
User's guide (**MC-55-239-001**)\*

\*Not Pictured

**Style # A942****Dimensions**

Height: 188mm (7.40 in.)  
Width: 228mm (8.98 in.)  
Depth: 105mm (4.13 in.)  
Weight: < 2.26 kg (< 5 lbs.)

**Output Power**

To 40 watts (0-10 watts in  
1/10th watt increments on  
both monopolar and bipolar)

**Line Voltage**

Universal Power Supply  
100 to 240 VAC

**Output Frequency:** 328-492kHz

**Line Frequency:** 50-60Hz

**Line Voltage:** 100-240 VAC

**Warranty:** **FOUR YEARS**

**▼ PRODUCTS INCLUDED**

A902	A804, A805, A806, A807, A806DE, A807DE	A837	A910	A910ST

## DERM 101/102 HIGH FREQUENCY DESICCATOR

The DERM 101 and DERM 102 are economical, 10 watt high frequency desiccators, thoughtfully designed to give practitioners the ability to perform in-office, minor skin procedures such as the removal of skin tags, benign lesions and premalignant lesions. The DERM 101 High Frequency Desiccator offers ten watts of monopolar power, adjustable in 1/10th watt increments, to desiccate and fulgurate tissue. The power can be adjusted using the 3-Button Power Control Handpiece (A902) or using the rotary dial on the unit. The unit features a large, bright, LED display for excellent visibility across the procedure room. For the practitioner utilizing bipolar electrosurgery, Bovie® offers the DERM 102 High Frequency Desiccator with 10 watts of bipolar power, adjustable in 1/10th watt increments.



### Style # DERM 101 & DERM 102

#### Dimensions

Height: 208mm (8.2 in.)  
 Width: 184mm (7.2 in.)  
 Depth: 69mm (2.7 in.)  
 Weight: < 1.4 kg (< 3 lbs.)

#### Output Power

Up to 10 watts (0-10 watts in 1/10th watt increments on both monopolar and bipolar)

#### Main Voltage

Universal Power Supply  
 DERM 101/102: 100 to 240 VAC±10%

#### Main Frequency: 50 – 60 Hz

#### Main Current: 320mA Max.

#### Warranty: **FOUR YEARS**



### ▼ PRODUCTS INCLUDED

A902	A804, A805	A806, A807
A806DE, A807DE	A910	A837

### ▼ Products Included

#### DERM 101 or DERM 102 unit

- A902** (1) three-button power control autoclavable handpiece
- A910** (2) non-sterile, disposable handpiece sheaths
- A804** (5) non-sterile single-use sharp dermal tips
- A805** (2) sterile single-use sharp dermal tips
- A806** (5) non-sterile single-use blunt dermal tips
- A807** (2) sterile single-use blunt dermal tips
- A806DE** (2) non-sterile DERM-Elite™ Premium blunt dermal tips
- A807DE** (2) sterile DERM-Elite Premium blunt dermal tips
- A837** (1) wall mount kit

110 VAC medical-grade power cord\*

User's guide CD (**MC-55-229-001**)\*


\* Not Pictured



## ELECTROSURGERY PENCILS

Most Bovie® electrosurgery pencils utilize a 3-pin industry-standard connection and are compatible with most electrosurgical generators. The disposable pencils feature standard push button, rocker switch or foot-control (1-pin connection) activation. The pencils come in a variety of kit configurations which include holsters, scratch pads, stainless steel electrodes or Resistick II™ Coated Electrodes. All disposable pencils are packaged individually, sterile. Additionally, Bovie offers two reusable pencils with 3-pin connection for high watt units. The ESPR2 fits standard 3/32" electrodes and is validated for up to 50 autoclave cycles. The ESPR-EL fits special 1/16" electrodes and is validated for up to 25 autoclave cycles.

### REUSABLE PENCILS



#### ESPR2 (Top 50)\*

Reusable Push Button Pencil (Top 50) - 1/box, Non-Sterile  
(electrode not included)

#### ESPR-EL (Top 25)\*

Reusable Push Button Pencil w/ 1/16" collet - 1/box, Non-Sterile  
(electrode not included)

\*Pencils use a standard 3-prong connector



### DISPOSABLE PENCILS



ESP1\*

#### ESP1 PUSH-BUTTON CONFIGURATIONS

**ESP1** pencil, blade - 50/box  
**ESP1H** pencil, blade, holster - 40/box  
**ESP1HS** pencil, blade, holster, scratch pad - 40/box  
**ESP1N** pencil, needle - 50/box  
**ESP1HN** pencil, needle, holster - 40/box

#### RESISTICK II™ PENCIL CONFIGURATIONS

**ESP1T** disposable pencil, ES01T coated blade - 50/box  
**ESP1TN** disposable pencil, ES02T coated needle - 50/box  
**ESP1TH** disposable pencil, ESPH holster, ES01 coated blade - 40/box



ESP1T\*



ESP1TN\*



ESP6\*

#### ESP6 ROCKER SWITCH CONFIGURATIONS

**ESP6** pencil, blade - 50/box  
**ESP6H** pencil, blade, holster - 40/box  
**ESP6HS** pencil, blade, holster, scratch pad - 40/box  
**ESP6N** pencil, needle - 50/box  
**ESP6HN** pencil, needle, holster - 40/box

#### RESISTICK II PENCIL CONFIGURATIONS


**ESP6T** disposable pencil, ES01T coated blade - 50/box  
**ESP6TN** disposable pencil, ES02T coated needle - 50/box



ESP6T\*



ESP6TN\*



ESP7†

#### ESP7 FOOT-CONTROL CONFIGURATIONS



**ESP7** pencil, blade - 50/box  
**ESP7H** pencil, blade, holster - 40/box

† Pencils use a single pin connector



**ELECTROSURGERY PENCILS** Continued...

For low watt generators, Bovie® offers the proprietary A901 and A902 autoclavable, 3-button, power-control pencils.

	<p><b>A901</b> autoclavable reusable 3-button pencil for A950, A940, &amp; A900 - comes with generator or 1/box</p>
	<p><b>A902</b> autoclavable reusable 3-button pencil for A952, A942, DERM 102, DERM 101 - comes with generator or 1/box</p>

**ADDITIONAL ACCESSORIES**

Disposable handpiece sheaths are used in conjunction with electrosurgery pencils to provide a cost effective alternative to continuous sterilization. The use of handpiece sheaths significantly increases the life of the reusable pencils.







	<p><b>A910</b> disposable handpiece sheath, non-sterile - 100/box</p>
	<p><b>A910ST</b> disposable handpiece and cord sheath, sterile, 38" (96.52cm) - 25/box</p>
	<p><b>ESPH</b> disposable pencil holster, sterile - 40/box</p>
	<p><b>ESSP</b> disposable scratch pad - 40/box</p>

**SUCTION COAGULATORS**

Suction coagulators allow fluid evacuation with or without simultaneous coagulation. Suction coagulators are compatible with Bovie generators, as well as, those of most other companies. The devices feature three different diameters, 6" flexible shaft and hand or foot activation. The single use devices are individually packaged sterile, ten per box.

**HANDSWITCH SUCTION COAGULATORS**

Catalog Number	Description	Tube Diameter (French Size)	Cable Length
<b>SCH08</b>	Disposable	08 Fr	3m
<b>SCH10</b>	Disposable	10 Fr	3m
<b>SCH12</b>	Disposable	12 Fr	3m

			
French Gauge	12	10	8
Diameter (mm)	4	3.3	2.7
Diameter (inches)	0.158	0.131	0.105
			

**FOOTSWITCH SUCTION COAGULATORS**






Catalog Number	Description	Tube Diameter (French Size)	Cable Length
<b>SCF08</b>	Disposable	08 Fr	3m
<b>SCF10</b>	Disposable	10 Fr	3m
<b>SCF12</b>	Disposable	12 Fr	3m




## DISPOSABLE ELECTRODES

All disposable electrodes utilize the standard 3/32" (2.3mm) stainless steel shaft. The molded electrodes feature safety grip insulators combining patient and user safety with easy insertion into and removal from the surgical pencil. Disposable electrodes are manufactured to the highest standards. They come individually packaged sterile\*.





\* non-sterile and not packaged individually.

	<b>A804*</b> Non-Sterile <b>A805</b> Sterile
<b>Sharp Dermal Tip for use with A952, A942, DERM 102, DERM 101, A950, A940, A900, &amp; A800.</b> Non-Sterile - 100/box Sterile - 50/box	
	<b>A806*</b> Non-Sterile <b>A807</b> Sterile
<b>Blunt Dermal Tip for use with A952, A942, DERM 102, DERM 101, A950, A940, A900, &amp; A800.</b> Non-Sterile - 100/box Sterile - 50/box	
	<b>A806DE*</b> Non-Sterile <b>A807DE</b> Sterile
<b>Derm-Elite Disposable Dermal Tip for use with A952, A942, DERM 102, DERM 101, A950, A940, A900, &amp; A800.</b> Non-Sterile - 100/box Sterile - 50/box	
	<b>H10012*</b> Non-Sterile <b>H10008</b> Sterile
<b>H-type Sharp Tip Disposable Electrode for use with A952, A942, DERM 102, DERM 101, A950, A940, A900, &amp; A800.</b> Non-Sterile - 100/box Sterile - 50/box	
	<b>H10112*</b> Non-Sterile <b>H10108</b> Sterile
<b>H-type Blunt Tip Disposable Electrode for use with A952, A942, DERM 102, DERM 101, A950, A940, A900, &amp; A800.</b> Non-Sterile - 100/box Sterile - 50/box	

	<b>ES01</b>
<b>Standard Blade Electrode</b> 2.75" (6.99cm) - 50/box	
	<b>ES37</b>
<b>Modified Blade Electrode</b> 2.75" (6.99cm) - 25/box	
	<b>ES18</b>
<b>Angled Blade 45° Electrode</b> 2.5" (6.35cm) - 5/box	
	<b>ES55</b>
<b>Extended Blade Electrode</b> 4" (10.16cm) - 25/box	
	<b>ES54</b>
<b>Extended Modified Blade Electrode</b> 4" (10.16cm) - 25/box	
	<b>ES04</b>
<b>Extended Blade Electrode</b> 6" (15.24cm) - 25/box	
	<b>ES39</b>
<b>Extended Modified Blade Electrode</b> 6" (15.24cm) - 25/box	

















	<b>ES50</b>
<b>Standard 1/8" Ball Electrode</b> 0.125" (3.18mm) - 25/box	
	<b>ES21</b>
<b>Standard 5/32" Ball Electrode</b> 0.156" (3.97mm) - 25/box	
	<b>ES20</b>
<b>Standard 3/16" Ball Electrode</b> 0.188" (4.78mm) - 25/box	
	<b>ES06</b>
<b>Extended Ball Electrode</b> 0.118" (3mm) - 5/box	
	<b>ES07</b>
<b>Extended Ball Electrode</b> 0.197" (5mm) - 5/box	


	<b>ES02</b>
<b>Standard Needle Electrode</b> 2.75" (6.99cm) - 25/box	
	<b>ES38</b>
<b>Modified Needle Electrode</b> 2.75" (6.99cm) - 25/box	
	<b>ES03</b>
<b>Extended Needle Electrode</b> 6" (15.24cm) - 25/box	
	<b>ES40</b>
<b>Extended Modified Needle Electrode</b> 6" (15.24cm) - 25/box	




## RESISTICK II™ COATED ELECTRODES




Resistick II coated blade, ball and needle electrodes are expertly coated with PTFE (polytetrafluoroethylene). Resisting eschar build-up the electrode moves more effectively through tissue increasing efficiency during surgical procedures. Simply use a damp sponge or cloth to almost effortlessly clean the debris off the Resistick II electrodes.

	ES58T
<b>Blade Electrode</b> 2.5" (6.35cm) - 12/box	
	ES59T
<b>Modified Blade Electrode</b> 2.5" (6.35cm) - 12/box	
	ES01T
<b>Standard Blade Electrode</b> 2.75" (6.99cm) - 12/box	
	ES37T
<b>Modified Blade Electrode</b> 2.75" (6.99cm) - 12/box	
	ES18T
<b>Angled Blade 45° Electrode</b> 2.75" (6.99cm) - 12/box	
	ES55T
<b>Extended Blade Electrode</b> 4" (10.16cm) - 12/box	
	ES54T
<b>Extended Modified Blade Electrode</b> 4" (10.16cm) - 12/box	
	ES04T
<b>Extended Blade Electrode</b> 6" (15.24cm) - 12/box	
	ES39T
<b>Extended Modified Blade Electrode</b> 6" (15.24cm) - 12/box	

	ES50T
<b>Standard 1/8" Ball Electrode</b> 0.125" (3.18mm) - 12/box	
	ES21T
<b>Standard 5/32" Ball Electrode</b> 0.156" (3.97mm) - 12/box	
	ES20T
<b>Standard 3/16" Ball Electrode</b> 0.188" (4.78mm) - 12/box	
	ES06T
<b>Extended Ball Electrode</b> 0.118" (3mm) - 12/box	
	ES07T
<b>Extended Ball Electrode</b> 0.197" (5mm) - 12/box	



	ES02T
<b>Standard Needle Electrode</b> 2.75" (6.99cm) - 12/box	
	ES38T
<b>Modified Needle Electrode</b> 2.75" (6.99cm) - 12/box	
	ES56T
<b>Extended Needle Electrode</b> 4" (10.16cm) - 12/box	
















	ES57T
<b>Extended Modified Needle Electrode</b> 4" (10.16cm) - 12/box	
	ES03T
<b>Extended Needle Electrode</b> 6" (15.24cm) - 12/box	
	ES40T
<b>Extended Modified Needle Electrode</b> 6" (15.24cm) - 12/box	



## LAPAROSCOPIC ELECTRODES



Bovie's® new Laparoscopic electrodes are used in minimally invasive procedures performed in outpatient facilities, surgery centers and operating rooms. The electrodes are 33cm long, with a 3/32" shaft, for insertion into standard electrosurgery pencils. Each electrode is also packaged with an optional 4mm shaft adapter for connection into 4mm active cables. The laparoscopic electrode tip styles come in a variety of popular configurations with standard stainless steel finish, as well as PTFE coated.



UNCOATED		UNCOATED		COATED		COATED	
	LJ07		LL06		LJ07BT		LL06BT
Laparoscopic J-Hook Wire Electrode with 4mm Adaptor, Sterile - 5/box		Laparoscopic Standard Blade Electrode with 4mm Adaptor, Sterile - 5/box		Coated Laparoscopic J-Hook Wire Electrode with 4mm Adaptor, Sterile - 5/box		Coated Laparoscopic Standard Blade Electrode with 4mm Adaptor, Sterile - 5/box	
	LH08		LC04		LH08BT		LC04BT
Laparoscopic L-Hook Wire Electrode with 4mm Adaptor, Sterile - 5/box		Laparoscopic Curved Blade Electrode with 4mm Adaptor, Sterile - 5/box		Coated Laparoscopic L-Hook Wire Electrode with 4mm Adaptor, Sterile - 5/box		Coated Laparoscopic Curved Blade Electrode with 4mm Adaptor, Sterile - 5/box	
	LH01		LB02		LH01BT		LB02BT
Laparoscopic L-Hook Heavy Wire Electrode with 4mm Adaptor, Sterile - 5/box		Laparoscopic Ball Electrode with 4mm Adaptor, Sterile - 5/box		Coated Laparoscopic L-Hook Heavy Wire Electrode with 4mm Adaptor, Sterile - 5/box		Coated Laparoscopic Ball Electrode with 4mm Adaptor, Sterile - 5/box	
	LN03		LS05		LN03BT		LS05BT
Laparoscopic Needle Electrode with 4mm Adaptor, Sterile - 5/box		Laparoscopic Spatula Electrode with 4mm Adaptor, Sterile - 5/box		Coated Laparoscopic Needle Electrode with 4mm Adaptor, Sterile - 5/box		Coated Laparoscopic Spatula Electrode with 4mm Adaptor, Sterile - 5/box	

## DISPOSABLE LOOP ELECTRODES

All Bovie® LLETZ Loop and Square electrodes feature tungsten wire for superior shape and integrity throughout the excision procedure. The electrodes use a standard 3/32" (2.3mm) shaft. All molded loops are packaged in a patented protective tray to prevent damage during shipping.



### LOOPS

**ES08**  
5mm x 5mm loop - 5/box

**ES09**  
10mm x 10mm loop - 5/box

**ES10**  
15mm x 8mm loop - 5/box

**ES11**  
15mm x 10mm loop - 5/box

**ES12**  
20mm x 8mm loop - 5/box

**ES13**  
20mm x 15mm loop - 5/box

**ES31**  
20mm x 20mm loop - 5/box

**ES42**  
20mm x 12mm loop - 5/box

**ES43**  
15mm x 15mm loop - 5/box

**ES44**  
15mm x 5mm loop - 5/box

**ES45**  
13mm x 13mm loop - 5/box

**ES46**  
10mm x 5mm loop - 5/box

**ES47**  
25mm x 10mm loop - 5/box

**ES49**  
20mm x 10mm loop - 5/box

**ES51**  
13mm x 8mm loop - 5/box

**ES52**  
10mm x 8mm loop - 5/box

**ES53**  
15mm x 12mm loop - 5/box

### SQUARES

**ES41**  
5mm x 10mm square - 5/box

**ES14**  
5mm x 5mm square - 5/box

**ES15**  
10mm x 4mm square - 5/box

**ES16**  
10mm x 8mm square - 5/box

## SPECIALTY ELECTRODES


Bovie® manufactures an expansive line of specialty electrodes for a variety of applications including; SuperCut™ Needle Electrodes for precise and delicate procedures, Dermal Loop Electrodes for dermatologic procedures, as well as Arthroscopic electrodes.

DISPOSABLE SHORT TUNGSTEN .008" WIRE LOOPS					DISPOSABLE SHORT TUNGSTEN .015" WIRE LOOPS				
<b>ES22-8</b>	1/4" Loop (6.35mm)	<b>ES23-8</b>	5/16" Loop (7.94mm)	<b>ES24-8</b>	3/8" Loop (9.53mm)	<b>ES25-8</b>	1/2" Loop (12.7mm)	<b>ES26-8</b>	3/4" Loop (19.05mm)
<b>ES22</b>	1/4" (6.35mm) short shaft loop, sterile - 5/box	<b>ES23</b>	5/16" (7.94mm) short shaft loop, sterile - 5/box	<b>ES24</b>	3/8" (9.53mm) short shaft loop, sterile - 5/box	<b>ES25</b>	1/2" (12.7mm) short shaft loop, sterile - 5/box	<b>ES26</b>	3/4" (19.05mm) short shaft loop, sterile - 5/box
<b>ES24-8</b>	3/8" (9.53mm) short shaft loop, sterile - 5/box	<b>ES25-8</b>	1/2" (12.7mm) short shaft loop, sterile - 5/box	<b>ES26-8</b>	3/4" (19.05mm) short shaft loop, sterile - 5/box				
<b>REUSABLE SHORT TUNGSTEN .015" WIRE LOOPS</b>									
<b>ES22R</b>	1/4" (6.35mm) short shaft loop, non-sterile - 1/box	<b>ES23R</b>	5/16" (7.94mm) short shaft loop, non-sterile - 1/box	<b>ES24R</b>	3/8" (9.53mm) short shaft loop, non-sterile - 1/box	<b>ES25R</b>	1/2" (12.7mm) short shaft loop, non-sterile - 1/box	<b>ES26R</b>	3/4" (19.05mm) short shaft loop, non-sterile - 1/box

SUPER CUT TUNGSTEN NEEDLES		ES62R	ES60R	ES63	ES62	ES61	ES60
<b>REUSABLE SUPER CUT TUNGSTEN NEEDLES</b>		<b>DISPOSABLE SUPER CUT TUNGSTEN NEEDLES</b>					
<b>ES60R</b>	needle super fine 3cm, non-sterile - 1/box	<b>ES60</b>	needle super fine 3cm, sterile - 5/box	<b>ES61</b>	modified needle super fine 3cm, sterile - 5/box	<b>ES62</b>	needle super fine 4.5cm, sterile - 5/box
<b>ES62R</b>	needle super fine 4.5cm, non-sterile - 1/box	<b>ES61HS</b>	modified fine needle with Heatshrink, sterile - 5/box	<b>ES63</b>	modified needle super fine 4.5cm, sterile - 5/box		
<b>DISPOSABLE SUPER CUT TUNGSTEN NEEDLES</b>							
<b>AR00</b>		hook 90, sterile - 5/box		<b>AR02</b>		Acromioplasty, sterile - 5/box	
<b>AR01</b>		hook 45, sterile - 5/box		<b>AR03</b>		Meniscectomy, sterile - 5/box	






## REUSABLE ELECTRODES

All reusable electrodes utilize the standard 3/32" (2.3mm) stainless steel shaft. The molded electrodes feature safety grip insulators combining patient and user safety with easy insertion and removal from the surgical pencil. Electrodes are 5 – 25 autoclave cycles, refer to the Instructions For Use in the package for proper reuse.



	A831
<b>Angled 3/32" Ball Electrode</b> 0.094" (2.38mm) - 1/box	
	A832M
<b>Standard 3/32" Ball Electrode</b> 0.094" (2.38mm) - 1/box	
	ES20R
<b>Standard 3/16" Ball Electrode</b> 0.188" (4.78mm) - 1/box	

	ES21R
<b>Standard 5/32" Ball Electrode</b> 0.156" (3.97mm) - 1/box	
	ES06R
<b>Extended Ball Electrode</b> 0.118" (3mm) - 1/box	
	ES07R
<b>Extended Ball Electrode</b> 0.197" (5mm) - 1/box	


	A833
<b>Straight Needle Electrode</b> - 1/box	
	A834
<b>Angled Fine Needle Electrode</b> - 1/box	
	A834T
<b>Angled Tungsten Fine Needle Electrode</b> - 1/box	
	ES02R
<b>Standard Needle Electrode</b> 2.75" (6.99cm) - 1/box	
	ES03R
<b>Extended Needle Electrode</b> 6" (15.24cm) - 1/box	

	A811
<b>Reusable Blunt Electrode</b> - 1/box	
	A830
<b>Angled Sharp Electrode</b> - 1/box	
	A836
<b>Extended Angled Sharp Electrode</b> - 1/box	
	ES01R
<b>Standard Blade Electrode</b> 2.75" (6.99cm) - 1/box	
	ES18R
<b>Angled Blade 45° Electrode</b> 2.75" (6.99cm) - 1/box	
	ES04R
<b>Extended Blade Electrode</b> 6" (15.24cm) - 1/box	



## REUSABLE LOOPS/SQUARES, ELECTRODES

All Bovie® LLETZ Loop and Square electrodes feature tungsten wire for superior shape and integrity throughout the excision procedure. The electrodes use a standard 3/32" (2.3mm) shaft. All molded loops are packaged in a patented protective tray to prevent damage during shipping. Electrodes are 5 – 25 autoclave cycles, refer to the Instructions For Use in the package for proper reuse.




### LOOPS

<b>ES08R</b> 5mm x 5mm loop - 1/box	<b>ES31R</b> 20mm x 20mm loop - 1/box	<b>ES47R</b> 25mm x 10mm loop - 1/box
<b>ES09R</b> 10mm x 10mm loop - 1/box	<b>ES42R</b> 20mm x 12mm loop - 1/box	<b>ES49R</b> 20mm x 10mm loop - 1/box
<b>ES10R</b> 15mm x 8mm loop - 1/box	<b>ES43R</b> 15mm x 15mm loop - 1/box	<b>ES51R</b> 13mm x 8mm loop - 1/box
<b>ES11R</b> 15mm x 10mm loop - 1/box	<b>ES44R</b> 15mm x 5mm loop - 1/box	<b>ES52R</b> 10mm x 8mm loop - 1/box
<b>ES12R</b> 20mm x 8mm loop - 1/box	<b>ES45R</b> 13mm x 13mm loop - 1/box	<b>ES53R</b> 15mm x 12mm loop - 1/box
<b>ES13R</b> 20mm x 15mm loop - 1/box	<b>ES46R</b> 10mm x 5mm loop - 1/box	

















### SQUARES

<b>ES41R</b> 5mm x 10mm square - 1/box	<b>ES14R</b> 5mm x 5mm square - 1/box	<b>ES15R</b> 10mm x 4mm square - 1/box	<b>ES16R</b> 10mm x 8mm square - 1/box
---	--	---	---



## REUSABLE BIPOLAR FORCEPS

Reusable Bipolar Forceps are available in micro (.5mm), 1.0mm, and 2.0 mm tip size and in lengths ranging from 3.5" to 8". The forceps are guaranteed for twenty autoclave cycles. Refer to the Instructions For Use in the package for proper reuse. Reusable bipolar cord, with quality silicone construction are available for all Bovie® generators.

				
<b>A820</b> 7.25" (184.15mm) bayonet, fine (1.0mm), smooth - 1/box	<b>A821</b> 8" (203.2mm) bayonet, medium (2.0mm), smooth - 1/box	<b>A822</b> 7" (177.8mm) straight, fine (1.0mm), smooth - 1/box		
				
<b>A823</b> 7" (177.8mm) curved, fine (1.0mm), smooth - 1/box	<b>A824</b> 4.5" (114.3mm) curved, micro (.5mm) - 1/box	<b>A825</b> 4.75" (120.65mm) straight, medium (2.0mm), smooth - 1/box	<b>A826</b> 4.5" (114.3mm) straight, micro (.5mm), smooth - 1/box	<b>A840</b> 4" (101.6mm) Tenzel, micro (.5mm), coated - 1/box
				
<b>A841</b> 3.5" (88.9mm) McPherson, curved, micro (.5mm), coated - 1/box	<b>A842</b> 3.5" (88.9mm) McPherson, straight, micro (.5mm), coated - 1/box	<b>A843</b> 4" (101.6mm) Tenzel, micro (.5mm), uncoated - 1/box	<b>A844</b> 3.5" (88.9mm) McPherson, curved, micro (.5mm), uncoated - 1/box	<b>A845</b> 3.5" (88.9mm) McPherson, straight, micro (.5mm), uncoated - 1/box
			<p>► Bipolar forceps require an electrosurgery footswitch for activation. See footswitch page for options.</p>	
<b>A827BP</b> forceps cord - 1/pouch fits A950, A940, & A900	<b>A827F</b> forceps cord - 1/pouch fits A1250U (only)	<b>A827V</b> forceps cord - 1/pouch fits IDS-310, A3350, IDS-210, A2350, DERM 102, DERM 101, A1250S, A952, A942, A3250, & A2250		








## STANDS, FOOTSWITCHES & WALL MOUNT KIT





Bovie® generators are compatible with the footswitches below, when footswitch activation is preferred. If mobility is a concern, select from a variety of mobile stands.





FOOTSWITCHES		STANDS AND WALL MOUNT			
	<b>A803</b> A952, A942, DERM 102, DERM 101, A940, & A900 footswitch		<b>A812-C</b> mobile stand for A952, A942, DERM 102, DERM 101, A950, A940, & A900 [includes bottom tray, top tray, clamp, and 12" (304.8 mm) height extension] 47" (1193.8 mm) high		<b>ESMS2</b> mobile stand for A1250S, A1250U, & A1250, 18.9" (483 mm) wide x 19.6" (500 mm) deep x 32" (814 mm) high
	<b>A1203 and A1203W (Waterproof)</b> A1200 & A950 footswitch		<b>A812</b> mobile stand for A952, A942, DERM 102, DERM 101, A950, A940, & A900, 35" (889.0 mm) high		<b>ESMS-B</b> bottom tray for ESMS 13" (330.2 mm) wide x 11.75" (298.5mm) deep
	<b>BV-1253B</b> IDS-310, A3350, IDS-210, A2350, A1250S, IDS-400, A3250, IDS-300, IDS-200, A2250, & A1250 footswitch		<b>A812-BT</b> bottom tray for A812 13" (330.2 mm) long x 11.75" (298.5 mm) deep		<b>BV-IDS-CS</b> mobile stand for IDS-310, A3350, IDS-210, A2350, A1250S, IDS-400, IDS-300, A3250, IDS-200, A2250, A1250U, & A1250, 16" wide x 20" deep x 35" high
	<b>BV-1253BGI</b> Icon Gi footswitch		<b>A808-T</b> top tray and clamp for A812 13.5" (342.9 mm) wide x 9.5" (241.3 mm) deep		<b>A813</b> table-top stainless steel stand for A952, A942, DERM 102, DERM 101, A1200, A950, A940, & A900
	<b>BV-1254B</b> IDS-310, A3350, IDS-210, A2350, IDS-400, A3250, IDS-300, IDS-200, & A2250 bipolar footswitch		<b>A812-E</b> 12" (304.8 mm) height extension for A812		<b>A837</b> wall mount kit for A952, A942, DERM 102, DERM 101, A950, A940, & A900

## RETURN ELECTRODE, ACTIVE CABLE, CORD AND ADAPTER

Return electrodes are utilized for a safe return of electrosurgical energy from the patient body. Bovie® disposable, single-use return electrodes are manufactured with a super adhesive hydrogel for secure patient application. Also available are reusable grounding plates used mostly in Animal Health procedures and reusable return electrode cables for connection of reusable and disposable pads to the ESU. In addition, Bovie provides active cables and adapter for the connection of a variety of compatible electrosurgery accessories.

REUSABLE RETURN ELECTRODE PLATES AND CORDS	
	<b>A802EU</b> metal return electrode plate for A942, DERM 102, DERM 101, A940, & A900 - 1/pouch
	<b>A1204</b> metal return electrode plate and cord for A950 - 1/pouch
	<b>A1204P</b> replacement metal return electrode plate for use with A1204C and A1254C - 1/pouch
	<b>A1204C</b> replacement cord for A950 - 1/pouch
	<b>A1254C</b> return electrode cord for IDS-310, A3350, IDS-210, A2350, A1250S, A952, A3250, A2250, A1250U, & A1250 - 1/pouch
	<b>A1202C</b> connecting cord for A1202 or ESRS return electrode for A950 - 1/pouch
	<b>A1252C</b> return electrode cord for ESRE and ESRS to IDS-310, A3350, IDS-210, A2350, A1250S, A952, A3250, A2250, A1250U, & A1250 - 1/pouch

PENCIL ADAPTERS	
	<b>A905B</b> pencil adapter for 1/32" (0.8mm) electrode — fits most ESU pencils
	<b>A905B-1</b> pencil adapter for 1/20" (1.25mm) electrode — fits most ESU pencils
FOOT-CONTROL ADAPTER	
	<b>A905A</b> adapter for connecting footswitch pencil for the A950, A940, and A900
	<b>A1255A</b> adapter plug for connecting footswitching pencil for the IDS-310, A3350, A2350, A1250S, A3250, A2250, A1250U, A1250, and most ESUs

DISPOSABLE RETURN ELECTRODE PADS	
	<b>A1202 or ESRS</b> solid return electrode A1202 - 5/pack ESRS - 50/box
	<b>ESRSC</b> solid adult return electrode w/2.8 m cable - 50/box
	<b>ESRE</b> split adult return electrode - 50/box
	<b>ESREC</b> split adult return electrode w/2.8 m cable - 50/box



## SMOKE SHARK™ II SMOKE EVACUATOR

The Bovie® Smoke Shark II is designed to be easy to use and extremely effective at smoke and particle removal at the operative site. While designed for use with all electro-surgical generators, the unit is also excellent for use in controlling laser and Harmonic Scalpel smoke plume.



### STYLE # SE02

#### ▼ Products Included

(SE02) unit with Long Life 35-hour Filter (SF35),  
pneumatic footswitch (FSSEP)  
power cord\*  
user's guide (MC-23088)\*

#### DIMENSIONS

**Height:** 15.2cm (6")  
**Width:** 27.9cm (11")  
**Depth:** 39.4cm (15.5")  
**Weight:** 4.3kg (10 lbs.)

#### LINE FREQUENCY

50 – 60 Hz

#### WARRANTY

two years

#### LINE VOLTAGE

SE02: 120 VAC

### TYPICAL SETUPS

#### ELECTROSURGERY

**SE02** Smoke Shark unit with pneumatic footswitch (FSSEP)  
**SF35** long-life 35-hour filter  
**786T** † 7/8" x 6' (22.3mm x 182.9cm) tube case of 24  
**SEAS**\* † arm stand

#### OB/GYN

**SE02** Smoke Shark unit with pneumatic footswitch (FSSEP)  
**SF35** long-life 35-hour filter  
**786T** † 7/8" x 6' (22.3mm X 182.9cm) tube case of 24 or  
**SERF**\* † reducer fitting 7/8" to 1/4" x 24" (22.225mm to 6.35mm x 609.6mm) [fits optional vaginal speculum]

#### LASER

**SE02** Smoke Shark unit with pneumatic footswitch (FSSEP)  
**SF35** long-life 35-hour filter  
**786T** † 7/8" x 6' (22.3mm X 182.9cm) tube case of 24 or  
**SELW**\* † 7/8" x 8" (22.3mm X 182.9cm) laser resistant wand 3/pak  
**SEAS**\* † arm stand

### ★ RECOMMENDED ACCESSORIES †

786T Tubing	SEPAT Pencil Adaptor	SERS2 Bovie EZ Link Remote Activation Switch

\* Not Pictured  
† Sold Separately

# ★ Smoke Evacuation Accessories

	<b>SF35</b> long-life 35-hour filter for SE02 1/box		<b>786T</b> <sup>†</sup> 7/8" x 6' tube, non-sterile 24/box		<b>SERS2</b> <sup>†</sup> Bovie® EZ Link remote switch activator 1/box
	<b>SF18</b> <sup>†</sup> 18-hour filter for SE01 model 1/box		<b>786TS</b> <sup>†</sup> 7/8" x 6' tube, sterile 24/box		<b>SERF</b> <sup>†</sup> reducer fitting 7/8" to 1/4" x 24" 10/box [fits optional vaginal speculum]
  Pencil Sold Separately <sup>†</sup>	<b>SEPA</b> <sup>†</sup> Smoke Evacuation Pencil Adapter 3/8" x 1' Hose with 3/8" to 3/8" Adapter 25 per case		<b>SETW</b> <sup>†</sup> 7/8" x 6' tube with wand & tip, non-sterile 24/box		<b>SERFS</b> <sup>†</sup> reducer fitting, sterile 10/box [fits optional vaginal speculum]
  Pencil Sold Separately <sup>†</sup>	<b>SEPAT</b> <sup>†</sup> Smoke Evacuation Pencil Adapter and Tubing- 3/8" x 4' Hose with 7/8" to 6' Tube 10 per case		<b>SETWS</b> <sup>†</sup> 7/8" x 6' tube with wand & tip, sterile 24/box		<b>SELW</b> <sup>†</sup> laser-resistant wand 7/8" x 8", non-sterile 3/box
	<b>FSSEP</b> Pneumatic footswitch included with SE01 & SE02 1/box		<b>SEVL</b> <sup>†</sup> 7/8" x 6' tubing reduced to 3/8" for ValleyLab & IEC pencil, sterile 5/box		<b>SEAS</b> <sup>†</sup> arm stand 1/box

## SMOKE EVACUATION

- An estimated 500,000 healthcare professionals are exposed to surgical smoke proven to carry viruses, CO and other chemical and biological hazards<sup>1</sup>. The Smoke Shark II minimizing exposure to toxic smoke plume during electrosurgery, laser and Harmonic® scalpel procedures.
- The Extended-Life Filter (SF35) offers up to 35 hours performance reducing procedure cost by over 30% when used at the lowest setting.
- The Bovie Smoke Evacuator is compatible with most electrosurgical generators. The unit can operate in conjunction with the optional Bovie Remote Activation Switch (SERS2) for maximized filter life and surgical efficiency.

<sup>1</sup> Barret, Dr. William L., Garber, Dr. Shawn M, "Surgical Smoke – A Review in Literature", Business Briefing: Global Surgery (2004)

HARMONIC® is a registered trademark of Ethicon

<sup>†</sup>Sold Separately

## COLPO-MASTER™ COLPOSCOPES

Bovie® introduces the new Colpo-Master line, a perfect addition to the already popular line of gynecologic products. The line of colposcopes offers superior, 3D, stereoscopic viewing with excellent clarity during gynecologic procedures where millimeters matter. The scopes are available with consistent and dependable zoom range including simple 3 point magnification (3.75x, 7.5x and 15x) or continuous zoom 4 – 27x magnification. For the EMR ready facility the Colpo-Master line is available with Trinocular camera-ready heads to provide quick and easy visual records of treatment. All Colpo-Master Colposcopes are backed by a 7 year limited warranty. Paired with the Bovie 1250-G Electrosurgery System with Smoke Evacuation for a complete GYN LEEP system.

### Colpo-Master I



Camera-Ready Trinocular head with camera mount

*Low center of gravity increases stability, while EZ-Roll locking casters provide smooth gliding movement between exam rooms.*

*CS-103LED includes 3-legged base as shown. Also available as CS-105LED with 5-legged base.*

Colpo-Master™ is a registered trademark of United Products and Instruments, Inc; DBA UNICO

### COLPO-MASTER I ZOOM SUSPENSION ARM COLPOSCOPES

The Colpo-Master is perfect for the GYN performing both exams and procedures. The swing-arm design allows easy access to the vaginal canal during surgical procedure.

Colpo-Master I features smooth, continuous zoom 4 – 27x magnification with 20x eyepieces and sharp focus for 3-D imaging. Available with Binocular or camera-ready Trinocular heads with 298mm working distance. The field of view (FOV) ranges from 58.5mm at 4x to 8.5mm at 27x. All models with your choice of 3 or 5 leg sturdy mobile bases.

Designed with your viewing comfort in mind; each model has dual diopter adjustments to independently adjust the focus of both Left and Right eyepieces to eliminate eyestrain. Camera-ready Trinocular models allow the user to view through both eyepieces while using the camera.

Colpo-Master I's Suspension Arm is rock-steady. To make head positioning easy it comes with locking casters and stainless steel handlebars, standard.

Complete with dust cover, power cord, manual, tool kit and 7 Year Limited Warranty, 1 year on electrical system (excludes bulbs and fuses). Optional accessories are available - please contact your Bovie Representative for additional information.

Item No.	Description
<b>CS-103LED</b>	Colpo-Master I Binocular Zoom Suspension-Arm Colposcope
<b>CS-103LED-E</b>	Colpo-Master I Binocular Zoom Suspension-Arm Colposcope, 220V
<b>CS-103T-LED</b>	Colpo-Master I Camera-Ready Zoom Swing-Arm Colposcope
<b>CS-103T-LED-E</b>	Colpo-Master I Camera-Ready Zoom Swing-Arm Colposcope, 220V
<b>CS-103T-DV</b>	Swing Arm Digital Video Package, includes CS-103T-LED and CS-8185
<b>CS-103T-SV</b>	Swing Arm S-Video Package, includes CS-103T-LED and CS-8021
<b>CS-103T-SVM</b>	Swing Arm S-Video Package with Monitor, includes CS-103T-LED, CS8021 and CS8022
<b>CS-103T-SDV</b>	Swing Arm S-Video Packaged with Digital Image Grabber and Monitor, includes CS-103T-LED, CS-8021, CS-8022 and CS-8023
<b>CS-105LED</b>	Same as CS-103LED above, with 5-legged base
<b>CS-105T-LED</b>	Same as CS-103T-LED above, with 5-legged base

\*These products may not be available in all markets.

## COLPO-MASTER™ II CENTER-POST MOUNT COLPOSCOPES\*

### Colpo-Master II

C205T-LED Camera-Ready Trinocular head with camera



Colpo-Master II is a compact Center-Post Colposcope with unique, continuous zoom magnification and sharp focusing for crisp 3-D imaging.

#### Features:

- **HEAD:** Full continuous zoom magnification range: 4 – 27x with 20x eyepieces. Binocular or video-ready Trinocular heads, with 298mm working distance. The field of view (FOV) ranges from 58.5mm at 4x to 8.5mm at 27x
- 3-D Stereoscopic Viewing
- Dual diopters allow user to independently focus each eyepiece to adjust for the natural vision differential between eyes for viewing comfort. Camera-ready Trinocular models allow the user to view through both eyepieces while using the camera.
- Sturdy, 5-legged, easy-to-manuever base slides in and under exam tables for optimal positioning
- Dual light intensity controls for both white and green LEDs
- Adjustable height center post with smooth suspension system
- S-Video and digital camera options available, trinocular version with built-in C-Mount Camera adapter

Complete with dust cover, power cord, manual, tool kit and 7 Year Limited Warranty, 1 year on electrical system (excludes bulbs and fuses). Optional accessories are available - please contact your Bovie® Representative for additional information.

Item No.	Description
<b>CS-205LED</b>	Colpo-Master II Binocular Zoom Center Post Colposcope
<b>CS-205LED-E</b>	Colpo-Master II Binocular Zoom Center Post Colposcope, 220V
<b>CS-205T-LED</b>	Colpo-Master II Camera-Ready Zoom Center Post Colposcope
<b>CS-205T-LED-E</b>	Colpo-Master II Camera-Ready Zoom Center Post Colposcope, 220V
<b>CS-205T-DV</b>	Center Post Digital Video Package, includes CS-205T-LED and CS-8185
<b>CS-205T-SV</b>	Center Post S-Video Package, includes CS-205T-LED and CS-8021
<b>CS-205T-SVM</b>	Center Post S-Video Package with Monitor, includes CS-205T-LED, CS-8021 and CS-8022
<b>CS-205T-SDV</b>	Center Post S-Video Packaged with Digital Image Grabber and Monitor, includes CS-205T-LED, CS-8021, CS-8022 and CS-8023

Colpo-Master™ is a registered trademark of United Products and Instruments, Inc; DBA UNICO

\*These products may not be available in all markets.

**Colpo-Master III****COLPO-MASTER™ III TRIPLE MAGNIFICATION COLPOSCOPES\***

Introducing the NEW Colpo-Master III Triple Magnification, LED Illuminated Center-Post Colposcopes. Select wheeled or pedestal base models.

Simple and easy to use, featuring brilliant, sharp optics and three magnification settings: 3.75x, 7.5x and 15x. Changing magnification is a simple rotation of the knob, while staying in focus!

The state of the art LED illumination provides bright, white light that will not burn out, or burn too hot!

**Features:**

- Three Magnification Settings: 3.75x, 7.5x and 15x
- Look straight through, no angled eyepieces
- Rheostat intensity controlled, bright LED illumination
- Paired 12.5x Widefield Eyepieces with rubber eye-guards
- 3-D Stereoscopic Viewing
- Dual Diopter, independent eye focusing controls
- Flip-down green filter, built into the head
- 300 mm Working Distance
- Extra-Wide 115 mm Interpupillary Distance Adjustment
- Center-Post stand with rolling, locking casters, or pedestal base
- Light weight, easy to maneuver
- Broad Range Universal Power Supply 90-240V
- Net Weight 23lbs, Shipping Weight 32lbs

Complete with dust cover, power cord, manual, tool kit and 7 Year Limited Warranty, 1 year on electrical system (excludes bulbs and fuses).

Item No.	Description
<b>CS-305</b>	Colpo-Master III Binocular Triple Magnification Colposcope with 5-Wheel Base
<b>CS-305E</b>	Colpo-Master III Binocular Triple Magnification Colposcope, 220V with 5-Wheel Base
<b>CS-8182</b>	Digital Eyepiece Camera, 2.0 megapixel with Image Capture software for Windows XP/7
<b>CS-8187</b>	HD Camera, 1080p with HDMI Output

\*These products may not be available in all markets.

Colpo-Master™ is a registered trademark of United Products and Instruments, Inc; DBA UNICO



## BOVIE® ACUDOP™ II DOPPLER SYSTEM\*

Fetal Doppler's have been used for years to hear and see the fetus heart rate. AcuDop Doppler system provides Unmatched Digital sound quality along with Breakthrough XD probe technology that provides up to a 56% increased sensitivity at shallow depths. Models available in Displayed, Non-Displayed and Rechargeable options. Interchangeable probes 2MHz, 3 MHz & waterproof. Convert your OB doppler to a vascular doppler by simply changing the probe to a 5MHz or a 8MHz.



- Unmatched Digital Sound Quality
- Interchangeable OB & Vascular Probes
- Breakthrough XD probes provide up to 56% increased sensitivity at shallow depths (2-5cm)
- Cost saving rechargeable models
- Sized to fit comfortably in your hand and pocket
- Only 5 year Warranty available on ALL Labor & ALL Parts



### VASCULAR PROBES

XD technology provides a clear high quality sound to determine patients vascular bloodflow.



### ACUDOP DOPPLER SYSTEMS

#### Style #

- AD-770\*\*** AcuDop II 770, Display Unit
- AD-770R** AcuDop II 770R, Rechargeable Unit
- AD-330\*\*** AcuDop II 330, Non-Display Unit
- AD-330R\*\*** AcuDop II 330, Rechargeable Non-Display Unit

\*\*Not Pictured

### OB PROBE CONFIGURATIONS

#### OB PROBES:

#### Style #

- A2** AcuDop II 2MHZ OB Probe
- A3** AcuDop II 3MHZ OB Probe

#### WATERPROOF OB PROBES:

#### Style #

- A2W** AcuDop II 2MHz Waterproof OB Probe
- A3W** AcuDop II 3MHz Waterproof OB Probe

### VASCULAR PROBE CONFIGURATIONS

#### Style #

- A5** AcuDop II 5MHZ Vascular Probe
- A8\*\*** AcuDop II 8MHZ Vascular Probe

\*\*Not Pictured

\*These products may not be available in all markets.

## CLEARSPEC® SPECULUM\*

The ClearSpec Speculum is the only speculum clinically proven to provide complete visualization of the cervix in 3 times as many patients.<sup>1</sup>




ClearSpec

### Style #

- 102** – Small Vaginal Speculum  
25/box
- 103** – Medium Vaginal Speculum  
25/box
- 104** – Large Vaginal Speculum  
25/box
- 105** – LED Illuminator  
12/box




Patent US 6,432,048 B1



105- LED Illuminator  
Patent Pending

- May help to improve the diagnosis by providing a more complete view of the cervix
- Improves workflow while decreasing procedure time and patient discomfort due to less need for speculum repositioning
- Reusable LED illuminator with 2-hours of light with an auto on/off option
- Illuminator conveniently fits into speculum handle providing light source with no cords or recharging batteries
- LED Illuminator is less expensive than traditional rechargeable lights that can accidentally be thrown away



<sup>1</sup>Hill et al., Sheathed versus standard speculum for visualization of the cervix. IJGO. 2014; 125:116-120

In this study the investigators were able to visualize the entire cervix in 42 (61.8%) patients when using the sheathed speculum [ClearSpec] compared with 11 (16.4%) patients undergoing standard speculum examination (p<0.0001).

\*These products may not be available in all markets.

## BOVIE® FREEZPOINT™

Cryosurgery is a popular treatment for simple skin lesion that can be performed by physicians, in-office. The Bovie Freezpoint is a safer, easier, and colder device, that delivers an improved, low cost per procedure, cryosurgical option. The Freezpoint device is 27% - 62% colder than common multi-use disposable systems. The slim tip design aids in clear visualization of the lesion for precise delivery of -89° C, non-flammable, Nitrous Oxide therapy. The outcome is a very safe and clinically effective treatment with reduced tissue damage to surrounding healthy skin.



### Style #

**CRY8** – 8 gram device

**CRY16** – 16 gram device

**CRY16-2\*** – 16 gram device (2 pack)  
\*not pictured

**Gas:** Nitrous Oxide (N<sub>2</sub>O)

**Temperature:** -89° C (-128° F)

**Shelf Life:** Two years and up to six months after activation

## BOVIE ICEOUT™ PORTABLE CRYOSURGICAL SYSTEM

Patient friendly, convenient to use and provides efficient and safe patient outcomes. Let IceOut help you expand your practice offerings.



### Features

- Results comparable to liquid nitrogen<sup>1</sup>
- Convenient and easy to use
- Effective cure rate
- Minimal treatment time for patients, 40 seconds or less
- Designed with the lowest (GWP) Global Warming Impact
- Reusable can with disposable applicators
- Estimated 50+ treatments per box system
- Low cost treatment with profitable reimbursement rates
- No Hazmat Fee

### Style #

**CRY2** – IceOut Portable Cryosurgical System, (2) 80ml aerosol dispensers and (60) 2mm applicators

**CRY5** – IceOut Portable Cryosurgical System, (2) 80ml aerosol dispensers and (50) 5mm applicators

**CRY25** – IceOut Portable Cryosurgical System, (2) 80ml aerosol dispensers and (60) 2mm & 5mm applicators

**Gas:** Dimethyl ether (DME) with propane / isobutane

**Temperature:** -55°C (-67°F)

**Shelf Life:** Three years – Refer to outer box for date

## DERMASCOUT™ II\*

The DermaScout II is an efficient and dependable instrument for visualization of dermal lesions. The dermatoscope is used to view, measure and monitor skin lesions, fungal infections tumors, warts and other condition that may affect the epidermis of the patient. The hand-held magnifier has 6x magnification, LED illumination, with a dual scale (Eng./Metric) measuring reticle as well as an independent diopter to allow the user to focus at differing focal points.



**Storage Case**



**LED illumination**

### Catalog #

**BVDS-6-02** – DermaScout II Dermatoscope

**DS-6-01** – Replacement Battery, Type CR2, 3













### Features

- Magnification with 40 mm top lens
- Diopter focus control for sharp focusing
- 40 mm field of view (flat field)
- 31 mm working distance
- LED illuminated with an array of 9 LED's
- Over 150 hours of battery power
- Reticule with dual measuring scales:  
millimeters and inches
- 20 mm: 10 left of center/10 right of center
- Measure in 0.1 mm increments
- Measure in 0.05 in. increments
- Durable locking transport and storage case
- Battery included (type: CR2, 3V)

\*These products may not be available in all markets.

## HIGH-TEMPERATURE CAUTERIES







Bovie® Medical has been the number one manufacturer of battery operated cauteries for decades. The new, safety-enhanced, snap-design cautery is the only battery-operated cautery in the world to address known safety hazards of improper disposal. The improved patent-pending cautery offers a safety-enhanced design, with a recessed button to reduce the potential for inadvertent activation. Bovie also modified the cautery units to allow the practitioner to “snap off” the end of the cautery for safe and effective removal of the battery power source.

HIGH-TEMPERATURE CAUTERIES				
 <p>← Recessed Button</p> <p>← Easy “Snap” Design</p> <p><b>■</b> New Style Cauteries</p>	 <p><b>AA01</b> high-temp fine tip, 4-yr. shelf life - 10/box</p>	 <p><b>AA17</b> high-temp fine tip extended 2” shaft, 4-yr. shelf life - 10/box</p>	 <p><b>■</b> Legacy Cauteries</p>	
	 <p><b>AA03</b> high-temp loop tip, 4-yr. shelf life - 10/box</p>	 <p><b>AA21</b> high-temp elongated tip, 4-yr. shelf life - 10/box</p>		
	 <p><b>AA05</b> high-temp loop tip extended 2” shaft, 4-yr. shelf life - 10/box</p>	 <p><b>AA25</b> high-temp fine tip micro, 4-yr. shelf life - 10/box</p>		
	 <p><b>AA09</b> high-temp loop tip extended 5” shaft, 4-yr. shelf life - 10/box</p>	 <p><b>AA27</b> high-temp loop tip micro, 4-yr. shelf life - 10/box</p>		
	 <p><b>AA11</b> high-temp adjustable fine tip, 4-yr. shelf life - 10/box</p>	 <p><b>AA29</b> high-temp elongated tip micro, 4-yr. shelf life - 10/box</p>		








## LOW-TEMPERATURE CAUTERIES

Bovie® Medical has been the number one manufacturing the battery operated cauteries for decades. The new, safety-enhanced, snap-design cautery has been enhanced based on valuable customer feedback. The improved cautery offers a new safety-enhanced design with a recessed button to reduce the potential for inadvertent activation. Bovie also modified the cautery units to allow the practitioner to “snap off” the end of the cautery for safe and effective removal of the battery power source. 10 per box.

LOW-TEMPERATURE CAUTERIES			
			
<span style="color: green;">■</span> <i>New Style Cauteries</i>		<span style="color: blue;">■</span> <i>Legacy Cauteries</i>	
<span style="color: green;">■</span>  <p><b>AA00</b> low-temp fine tip, 4-yr. shelf life - 10/box</p>	<span style="color: blue;">■</span>  <p><b>AA02</b> low-temp adjustable fine tip, 4-yr. shelf life - 10/box</p>	<span style="color: green;">■</span>  <p><b>AA04</b> low-temp elongated fine tip, 4-yr. shelf life - 10/box</p>	<span style="color: green;">■</span>  <p><b>AA90</b> low-temp micro fine tip, 4-yr. shelf life - 10/box</p>

## CHANGE-A-TIP™ CAUTERIES

These replaceable battery & replaceable tip cauteries are excellent battery operated cauteries where reusability is a must. These products are frequently utilized for the office, ER, and clinical procedures.









<span style="color: green;">■</span> <i>New Style Cauteries</i> 	<p><b>HIT0</b> Change-A-Tip low-temp handle (includes one H100 non-sterile tip and one “AA” alkaline battery) - 1/box</p>	<span style="color: green;">■</span> <i>New Style Cauteries</i> 	<p><b>HIT1</b> Change-A-Tip high-temp handle (includes one H101 non-sterile tip and two “AA” alkaline batteries) - 1/box</p>
<span style="color: green;">■</span> <i>New Style Cauteries</i> 	<p><b>DELO</b> Change-A-Tip deluxe <b>low-temp</b> cautery kit (includes one low-temperature handle, two H100 non-sterile tips, two H104 non-sterile tips, one “AA” battery and a foam-lined case) - 1/box</p>	<span style="color: green;">■</span> <i>New Style Cauteries</i> 	<p><b>DEL1</b> Change-A-Tip deluxe <b>high-temp</b> cautery kit (includes one high-temperature handle, two H101 non-sterile tips, two H103 non-sterile tips, two “AA” alkaline batteries and a foam-lined case) - 1/box</p>
<span style="color: green;">■</span> <i>New Style Cauteries</i> 	<p><b>DEL2</b> Change-A-Tip deluxe HI-LO cautery kit (includes <b>one low-temp</b> handle, <b>one high-temp</b> handle, one sterile H100 tip, one sterile H101 tip, one sterile H103 tip, one sterile H104 tip, six “AA” alkaline batteries and a foam-lined case) - 1/box</p>		






## CAUTERY REPLACEMENT TIPS

Sterile cautery tips include a sterile drape to cover the non-sterile handle. Tips are individually packaged 10 per box. The replacement tips are used in conjunction with the Change-A-Tip™ cautery handles as well as the deluxe kits.



Protective Tray

	<b>H100</b> low-temp fine tip - 10/box		<b>H105</b> high-temp 2" (50.8mm) flex loop tip - 10/box
	<b>H101</b> high-temp fine tip - 10/box		<b>H106</b> high-temp 2" (50.8mm) flex fine tip - 10/box
	<b>H103</b> high-temp loop tip - 10/box		<b>H109</b> high-temp 5" (127.0mm) loop tip - 10/box
	<b>H104</b> low-temp elongated fine tip - 10/box		<b>H112</b> high-temp 5" (127.0mm) fine tip - 10/box

	<b>H110</b> contra angle fine tip - 10/box		<b>H121</b> high-temp elongated tip - 10/box	 <i>Change-A-Tip sold separately</i>	<b>HISL</b> Sterile sheath for replaceable cauteries - 10/box
	<b>H111</b> high-temp circular tip - 10/box		<b>H101-ADH</b> fine tip w/sheath seal - 10/box		

Cauteries

## FLEXIBLE LIGHTS AND STYLETS

Bovie® manufactures two groups of flexible lighting products. The first is available in three different sizes 5" (127.0mm), 10" (254.0mm), and 15" (381.0mm) (refers to the flexible shaft portion of the light). These lights are for use in procedures where additional light is needed. The second group of flexible lights is our lighted stylets, for difficult intubations.

### FLEXIBLE LIGHTS



**NS10** Surch-Lite™ 10" (254.0mm) non-sterile - 1 ea



**ST05** Surch-Lite 5" (127.0mm) sterile - 10/box



**ST10** Surch-Lite 10" (254.0mm) sterile - 3/box



**ST15** Surch-Lite 15" (381.0mm) sterile - 3/box

### LIGHTED STYLETS



**NLOT** orotracheal stylet, non-sterile - 1/box



**SLOT** orotracheal stylet, sterile - 1/box

## UV WOODS LIGHT, SPECIALTY LIGHTS, PENLIGHTS

### SPECIALTY LIGHTS



**6802**  
Fluoro-Dot™ - 1/box

Specialty lights are used in conjunction with fluorescein dye to find foreign matter in the eye.

### WOODS LIGHT



**UV59**  
Hand-held woods light, excellent for detecting certain skin disorders and rape detection procedures - 1/box

### PENLIGHT PACKS

**6303**

cobalt penlight - 3/pack



**6666**

doctor's disposable penlight - 6/pack



**66RN**

penlight w/pupil gauge - 6/pack



**1125**

doctor's disposable penlight - 1/box










**7156**

cobalt filter only - 1/pack





## OPHTHALMIC PRODUCTS

Bovie® Medical manufactures a diverse line of ophthalmic products, including an ophthalmic burr, and a clear protective eye bubble for ophthalmic and neurological procedures.

OPHTHALMIC BURR HANDLE		EYE BUBBLE	
			
<p><b>0010</b> Bovie Burr 2 (includes one disposable 1 mm burr tip [0001] for corneal rust ring removal, "AA" alkaline battery and a foam-lined case, a diamond burr tip [0011] may be purchased separately for use in polishing the pterygium bed after surgical removal and for lid margin lesions) - 1/box</p>		<p><b>0002</b> eye bubble - 10/box</p>	
OPHTHALMIC BURR TIPS			
			
<p><b>0011</b> diamond burr 2.5mm for use in conjunction with the 0010 or equivalent - 1/box</p>		<p><b>0012</b> diamond burr 5mm for use in conjunction with the 0010 or equivalent - 1/box</p>	
			
<p><b>0001</b> 1 mm burr tip for use with the 0010 or equivalent for corneal rust ring removal, sterile - 10/box</p>	<p><b>AB01</b> 1mm burr tip for use with AB03, AB04 or equivalent, sterile - 10/box</p>	<p><b>AB05</b> 0.5mm burr tip for use with AB03, AB04 or equivalent, sterile - 10/box</p>	

## NERVE LOCATORS

Neuro-Pulse™ nerve locators were designed for identifying exposed motor nerves while at the same time reducing the possibility of accidental nerve damage or severance. These self-powered units are portable and battery operated. Ideal for head, neck, hand, plastic, and facial applications.

NEURO-PULSE


<p><b>0003Y (yellow LED, for use in Europe)</b> Neuro-Pulse has three output settings - 0.5 mA, 1 mA, and 2 mA at 3 volts DC - 10/box</p>



# NOTES

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

# NOTES

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---



UK Distributor

T 01702 602050

E [info@surgicalholdings.co.uk](mailto:info@surgicalholdings.co.uk)

# Bovie®

## **Bovie Medical**

Clearwater, FL 33760, USA

Int'l. Phone +1 727 384 2323 • Fax +1 727 347 9144

[Sales@BovieMed.com](mailto:Sales@BovieMed.com) • [BovieMedical.com](http://BovieMedical.com)

©2017 Bovie Medical Corporation, All Rights Reserved.

Bovie® is a registered trademark of Bovie Medical Corporation.

MC-55-275-001\_r0, 10/17