


CHIROBLOC

3 models with different ergonomies

3 lead hands without lead.

Made out of a composite bio-compatible material.



CHIROBLOC3 - ref. CHIRO3


Weight: 0.27 kg
Warranty: 1 year

Supplied with 5 ForgetMeNot
see accessories




CHIROBLOC Large - ref. CBM

Weight: 0.5 kg
Warranty: 1 year



CHIROBLOC Small - ref. CBP

Weight: 0.3 kg
Warranty: 1 an



DISTRIBUTED BY



Manufacturer:
AREX®
3 allée du Clos Tonnerre
91120 Palaiseau - FRANCE


www.arex.fr

Design: e-monneaux.fr / 254 éd.08 - © AREX® - 2019

www.arex.fr

CHIROBLOC3 - ref. CHIRO3

Made out of a composite bio-compatible material.
Radio-lucency without lead weight.
Will not stain surgical gloves.



One wrist strap: two possible adjustments



5 ForgetMeNot (FMNY2B/ FMNB2B) 2 functions:

- Reusable finger tourniquet
- ForgetMeNot units may be used as an attachment to any CHIROBLOC

ACCESSORIES



CHIROBLOC Ballast - ref. CBL

It is an independent spare part made out of stainless steel. It may be attached to CHIROBLOC to weigh it.
Weight: 0.725 kg
Warranty: 5 years



ForgetMeNot

Yellow color - ref. FMNY2B
Blue color - ref. FMNB2B
Sold by box of 10 units.

FMN is a reusable finger tourniquet, it is sold non sterile.
FMN may also be used as an attachment to any CHIROBLOC.

Warranty: 1 year



ForgetMeNot is specially useful to replace elastic bands made with surgical gloves.