


STILLE

Micro Instruments

Solid Performance - Outstanding Precision



STILLE offers the highest quality handmade instruments for the world's best surgeons

The line of STILLE Micro Instruments includes a wide range of stainless steel and titanium patterns. Our instruments provide the surgeon with the fingertip feel, perfect balance and reliability they expect during long and difficult cases. STILLE enhances your investment in quality by providing you with a 30-year instrument warranty.



SYMBOLS

SuperCut



Diamond Jaws



TC-inserts

STILLE MICRO SCISSORS

The STILLE Micro scissors consist of both straight and curved models in different lengths, blade sizes and angles. All models are made with the original STILLE Supercut blades for consistent cutting performance in soft tissue.

Our joint screw design allows for easy removal during instrument service and repair. After service, the instrument will perform with the same feel and reliability you come to expect from a STILLE instrument.

STILLE MICRO FORCEPS


All STILLE Micro Forceps are designed to give the same feel and balance as our standard line of forceps. Our attention to detail by our skilled artisans in Sweden, combines both perfect spring tension and tip alignment for reliable grasping of tissue.

They are available in Titanium or stainless steel that include tying platform, 1x2 teeth, ring tips and DeBakey. The tips are available in stainless steel or STILLE Diamond jaws for added durability.

STILLE MICRO NEEDLE HOLDERS

The STILLE Micro Needle Holders are designed to hold the finest sutures in a tight and secure grip. Only minimal pressure will need to be applied by surgeon's fingers.


All models are made with STILLE Diamond Jaws and the heavier models with TC-inserts. They can be delivered as straight or curved, with or without ratchet and in several different lengths. Regular box lock for 7-0 sutures and finer, Streamline box lock for 5-0 sutures and finer. (heavy models for 2-0 to 4-0).



Anatomically adapted
to the surgeon's hand


Satin finish prevents
glare, keeping the
surgeon's eyes
comfortable during
long operations





The edges, spring unit, pivot and grip are designed to prevent glove damage

Grooved handle in order to reduce the weight of the instrument



Micro Instruments

IN STAINLESS STEEL AND TITANIUM

- 8 **Micro Scissors** in Stainless Steel
- 10 **Micro Forceps** in Stainless Steel
- 13 **Micro Forceps** in Titanium
- 15 **Micro Needle Holders** in Stainless Steel
- 18 **Micro Needle Holders** in Titanium







STILLE TITANIUM MICRO INSTRUMENTS

Titanium is a non-magnetic metal that offers the surgeon a lighter, stronger instrument alternative. It minimizes the frustration over magnetized needles and is very resistant to corrosion. That provides you with an excellent long-term investment in quality.

STILLE Titanium Instruments come from only the finest, U.S. origin alloy (Ti6Al4V). The end result is improved durability and better fracture toughness than other titanium options.

STILLE Titanium Forceps and Needle Holders are clearly distinguishable by our blue coated round handles. The forceps (except for DeBakey models) and needle holders are all equipped with STILLE Diamond Jaws coating for improved durability and proper gripping of fine sutures.






Art. Number: 117-A110-18


STILLE SUPERCUT MICRO SCISSORS

Round handle, STILLE original 14 mm SuperCut blades.

ART. NUMBER	TIPS	BLADES	LENGTH	
117-A110-15	Sharp	Straight	15 cm 6"	•
117-A110-18	Sharp	Straight	18 cm 7"	•
117-A110-21	Sharp	Straight	21 cm 8.25"	•
117-A111-15	Blunt	Curved	15 cm 6"	•
117-A111-18	Blunt	Curved	18 cm 7"	•
117-A111-21	Blunt	Curved	21 cm 8.25"	•
117-A112-15	Blunt	Straight	15 cm 6"	•
117-A112-18	Blunt	Straight	18 cm 7"	•
117-A112-21	Blunt	Straight	21 cm 8.25"	•
117-A113-15	Sharp	Curved	15 cm 6"	•
117-A113-18	Sharp	Curved	18 cm 7"	•
117-A113-21	Sharp	Curved	21 cm 8.25"	•

STILLE SUPERCUT MICRO VASCULAR SCISSORS

Round handle, STILLE original fine 11 mm SuperCut blades.

ART. NUMBER	ANGLED	LENGTH	
117-A120-25	25°	18 cm 7"	•
117-A120-45	45°	18 cm 7"	•
117-A120-60	60°	18 cm 7"	•
117-A120-90	90°	18 cm 7"	•
117-A120-125	125°	18 cm 7"	•
117-A121-25	25°	21 cm 8.25"	•
117-A121-45	45°	21 cm 8.25"	•
117-A121-60	60°	21 cm 8.25"	•
117-A121-90	90°	21 cm 8.25"	•
117-A121-125	125°	21 cm 8.25"	•




Art. Number: 117-A120-45



Art. Number: 117-A125-60

STILLE SUPERCUT FINE MICRO VASCULAR SCISSORS

Round handle, STILLE original extra fine 7.5 mm SuperCut blades.

ART. NUMBER	ANGLED	LENGTH	
117-A125-25	25°	18 cm 7"	•
117-A125-45	45°	18 cm 7"	•
117-A125-60	60°	18 cm 7"	•
117-A125-90	90°	18 cm 7"	•
117-A125-125	125°	18 cm 7"	•

“STILLE stands for impeccable quality, they’re dependable and the durability is just outstanding.”

- Dr. Mark L. Jewell, M.D.
Plastic Surgeon

STILLE MICRO FORCEPS

Well balanced and ergonomic micro forceps of ultra high quality. With tying platform and 0.4 mm tips.


ART. NUMBER	JAWS	LENGTH
117-B110-15	Straight	15 cm 6"
117-B110-18	Straight	18 cm 7"
117-B110-21	Straight	21 cm 8.25"
117-B110-23	Straight	23 cm 9"
117-B111-15	Curved	15 cm 6"
117-B111-18	Curved	18 cm 7"
117-B111-21	Curved	21 cm 8.25"



Art. Number: 117-B110-18

STILLE DIAMOND JAWS MICRO FORCEPS

Fine micro forceps, with Diamond Jaws, for increased tip durability and secure grip of fine needles. 0.4 mm tips.


ART. NUMBER	JAWS	LENGTH	
117-B120-15	Straight	15 cm 6"	•
117-B120-18	Straight	18 cm 7"	•
117-B120-21	Straight	21 cm 8.25"	•
117-B120-23	Straight	23 cm 9"	•
117-B121-15	Curved	15 cm 6"	•
117-B121-18	Curved	18 cm 7"	•
117-B121-21	Curved	21 cm 8.25"	•



Art. Number: 117-B121-15

STILLE DIAMOND JAWS COUNTER-BALANCED MICRO FORCEPS

Micro forceps with Diamond Jaws and Counter-balance.
0.4 mm tips.

ART. NUMBER	JAWS	LENGTH	
117-B130-18	Straight	18 cm 7"	•
117-B130-21	Straight	21 cm 8.25"	•
117-B130-23	Straight	23 cm 9"	•
117-B131-18	Curved	18 cm 7"	•
117-B131-21	Curved	21 cm 8.25"	•



Art. Number: 117-B130-18

STILLE MICRO RING FORCEPS

Counter-balanced micro forceps with fine 1.0x0.5 mm Diamond Jaws ring tips.

ART. NUMBER	RING TIPS	JAWS	LENGTH	
117-B140-18	1.0x0.5 mm	Straight	18 cm 7"	•
117-B140-21	1.0x0.5 mm	Straight	21 cm 8.25"	•
117-B140-23	1.0x0.5 mm	Straight	23 cm 9"	•



Art. Number: 117-B140-18

STILLE MICRO TISSUE FORCEPS

Round handle, with 1x2 fine teeth on 0.4 mm tips.

ART. NUMBER	TEETH	JAWS	LENGTH
117-B150-15	1x2	Straight	15 cm 6"
117-B150-18	1x2	Straight	18 cm 7"
117-B150-21	1x2	Straight	21 cm 8.25"



Art. Number: 117-B150-15

STILLE MICRO DEBAKEY FORCEPS

Micro forceps with round handle and 1.2 mm fine atraumatic DeBakey tips.

ART. NUMBER	TIPS	JAWS	LENGTH
117-B160-18	1.2 mm Debakey	Straight	18 cm 7"
117-B160-21	1.2 mm Debakey	Straight	21 cm 8.25"



Art. Number: 117-B160-18

STILLE MICRO DEBAKEY COUNTER-BALANCED FORCEPS

Counter-balanced micro forceps with round handle and 1.2 mm fine atraumatic DeBakey tips.


ART. NUMBER	TIPS	JAWS	LENGTH
117-B165-18	1.2 mm Debakey	Straight	18 cm 7"
117-B165-21	1.2 mm Debakey	Straight	21 cm 8.25"



Art. Number: 117-B165-18

STILLE TITANIUM MICRO FORCEPS

Round handle, with tying platform, straight or curved 0.4 mm tips.


ART. NUMBER	JAWS	LENGTH	
117-B120-18-TI	Straight	18 cm 7"	•
117-B120-21-TI	Straight	21 cm 8.25"	•
117-B120-23-TI	Straight	23 cm 9"	•
117-B120-25-TI	Straight	25 cm 10"	•
117-B121-18-TI	Curved	18 cm 7"	•
117-B121-21-TI	Curved	21 cm 8.25"	•
117-B121-23-TI	Curved	23 cm 9"	•
117-B121-25-TI	Curved	25 cm 10"	•



Art. Number: 117-B120-18-TI

STILLE TITANIUM COUNTER-BALANCED MICRO FORCEPS

Round handle, with tying platform, counter-balanced, straight or curved 0.4 mm tips.


ART. NUMBER	JAWS	LENGTH	
117-B130-18-TI	Straight	18 cm 7"	•
117-B130-21-TI	Straight	21 cm 8.25"	•
117-B131-18-TI	Curved	18 cm 7"	•
117-B131-21-TI	Curved	21 cm 8.25"	•



Art. Number: 117-B130-21-TI

STILLE TITANIUM MICRO RING FORCEPS

Round handle, with fine ring tips, counter-balanced.

ART. NUMBER	RING TIPS	JAWS	LENGTH	
117-B140-18-TI	1.0x0.5 mm	Straight	18 cm 7"	•
117-B140-21-TI	1.0x0.5 mm	Straight	21 cm 8.25"	•
117-B140-23-TI	1.0x0.5 mm	Straight	23 cm 9"	•



Art. Number: 117-B140-21-TI

STILLE TITANIUM MICRO DEBAKEY FORCEPS

Round handle, with 1.2 mm fine atraumatic DeBakey tips.

ART. NUMBER	TIPS	JAWS	LENGTH
117-B160-18-TI	1.2 mm DeBakey	Straight	18 cm 7"
117-B160-21-TI	1.2 mm DeBakey	Straight	21 cm 8.25"




Art. Number: 117-B160-18-TI



Art. Number: 117-C110-15


STILLE MICRO NEEDLE HOLDER REGULAR BOX LOCK WITHOUT RATCHET

Round handle, fine 0.8 mm Diamond Jaws and Regular Box Lock for sutures 7-0 and smaller. Without ratchet.

ART. NUMBER	JAWS	LENGTH	
117-C110-15	Straight	15 cm 6"	•
117-C110-18	Straight	18 cm 7"	•
117-C110-21	Straight	21 cm 8.25"	•
117-C110-23	Straight	23 cm 9"	•
117-C111-15	Curved	15 cm 6"	•
117-C111-18	Curved	18 cm 7"	•
117-C111-21	Curved	21 cm 8.25"	•
117-C111-23	Curved	23 cm 9"	•

STILLE MICRO NEEDLE HOLDER REGULAR BOX LOCK WITH RATCHET

Round handle, fine 0.8 mm Diamond Jaws and Regular Box Lock for sutures 7-0 and smaller. With ratchet.

ART. NUMBER	JAWS	LENGTH	
117-C120-15	Straight	15 cm 6"	•
117-C120-18	Straight	18 cm 7"	•
117-C120-21	Straight	21 cm 8.25"	•
117-C120-23	Straight	23 cm 9"	•
117-C121-15	Curved	15 cm 6"	•
117-C121-18	Curved	18 cm 7"	•
117-C121-21	Curved	21 cm 8.25"	•
117-C121-23	Curved	23 cm 9"	•




Art. Number: 117-C121-18



Art. Number: 117-C115-15


STILLE MICRO NEEDLE HOLDER STREAMLINE BOX LOCK WITHOUT RATCHET

Round handle, 1.2 mm Diamond Jaws and Streamline Box Lock for sutures 5-0 and smaller. Without ratchet.

ART. NUMBER	JAWS	LENGTH	
117-C115-15	Straight	15 cm 6"	•
117-C115-18	Straight	18 cm 7"	•
117-C115-21	Straight	21 cm 8.25"	•
117-C115-23	Straight	23 cm 9"	•
117-C116-15	Curved	15 cm 6"	•
117-C116-18	Curved	18 cm 7"	•
117-C116-21	Curved	21 cm 8.25"	•
117-C116-23	Curved	23 cm 9"	•

STILLE MICRO NEEDLE HOLDER STREAMLINE BOX LOCK WITH RATCHET

Round handle, 1.2 mm Diamond Jaws and Streamline Box Lock for sutures 5-0 and smaller. With ratchet.

ART. NUMBER	JAWS	LENGTH	
117-C125-15	Straight	15 cm 6"	•
117-C125-18	Straight	18 cm 7"	•
117-C125-21	Straight	21 cm 8.25"	•
117-C125-23	Straight	23 cm 9"	•
117-C126-15	Curved	15 cm 6"	•
117-C126-18	Curved	18 cm 7"	•
117-C126-21	Curved	21 cm 8.25"	•
117-C126-23	Curved	23 cm 9"	•



Art. Number: 117-C125-15



Art. Number: 117-C140-25

STILLE MICRO NEEDLE HOLDER HEAVY MODEL

Extra heavy model for handling of larger needles with TC-inserts. Intended for 2/0 suture and finer.


ART. NUMBER	RATCHET	JAWS	LENGTH	TC
117-C130-19	Without	Straight	19 cm 7.5"	•
117-C130-23	Without	Straight	23 cm 9"	•
117-C130-25	Without	Straight	25 cm 10"	•
117-C140-19	With	Straight	19 cm 7.5"	•
117-C140-23	With	Straight	23 cm 9"	•
117-C140-25	With	Straight	25 cm 10"	•

“An investment in a STILLE instrument is an investment that will last throughout your career.”

- Dr. Jordan Deschamps-Braly, M.D.
Plastic Surgeon

**STILLE TITANIUM MICRO NEEDLE HOLDER
REGULAR BOX LOCK**

Round handle, fine 0.8 mm Diamond Jaws and Regular Box Lock for sutures 7-0 and smaller. With or without ratchet.


ART. NUMBER	RATCHET	JAWS	LENGTH	
117-C110-18-TI	Without	Straight	18 cm 7"	•
117-C110-21-TI	Without	Straight	21 cm 8.25"	•
117-C110-23-TI	Without	Straight	23 cm 9"	•
117-C111-15-TI	Without	Curved	15 cm 6"	•
117-C111-18-TI	Without	Curved	18 cm 7"	•
117-C111-21-TI	Without	Curved	21 cm 8.25"	•
117-C120-18-TI	With	Straight	18 cm 7"	•
117-C120-21-TI	With	Straight	21 cm 8.25"	•
117-C120-23-TI	With	Straight	23 cm 9"	•
117-C121-18-TI	With	Curved	18 cm 7"	•
117-C121-21-TI	With	Curved	21 cm 8.25"	•
117-C121-23-TI	With	Curved	23 cm 9"	•



Art. Number: 117-C110-18-TI



**STILLE TITANIUM MICRO NEEDLE HOLDER
STREAMLINE BOX LOCK**

Round handle, 1.2 mm Diamond Jaws and Streamline Box Lock for sutures 5-0 and smaller. With or without ratchet.

ART. NUMBER	RATCHET	JAWS	LENGTH	
117-C115-18-TI	Without	Straight	18 cm 7"	•
117-C115-21-TI	Without	Straight	21 cm 8.25"	•
117-C115-23-TI	Without	Straight	23 cm 9"	•
117-C125-18-TI	With	Straight	18 cm 7"	•
117-C125-21-TI	With	Straight	21 cm 8.25"	•
117-C125-23-TI	With	Straight	23 cm 9"	•
117-C126-18-TI	With	Curved	18 cm 7"	•
117-C126-21-TI	With	Curved	21 cm 8.25"	•
117-C126-23-TI	With	Curved	23 cm 9"	•



Art. Number: 117-C125-21-TI



“My patients deserve the best and STILLE is the natural instrument choice which enables me to perform at my best.”

- Dr. Mark L. Jewell, M.D.
Plastic Surgeon



**Do you want to
learn more?**

Scan the QR code to
visit our website



950-G117-0001-EN-01

© 2023-10-31 STILLE. All rights reserved.

Surgical perfection. For life.



Surgical perfection. For life.

STILLE AB

Ekbacken 11
SE-644 30 Torshälla, Sweden
☎ +46 8 588 580 00
✉ info@stille.se
🌐 www.stille.se

STILLE SURGICAL INC.

840 Parkview Blvd
Lombard, IL 60148
United States
Phone: 1.800.6651614

DISTRIBUTED BY:

Surgical Holdings
Parkside Centre
Temple Farm Ind. Est.
Southend-on-Sea
Essex SS2 5SJ
T: 01702 602050
E: info@surgicalholdings.co.uk