



# Hip Instruments Catalogue

Version 2



 Proud to be a family-run British company

  
Promoting education;  
sharing knowledge and  
offering training opportunities

## Contents

About Surgical Holdings	2
Curettes	3
Elevators	4
Extractor	6
Hip Set	6
Hooks	7
Impactor	7
Levers	7
Mallets	8
Nail	9
Osteotomes, Chisels and Gouges	10
Pressurisers	12
Pushers	13
Rasp	14
Restrictor	14
Retractors	15
Skids	18

**We're proud to present a selection of our reusable instrumentation which may be suitable for the compilation of your hip surgery sets. The sets can be tailored to your requirements.**



**We have over 50 years of knowledge and experience in our workshop. Our team of highly skilled technicians undertake a minimum apprenticeship of 5 years and are trained to manufacture instruments, rather than just repair them. This gives us 360 degree understanding of the challenge and the importance of correctly functioning medical devices.**

We care immensely about the work we do, never forgetting that instruments we create and repair are responsible for saving lives.

We insist on controlling the quality of every item we manufacture, from the materials used to the rigorous procedures we apply to testing. Our range of supporting products are serviced and repaired in-house, ensuring we are always able to exceed our customers' expectations.

Customers can receive training, participate in factory visits or simply call their local representative for guidance with instrumentation or management.

Skill-sharing has been at the heart of our business for four generations and we frequently arrange workshops and attend events at which we can educate others.



## How To Order

- 1 Choose your desired product from our extensive range or send your detailed tray list to us and we'll do the work.
- 2 Email [info@surgicalholdings.co.uk](mailto:info@surgicalholdings.co.uk) or call +44 (0)1702 602050 with the product codes and we will provide you with a quote.
- 3 Upon approval of the quote we will process your order.



**NEED A BESPOKE INSTRUMENT?**  
CONTACT OUR FRIENDLY TEAM  
TO DISCUSS YOUR REQUIREMENTS

## How To Pay

30 days credit is available (subject to status checks) or we accept Visa, Mastercard, Amex cards and cheque payments.


**30 DAYS CREDIT**



# Curettes

## Charnley Bone Curette


With T handle, 380mm long



Code	Description
SC-OR/N040	Large, 20 x 14mm
SC-OR/N041	Medium, 17 x 12mm
SC-OR/N042	Small, 11 x 8mm

## Charnley Ring Curette

With T handle



Code
SC-OR/N042C



## Refurbishment of orthopaedic sets

Refurbishing  
instruments  
can save up to  
**80%**

of the cost of buying new,  
with compelling evidence to  
demonstrate effectiveness

Independent research<sup>1</sup> showed that instruments  
refurbished by Surgical Holdings had a 4x reduction in  
the amount of residual protein.

The smooth and polished surface has an increased cleaning  
efficacy, and all refurbishments come with an 18 month  
warranty as standard. Standard turnaround times are  
1 week and an express 2 day service is also available.

Reference: <sup>1</sup>Data available on file.

Book your orthopaedic set refurbishment  
on [info@surgicalholdings.co.uk](mailto:info@surgicalholdings.co.uk)




Before  
Refurb

After  
Refurb

## Elevators

### Hohmann-Aldinger Bone Elevator


Curved



Code	Description
SC-OR/304	270mm long x 24mm wide

### Hohmann-Muller Spike Elevator

With narrow spade



Code	Description
SC-OR/307f	220mm long, 10mm wide
SC-OR/307	230mm long, 12mm wide
SC-OR/308	240mm long, 18mm wide

### Hohmann-Muller Spike Bone Elevator


Curved, 235mm long



Code	Description
SC-OR/308D	45mm wide

### Hohmann Bone Elevator

Curved



Code	Description
SC-OR/305	285mm long x 34mm wide


## Sterile Arthrodesis & Kirschner Wires

- ✓ Highest grade implant quality stainless steel to ISO 5832-1
- ✓ Reduces risk of sharps injuries to Sterile Services personnel
- ✓ All popular sizes in stock and can manufacture any custom size, including threaded wires



### Hohmann Bone Lever

Curved, Leicester Royal Infirmary pattern



Code	Description
SC-OR/LGH-HOH	340mm long x 23mm wide

## Elevators

### Lanes Elevator

355mm long



#### Code

SC-OR/N049

### Lange-Hohmann Bone Elevator

250mm long




#### Code

SC-OR/307i

#### Description

24mm wide

### Muller Soft Tissue Retractor



#### Code

SC-OR/307aNOR

#### Description


410mm long, 40mm wide

SC-OR/307bNOR

300mm long, 40mm wide

### Schumacher Bone Elevator

290mm long



#### Code


SC-OR/307k

#### Description

35mm wide, 18mm at edge

### Schumacher Bone Elevator

345mm long



#### Code

SC-OR/307j

#### Description

37mm wide, 18mm at edge

### Wagner Bone Elevator

220mm long



#### Code

SC-OR/307h

#### Description


17mm wide

## Extractor

### Moore Femoral Head Extractor

With T handle

Code	Description
SC-OR/N005	150mm long
SC-OR/N006	205mm long



## Hip Set

### Curved Osteotome

260mm



Code	Description
SC-OR/N059D	With blue acetal suregrip handle - 2

### Muller Curette

270mm



Code	Description
SC-OR/N059C	Curette 9mm with blue acetal suregrip handle - small

### Femoral Dislocation Lever

285mm



Code	Description
SC-OR/N059A	Spoon elevator 25mm with blue acetal suregrip handle - large

### Small Gouge

240mm



Code	Description
SC-OR/N059E	8mm cranked blue acetal suregrip handle

### Muller Curette

240mm




Code	Description
SC-OR/N059B	Curette 27mm with blue acetal suregrip handle - large

## Hooks


### Dislocation Hip Hook

Blue acetal "Suregrip" handle



Code	Description
SC-OR/420	185mm long

### Dislocation Hook




Code	Description
SC-OR/SH-5703	253mm long

## Impactor

### Judet Impactor

Stainless steel handle with nylon head




Code	Description
SC-OR/N009	140mm long

## Levers

### Bennett Cobra Bone Lever

240mm long



Code
SC-OR/425



### Charnley Femoral Lever

300mm long



Code
SC-OR-694

### Hohmann Bone Lever

320mm long



Code
SC-OR/646

### Tager Bone Lever Angle Joint

175mm long




Code	Description
SH009129	12mm wide



## Mallets

### Heath Mallet


Stainless steel, 135mm long



Code	Description
SC-OR/N044c	650 grams
SC-OR/N044b	920 grams
SC-OR/N044	1150 grams
SC-OR/N044a	1340 grams

### Mallet


Blue acetal "Suregrip" handle, 275mm long



Code
SC-OR/270A-BH

### Mallet


Blue acetal, 255mm long



Code	Description
SC-OR/22a	520 grams, vertical end
SC-OR/22	650 grams, horizontal end

### Mallet


Stainless steel, 270mm long



Code	Description
SC-OR/270	640 grams
SC-OR/271	1020 grams

### Orthopaedic Mallet


Surgical Holdings pattern, stainless steel, 170mm long



Code	Description
SC-OR/N055b	800 grams
SC-OR/N055	880 grams
SC-OR/N055a	1020 grams

### Surgical Holdings Mallet

Stainless steel with blue acetal "Suregrip" handle domed faced head, 200mm long




Code	Description
SC-OR/418BH	650 grams



## Mallets

### Surgical Holdings Orthopaedic Mallet


Stainless steel with blue acetal "Suregrip" handle  
200mm long



Code	Description
SC-OR/414BH	660 grams
SC-OR/415BH	790 grams
SC-OR/416BH	1020 grams
SC-OR/417BH	1400 grams

### Surgical Holdings Orthopaedic Mallet

Stainless steel with blue acetal "Suregrip" handle  
270mm long



Code	Description
SC-OR/410	720 grams
SC-OR/411	800 grams
SC-OR/412	900 grams
SC-OR/413	970 grams

## Nail

### Square Headed Nail

150mm long



Code	Description
SC-OR/N043B	6mm wide

Did you know...

**We repair power tools,  
pneumatic tools and hoses**




- ✓ Stryker, Synthes, Hall, Microaire, De Soutter, and all other common brands repaired
- ✓ All dental tools can be repaired, including bearings and rotor replacement
- ✓ We also repair Phaco Handpieces and Ultrasound Handpieces

[info@surgicalholdings.co.uk](mailto:info@surgicalholdings.co.uk)

## Osteotomes, Chisels and Gouges

### Capener Osteotome

Blue acetal "Suregrip" handle, straight, 255mm long



Code	Description	Code	Description
SC-OR/374	Width 6mm	SC-OR/379	Width 15mm
SC-OR/375	Width 8mm	SC-OR/380	Width 18mm
SC-OR/376	Width 10mm	SC-OR/381	Width 20mm
SC-OR/377	Width 12mm	SC-OR/382	Width 25mm
SC-OR/378	Width 13mm	SC-OR/383	Width 30mm

### Capener Osteotome

Blue acetal "Suregrip" handle, 250mm long



#### Shallow curved

Code	Description
SC-OR/384	Width 6mm
SC-OR/385	Width 10mm
SC-OR/386	Width 12mm
SC-OR/387	Width 15mm
SC-OR/388	Width 20mm
SC-OR/389	Width 25mm
SC-OR/390	Width 30mm

#### Acute curved

Code	Description
SC-OR/390a	Width 6mm
SC-OR/390b	Width 10mm
SC-OR/390c	Width 15mm
SC-OR/390d	Width 20mm
SC-OR/390e	Width 25mm

### Lexer Osteotome

Blue acetal "Suregrip" handle, 180mm long



Code	Description
SC-OR/SH5347a	4mm

### Moore Hollow Osteotome

Rectangular, for curved standard stems



Code	Description
SC-OR/N008	135mm long
SC-OR/N008a	255mm long

### Box Chisel

Blue acetal "Suregrip" handle, 285mm long



Code
SC-OR/N008b

### Capener Chisel

Blue acetal "Suregrip" handle straight, 250mm long




Code	Description	Code	Description
SC-OR/368	Width 6mm	SC-OR/371	Width 15mm
SC-OR/369	Width 10mm	SC-OR/372	Width 20mm
SC-OR/370	Width 12mm	SC-OR/373	Width 25mm

# Osteotomes, Chisels and Gouges

## Capener Lamina Gouge

Cranked, blue acetal "Suregrip" handle, 250mm long




### With reversed edge

Code	Description	Code	Description
SC-OR/373c	Width 6mm	SC-OR/373cre	Width 6mm
SC-OR/373d	Width 8mm	SC-OR/373dre	Width 8mm
SC-OR/373j	Width 10mm	SC-OR/373jre	Width 10mm
SC-OR/373e	Width 11mm	SC-OR/373ere	Width 11mm
SC-OR/373f	Width 16mm	SC-OR/373fre	Width 16mm
SC-OR/373k	Width 17mm	SC-OR/373kre	Width 17mm

OTHER LENGTHS AND SIZES AVAILABLE

## Capener Lamina Gouge

Swan neck, blue acetal "Suregrip" handle, 250mm long



Code	Description
SC-OR/373g	Width 8mm

## Capener Lamina Gouge


Blue acetal "Suregrip" handle, straight, 265mm long



Code	Description	Code	Description
SC-OR/362	Width 6mm	SC-OR/365	Width 15mm
SC-OR/363	Width 10mm	SC-OR/366	Width 20mm
SC-OR/364	Width 12mm	SC-OR/367	Width 25mm

## Capener Lamina Gouge

Reverse cranked, blue acetal "Suregrip" handle, 245mm long




Code	Description	Code	Description
SC-OR/373ar	Width 11mm	SC-OR/373or	Width 6mm
SC-OR/373nr	Width 16mm	SC-OR/373lr	Width 8mm
SC-OR/373jr	Width 17mm	SC-OR/373mr	Width 10mm

## Pressurisers

### Acetabular Pressuriser


260mm long, with T handle



Code	Description
SH-5636A	25mm x 6mm end

### Acetabular Pressuriser


370mm long, domed head, with T handle



Code
SH-5636-A

### Acetabular Pressuriser

370mm, long flat end 40mm x 14mm, with T handle



Code
SH-5636

# Pushers

## Charnley Cup Pusher


With nylon head



Code	Description
SC-OR/N020	295mm long

## Charnley Prosthesis Pusher


Stainless steel handle with nylon insert



Code	Description
SC-OR/N021	180mm long

## Femoral Pusher

Head impactor, white Pomalux




Code	Description
SC-OR/N053	31mm ø x 180mm long

## Rasp

### Bone File and Rasp

216mm long

Code	Description
SH002424	18mm wide



## Restrictor

### Hardinge Cement Restrictor

With T handle

Code	Description
SC-OR/HR001	335mm long



## Arthroscope Repair

Surgical Holdings' highly skilled and experienced technicians will return scopes to the manufacturer's original specification. All parts are brand new and have full traceability. Surgical Holdings also manufacture their own fibre optic light bundles.



18 month guarantee on the repair



ISO 13485 accredited as recommended by MHRA\* advice around endoscope storage and care



OEM optical specification followed where possible



Join a free educational visit to our scope workshop and see all of our repair facilities



Personal collection and delivery



Laser welding for sheaths

\*MHRA Managing Medical Devices – April 2016



## Retractors

### Anterior Retractor

410mm long

Code
SC-OR/5914




### Anterior Hip Retractor

400mm long

Code	Description
SC-OR/5913	22mm wide



### Charnley Initial Incision Retractor



Code	Description
SC-OR/N038a	Frame only
SC-OR/N038b	Weight only, 75mm x 50mm ø
SC-OR/N038s	Goal post pins and chain

#### Short blades

Code	Description
SC-OR/N038r	Short blade, 25mm deep
SH006084	Short blade, 35mm deep
SC-OR/N038c	Short blade, 50mm deep
SC-OR/N038l	Short blade, 65mm
SC-OR/N038e	Short blade, 75mm deep
SC-OR/N038q	Short blade, 80mm deep
SC-OR/N038g	Short blade, 100mm deep
SH003879	Short blade, 125mm
SH003880	Short blade, 150mm


#### Long blades

Code	Description
SH006085	Long blade, 35mm deep
SC-OR/N038d	Long blade, 50mm deep
SC-OR/N038k	Long blade, 65mm
SC-OR/N038p	Long blade, 70mm deep
SC-OR/N038f	Long blade, 75mm deep
SC-OR/N038h	Long blade, 100mm deep
SH003877	Long blade, 125mm
SH003878	Long blade, 150mm
SC-OR/N038m	Long round blade, 50mm deep

## Retractors

### Charnley Pin Retractor

With goal posts



**Code**

SC-OR/N038I

### Cobra Retractor

350mm long, double prong 51mm



**Code**

SH007226

### Exeter Inferior Acetabular Retractor



**Code**

SC-OR/5913A

**Description**


400mm long, 22mm wide

SC-OR/5913B

345mm long, 42mm wide

### Femoral Neck Delivery Retractor

330mm long




**Code**

SC-OR/642

### Gluteus Medius Retractor

175mm long



**Code**

SC-OR/695

### Hohmann Angled Knee Retractor

225mm long



**Code**


SH/5252



## Retractors

### Hohmann Retractor


360mm long



Code	Description
SH003859	34mm wide blade

### Large Hip Retractor

330mm long



Code
SC-OR/307L

### Lateral Patellar Retractor

265mm long double prong, 20mm wide



Code
SC-OR/710

### Muller Femoral Retractor

350mm long, with blue acetal handle



Code
SC-OR/705

### Muller Soft Tissue Retractor

410mm long



Code	Description
SC-OR/696	30mm wide

### PCL Retractor

240mm long



Code
SC-OR/608

### PCL Tibial Retractor



Code
SC-OR/711

## Skids

### Murphy Skids

Double ended, 355mm long



Code	Description
SC-OR/N007	Uncoated
SC-DIA/111	Coated

### Murphy Skids

Single ended, 330mm long



Code	Description
SC-OR/N007a	Coated

## Index

### A

Acetabular Pressuriser	12
Anterior Hip Retractor	15
Anterior Retractor	15

### B

Bennett Cobra Bone Lever	7
Bone File and Rasp	14
Box Chisel	10

### C

Capener Chisel	10
Capener Lamina Gouge	11
Capener Osteotome	10
Charnley Bone Curette	3

Charnley Cup Pusher	13
Charnley Femoral Lever	7
Charnley Initial Incision Retractor	15
Charnley Pin Retractor	16
Charnley Prosthesis Pusher	13
Charnley Ring Curette	3
Cobra Retractor	16
Curved Osteotome	6

### D









Dislocation Hip Hook	7
Dislocation Hook	7

### E

Exeter Inferior Acetabular Retractor	16
--------------------------------------	----

<b>F</b>		
Femoral Dislocation Lever	6	
Femoral Neck Delivery Retractor	16	
Femoral Pusher	13	
<b>G</b>		
Gluteus Medius Retractor	16	
<b>H</b>		
Hardinge Cement Restrictor	14	
Heath Mallet	8	
Hohmann-Aldinger Bone Elevator	4	
Hohmann Angled Knee Retractor	16	
Hohmann Bone Elevator	4	
Hohmann Bone Lever	4, 7	
Hohmann-Muller Spike Bone Elevator	4	
Hohmann-Muller Spike Elevator	4	
Hohmann Retractor	17	
<b>J</b>		
Judet Impactor	7	
<b>L</b>		
Lanes Elevator	5	
Lange-Hohmann Bone Elevator	5	
Large Hip Retractor	17	
Lateral Patellar Retractor	17	
Lexer Osteotome	10	
<b>M</b>		
Mallet	8	
Moore Femoral Head Extractor	6	
Moore Hollow Osteotome	10	
Muller Curette	6	
Muller Femoral Retractor	17	
Muller Soft Tissue Retractor	5, 17	
Murphy Skids	18	
<b>O</b>		
Orthopaedic Mallet	8	
<b>P</b>		
PCL Retractor	17	
PCL Tibial Retractor	17	
<b>S</b>		
Schumacher Bone Elevator	5	
Small Gouge	6	
Square Headed Nail	9	
Surgical Holdings Mallet	8	
Surgical Holdings Orthopaedic Mallet	9	
<b>T</b>		
Tager Bone Lever Angle Joint	7	
<b>W</b>		
Wagner Bone Elevator	5	

## Instrument Sales

-  Surgical Holdings Instruments
-  Symmetry Surgical
-  Hermann Laparoscopic
-  Bissinger Bipolar
-  Light Leads
-  DIN Baskets & Containers
-  Sterile Packed K-Wires
-  Custom-made Instruments

## Refurbish & Repair

-  Surgical Instruments
-  Laparoscopic Instrumentation
-  Electrosurgery Devices
-  Rigid Endoscopes



**Surgical Holdings** Parkside Centre, Temple Farm Ind. Est., Southend-on-Sea, Essex, SS2 5SJ  
E. [info@surgicalholdings.co.uk](mailto:info@surgicalholdings.co.uk) T. +44 (0)1702 602050  
[surgicalholdings.co.uk](http://surgicalholdings.co.uk)

