




General Surgical Instruments Catalogue





Working in close partnership with our customers, we manufacture innovative and bespoke instrumentation.

We have over 50 years of knowledge and experience in our workshop. Our team of highly skilled technicians undertake a minimum apprenticeship of 5 years and are trained to manufacture instruments, rather than just repair them. This gives us 360 degree understanding of the challenge and the importance of correctly functioning medical devices.



We care immensely about the work we do, never forgetting that instruments we create and repair are responsible for saving lives.

We insist on controlling the quality of every item we manufacture, from the materials used to the rigorous procedures we apply to testing. Our range of supporting products are serviced and repaired in-house, ensuring we are always able to exceed our customers' expectations. Customers can receive training, participate in factory visits or simply call their local representative for guidance with instrumentation or management.

Skill-sharing has been at the heart of our business for four generations and we frequently arrange workshops and attend events at which we can educate others.

Quality Assurance

We are an ISO 13485 accredited company.

Our comprehensive quality management system is in accordance with the requirements of the Medical Device Directive 93/42/EEC - Annex I & II, for temporary orthopaedic surgical implants - skeletal pins, wires, staples, guide wires - and all associated instrumentation.

All Surgical Holdings products undergo rigorous testing and inspection to applicable harmonised standards before leaving the factory. We offer CE marking on all Class I instrumentation, Class IIa and Ilb implants, and can provide certificates of conformity for all materials used. Complete traceability throughout all aspects of the production process provides our customers with a fully documented quality control system.

All instrumentation carries a lifetime guarantee from the date of purchase on any fault found as a result of poor materials. This warranty is void, however, if repairs are carried out by anyone other than Surgical Holdings repair personnel, or if it is deemed to have been misused.

Please contact us or visit our website for current reprocessing instruction.





Visit us for a Workshop Tour



Want to see our strong family ethos in action? Our 24 hour promise means you can drop in on us at short notice for a guided tour. No grand clean up or special planning will be needed - we're confident for you to see us as we really are.

Pop in any day and you'll see:

- Instrument manufacturing
- Ortho implants manufacturing Instrument repair
- Laparoscopic instrument repair
- Power tool repair
 - Scope repair

We've invested in new machinery to help improve quality and reduce turnaround times. We'll show you our recent upgrades including:

- £100k worth of ultrasonic cleaners, taking the cleaning and passivating process from 30 minutes to just 3 minutes!
- Laser welders
- Marker machines



We're in an excellent location for export customers, just 40 miles from London, at the mouth of the Thames Estuary.

- By road: A13 or A127 to M25 Dartford Tunnel and south coast ports
- By rail: Direct links to London Liverpool Street, London Fenchurch Street & St Pancras International via Stratford
- **By air:** International flights from Heathrow or Gatwick 90 minutes away, Stansted 60 minutes, and Southend airport is just 10 minutes away with direct flights to major cities across Europe



Surgical Holdings, Parkside Centre, Temple Farm Industrial Estate. Southend-on-Sea, Essex, SS2 5SJ

Contents

Introduction	2
About Surgical Holdings	4
Our Services	6
CPD Training Procedures	10
CPD Accredited Webinars	11
Case Study	12
General Instruments	13
Artery Forceps	14
Clamps	17
Spring Forceps	19
Tissue Forceps	20
Sponge Holding and	22
Sterilising Forceps	
Other Forceps	23
Scissors	24
Needle Holders	29
Combined Scissor and	33
Needle Holders	
Bag and Towel Clips	34
Retractors	35
Dissectors and Probes	43
Miscellaneous	44
Hollow Ware	45
Silicone Matting	46
Sterilising Baskets	46
Optical Fibre Cables	48
Index	50
Our Manufacturing Partners	54

How To Order

- 1 Choose your desired product from our extensive range or send your detailed tray list to us and we'll do the work.
- Email info@surgicalholdings.co.uk or call +44 (0)1702 602050 with the product codes and we will provide , you with a quote.

Upon approval of the quote we will process your order.



NEED A BESPOKE **INSTRUMENT?** CONTACT OUR FRIENDLY TEAM TO DISCUSS YOUR REQUIREMENTS

How To Pay

30 days credit is available (subject to status checks) or we accept Visa, Mastercard, Amex cards and cheque payments.








Surgical Holdings is a family run business, established for 30 years and counting.

Based in Southend-on-Sea, we are an award-winning British surgical instrument manufacturer and repairer, facilitating sustainability in surgery; through our range of reusable instruments and laparoscopic instruments, of which can be serviced in-house by our experienced and specialist technicians.

Our SMART REPAIR service for rigid endoscopes and our 'as new' full surgical tray refurbishments, save 50% in CO² vs buying new. We also distribute key industry brands, including the Bookwalter retractor.



Our Company History



Mission, Vision and Values

Our **Mission** is to sustainably manufacture, service and distribute medical devices for the UK healthcare sectors and export partners worldwide. Expertly working to our specialisms, whilst promoting education and continual development of our employees, we maintain our dedication to our safety, unique skill, products and services.

Surgical Holdings' Vision is to be recognised by our stakeholders as the most passionate, authentic and knowledgeable surgical instrument manufacturer and service provider in the UK.

Values:





TRANSPARENCY F

PASSION



2017

Rigid Endoscope Repair Workshop opens, appointed with innovative technical equipment to facilitate the highest quality repairs. Our export base continues to grow, with appointment of a key partner in Malaysia.

We start working with one of our key brands Symmetry Surgical, bringing products like the Bookwalter Retractor and the Thomas Codman instruments to our customers.



2019

Surgifix is made available to customers for the first time, as our repair service is merged and identified as SMART REPAIR, giving customers unique dashboard access.

Version 4 of our catalogue is published.

We exhibit at Medica, this time promoting our rigid endoscope repair service for our export marketplace.



2021 Orthopaedic Power Tool Repair

Workshop opens, offering repair of Large Bone Battery Hand Pieces.

Surgical Holdings Accredits their first CPD approved workshop for continuing professional development in Surgical Instruments and medical devices.

Employment of our first marketing executive, to support customers and our sales team.

<u>e</u>

Promoting education; haring knowledge and training opportunities

2022 We establish our modern day driven themes and company values, putting education at the forefront to trainees, apprentices, graduates, and industry enthusiasts.

We also expand our premises into neighbouring unit 10 following Southend gaining a city status. Investment in new state-of-the-art

cleaning line, and laser marking machines. Appointment of additional sales team member in South East & London and Scotland.

TrustLinks

Established a partnership with TrustLinks, South East Essex's local independent charity for mental health and wellbeing.

Successful recertification audit for ISO 13485 and CE Marking for orthopaedic implants.

We have successfully launched our CPD Training Procedures including our Think Twice Educational Programme, enabling our qualified Representatives to host Rigid Endoscope Fault Identification training courses and conduct Scope Audits.

ļ.

2016 We begin working with operating rooms in Northern Ireland, extending our sales coverage across the whole of the UK.



2018

2018 Our focus on quality continues, as we appoint a new quality manager to drive the business forward.

SMXRT REPAIR

With success from our Surgifix platform, SMART REPAIR is put to the forefront, tying in with our company focus of sustainability and ethical duty. This unique repair service promises to replace an instrument with a new one if we can't repair an instrument to 'as new' condition.

THICAL DUT

We complete the 5000th in-house rigid scope repair.

2021

At the present day, our once apprentices are now established as Operations Manager, Assistant Factory Manager, and our most specialist Repair Technicians.



And now... Dan, Luke, Reece, James and Dean



Looking ahead We continue to invest and train our

We continue to invest and train or instrument makers of the future.

Our team continues to grow, with us now employing 50 people in our Southend city base, enabling us to support our customers' needs. Please add and further invest in the company's Social Value.

Customers are valuing British made instrumentation, as we strive to support NHS and private hospitals to recover from waiting lists as a result of Covid.

All our company vehicles for representatives are either fully electric or hybrid.

Proud to be a family-run British company

Instruments are repaired to BS 5194 where possible, OEM specification in accordance with the recommendations laid out in MHRA 'Managing Medical Devices,' ensuring CE marking remains valid.

REPAIRED

High quality materials are used, providing full traceability as every part is recorded. Replacing instruments is avoided wherever possible.

Independent testing helps benchmark repairs against industry standards.

CONCUPERT

If an instrument is not to our high standards or we can't repair it, we will replace it with a brand new one for the cost of a repair*. This ensures continual upcycling of your trays.

Instruments can be laser marked, encoded and GS1 2D matrix marked at the same time as a repair.

*Applies to selected instruments only, please get in touch for full details

Repairs are ISO 13485 accredited.

📵 GUARANTEED

Customers receive an industry leading 18 month warranty as standard on all instrument repairs.

We carry full product, professional liability and transit insurance.



Standard turnaround times are 1 week, an express 2 day service is available.

Our instrument technicians carry out a minimum apprenticeship of 5 years and are trained to manufacture instruments rather than simply repair them. This gives us 360 degree understanding of the challenge and the importance of correctly functioning medical devices.

The Repair Process

Request a collection by contacting Surgical Holdings or by using the Surgifix Repair Portal

> Complete paperwork, pack the decontaminated instruments and call Surgical Holdings to arrange a courier or sales representative collection

> > Quotation is raised

Provide purchase order details to Surgical Holdings

Instrument is repaired, polished, laser marked for traceability and passivated

> Returned by courier within 7 days, with a repair report (for scopes and air tools)

+

Our expertise means that instruments can be modified or remanufactured if required. Customers can also receive training, participate in factory visits or simply call their local representative for guidance with instrumentation or management.

Need advice? Our friendly team of experts will be happy to assist you

Instrument and Power Tool Repair

We repair general instruments, laparoscopic and electrosurgery instruments

- Instruments are repaired by an experienced instrument 'maker', not just trained on repairs
- All market brands of leading power, battery and pneumatic tools repaired
- Brand new parts with full traceability
- Full tracking of devices from handover to return through Surgifix
- Option for free personal pick up and delivery through our team







SMXRT REPAIR





SMART REPAIR is our unique repair service for surgical instruments, offering the promise of long-term value for money in choosing the sustainable option, and a product replacement if we are unable to return the instrument to an 'as new' condition*.

In choosing our innovative Smart Repair service, customers benefit from dashboard access to our Surgifix Online Portal.

*For the price of the repair. Applies to the 50 most used and repaired products, correct at the time of going to press,

Surgifix Repair Portal

Designed to make managing repairs easy and quick

Manage and monitor repairs online, and communicate directly with the repairs team.

An easy to use system with up to date information on all repairs

225

5720

On the Surgifix Repair Portal customers can:

View pricing, create guotes, and view £... the progress of repairs. These can be printed or saved.



Historical endoscope repair reports are available.




Instant notifications allow the repairs team and customers to communicate seamlessly regarding what needs a repair or refurbishment.



Training is available and can be done on a video call if required.

Endoscope Repair

All major brands of rigid endoscope can be repaired at Surgical Holdings workshop in Southend-on-Sea. Using the latest optical technology, Surgical Holdings' highly skilled and experienced technicians will return scopes to the manufacturer's original specification. All parts are brand new and have full traceability. A hassle-free and easy endoscope repair service



*MHRA Managing Medical Devices - April 2016

Instrument Refurbishment

Refurbishing instruments can save up to

BU/0 of the cost of buying new,

with compelling evidence to demonstrate effectiveness

Independent research¹ showed that instruments refurbished by Surgical Holdings had a 4x reduction in the amount of residual protein.

The smooth and polished surface has an increased cleaning efficacy, and all refurbishments come with an 18 month warranty as standard. Standard turnaround times are 1 week, and an express 2 day service available.

Reference: 'Data available on file.



Refurbishment

After Refurbishment

Endoscope Audits

Is your scope suitable for use? Think TW CE

[A] **<u>T</u>IP**

Check the tip of the endoscope for sharpness and any type of damage to the tube such as dents and cracks

WINDOWS

Make sure the windows at the tip and evepiece are free from dirt and cracks



ILLUMINATION

Hold the scope lightpost to a light source to see if the output is bright and clear at the tip

○ CLARITY

Look through the eyepiece, is the image clear and free from shadows, dots and lines

B ENVIRONMENT

If you think twice and the endoscope isn't suitable for use then think of the environment and repair instead of replace

We offer a comprehensive scope auditing service

Can we help with your auditing? Please contact your local Surgical Holdings representative to book.

CPD Training Procedures





At Surgical Holdings we have 500 years worth of combined medical device manufacture experience. We're passionate about what we do and are keen to pass that wealth of knowledge onto our customers. We believe it is important to educate and train others, to help them get the most out of their instruments and so they last for many, many years.

We offer Rigid Endoscopy Fault Identification Workshops where we visit you and:



Show you the difference between poor quality and high quality instruments to help you make the right choices when purchasing



Train you on cleaning and maintenance to help keep your instruments in good condition



Walk you through our repair and refurbishment processes, so when the time comes you'll know exactly how your instruments will be cared for



Need advice? Our friendly team of experts will be happy to assist you

We also offer CPD accredited webinars, to boost your ongoing learning and help with your revalidation.

This includes our new module on the **Quality Checking of Rigid Endoscopes**, without the need for expensive checking apparatus. We introduced this module following feedback from stakeholders in sterile services and theatres.

Our team have designed a training SOP that can be integrated into your internal QMS systems, covering the different faults we see and how they can be identified. Once the training module is completed, you will then be issued a certificate directly via CPD. We'll leave you competently trained and externally accredited, so you can then train other staff members.

We are looking forward to working with you and sharing with trusted training model so please get in touch for more information!

Contact your local rep to arrange a training webinar



We also host a series of educational webinars featuring guest speakers and industry professionals

Surgical Instrument Set Refurbishment at North Devon District Hospital

Surgical Holdings have been recently managing a refurbishment programme for North Devon District Hospital. Here is their report on the implementation and success of their refurbishment programme so far.

Background

North Devon District Hospital is part of the Northern Devon Healthcare NHS Trust and the theatre complex has an inventory in excess of 1000 sets. The majority of sets are owned by the trust but this is supplemented with "Loan and Consignment" Orthopaedic sets.

Need for Change

With the need for clinical efficiency and productivity in the theatres instrumentation is required to be efficient and fit for purpose.

Limited sets are being used with a high turnover affecting compliance, sharpness, tension, etc. New instruments are available to replace old, well used stock but the capital required for the instrumentation is hard to obtain.

Implementation

To offset the need for capital to fund instrument replacements the CSSD and Theatres took the decision to try the Surgical Holdings option of refurbishing the instruments sets rather than disposing and purchasing new instruments. Instruments will either be refurbished back to new quality or replaced with new items for an inclusive price.

A "Basic" set was delivered to Surgical Holdings on the understanding the returned set would be assessed by the theatre staff to ensure the returned items met the requirements of the clinical team; e.g. sharpness, tension, colour, glare (lack of), etc.

The set passed the test and the remaining 22 "Basic" sets were refurbished on a rolling plan, the benefit of sending one set away at a time ensured a maximum stock level was available for clinical use. Following the success of the 22 "Basic" sets a further 10 Oral Surgery Sets have been refurbished to the same standard.

Currently we have a Large Ortho Basic back after refurbishment awaiting theatre approval.

Implications

Finance

- No capital required to fund purchase
- Revenue spend to fund refurbishment

CSSD

- Minimum of sets out of service due to rolling programme
- Less incomplete sets as items are like new on refurbished sets, not needing further repairs

Theatre Staff

- Increased Clinical confidence instruments like new
- Not subject to shortage of sets during process
- Reduction in procurement process process agreed then rolling work plan
- · Clinical staff sign-off regarding quality

"

The refurbishment of the theatre sets has allowed high usage instrumentation to be brought back up to a standard equal to new, this has been met with appreciation from theatre staff and has ensured clinical requirement for sets has not been compromised.


Fraser Bishop, Facilities Clinical Services Manager Trust Decontamination Lead, North Devon District Hospital

General Instruments

Spring Forceps19Tissue Forceps20Sponge Holding and Sterilising Forceps21Other Forceps22Scissors22Needle Holders23Combined Scissor and Needle Holders33Bag and Towel Clips34Retractors35Dissectors and Probes44Miscellaneous44Hollow Ware44Silicone Matting44Sterilising Baskets44	Artery Forceps	14
Tissue Forceps20Sponge Holding and Sterilising Forceps20Other Forceps21Scissors22Needle Holders21Combined Scissor and Needle Holders31Bag and Towel Clips34Retractors31Dissectors and Probes44Miscellaneous44Hollow Ware44Silicone Matting44Sterilising Baskets44	Clamps	17
Sponge Holding and Sterilising Forceps2Other Forceps2Scissors2Needle Holders2Combined Scissor and Needle Holders3Bag and Towel Clips3Retractors3Dissectors and Probes4Miscellaneous4Hollow Ware4Silicone Matting4Sterilising Baskets4	Spring Forceps	19
Other Forceps2Scissors2Needle Holders2Combined Scissor and Needle Holders3Bag and Towel Clips3Retractors3Dissectors and Probes4Miscellaneous4Hollow Ware4Silicone Matting4Sterilising Baskets4	Tissue Forceps	20
Scissors2Needle Holders2Combined Scissor and Needle Holders3Bag and Towel Clips3Retractors3Dissectors and Probes4Miscellaneous4Hollow Ware4Silicone Matting4Sterilising Baskets4	Sponge Holding and Sterilising Forceps	22
Needle Holders2Combined Scissor and Needle Holders3Bag and Towel Clips3Retractors3Dissectors and Probes4Miscellaneous4Hollow Ware4Silicone Matting4Sterilising Baskets4	Other Forceps	23
Combined Scissor and Needle Holders3Bag and Towel Clips3Retractors3Dissectors and Probes4Miscellaneous4Hollow Ware4Silicone Matting4Sterilising Baskets4	Scissors	24
Bag and Towel Clips3-Retractors3-Dissectors and Probes4-Miscellaneous4-Hollow Ware4-Silicone Matting4-Sterilising Baskets4-	Needle Holders	29
Retractors33Dissectors and Probes4Miscellaneous4Hollow Ware4Silicone Matting4Sterilising Baskets4	Combined Scissor and Needle Holders	33
Dissectors and Probes 4 Miscellaneous 4 Hollow Ware 4 Silicone Matting 4 Sterilising Baskets 4	Bag and Towel Clips	34
Miscellaneous4Hollow Ware4Silicone Matting4Sterilising Baskets4	Retractors	35
Hollow Ware 4 Silicone Matting 4 Sterilising Baskets 4	Dissectors and Probes	43
Silicone Matting 40 Sterilising Baskets 40	Miscellaneous	44
Sterilising Baskets 4	Hollow Ware	45
	Silicone Matting	46
Optical Fibre Cables4	Sterilising Baskets	46
	Optical Fibre Cables	48

Artery Forceps

Adson Artery Forceps




Serrated jaws

Code	Description
SC-GI/034	Straight
SC-GI/034a	Curved

1	I-2 teeth	
	Code	Description
	SC-GI/034b	Straight
	SC-GI/034c	Curved

Cushing Artery Forceps


Box joint, fine, half serrated jaws, 140mm long





Code	Description	
SC-GI/001	Straight	
SC-GI/002	Curved	

Cairns Artery Forceps




Dunhill Artery Forceps

Box joint, fine serrations




Straight Curved				
Code	Description		Code	Description
SC-GI/065	135mm long		SC-GI/043	135mm long
SC-GI/065a	190mm long		SC-GI/468	180mm long

Crile Artery Forceps



Lawson Tait Artery Forcep

Box joint, straight jaw



Halstead Mosquito Forceps

Standard pattern, haemostatic, box joint





Descri

SC-GI/016a Straight

SC-GI/017b Curved

1-2 teeth

1-2 teeth		
Code	Description	
SC-GI/010	Straight 125mm long	
SC-GI/012	Straight 190mm long	
SC-GI/011	Curved 125mm long	
SC-GI/013	Curved 190mm long	

Heiss Artery Forceps

Box joint, 205mm long

Descriptio

Straight SC-GI/017 Curved

Serrated jaws

Code

SC-GI/016


Micro, 100mm long			
Code	Description		
SC-GI/003	Straight, serrated		
SC-GI/004	Curved, serrated		
SC-GI/005	Straight, 1-2 teeth		
SC-GI/006	Curved, 1-2 teeth		

Kelly Artery Forceps

Box joint, half serrated jaws

Straight			Curved	
Code	Description		Code	Description
SC-GI/022	140mm long		SC-GI/024	140mm long
SC-GI/023	160mm long		SC-GI/025	160mm long

Fraser Kelly Artery Forceps Curved



Kocher Artery Forceps

Box joint, 1-2 teeth



8.....

Straight

otrangine		
Code	Description	Code
SC-GI/052	130mm long	SC-GI/05
SC-GI/053	150mm long	SC-GI/05
SC-GI/054	180mm long	SC-GI/05
SC-GI/055	200mm long	SC-GI/06
SC-GI/056	230mm long	

Curved

Code	Description
SC-GI/057	150mm long
SC-GI/058	180mm long
SC-GI/059	200mm long
SC-GI/060	230mm long


General Instruments

Moynihan Artery Forceps

Box joint, half-serrated jaws




Rochester-Pean Artery Forceps



Code	Description	Code	Description
SC-GI/044	165mm long	SC-GI/048	165mm long
SC-GI/045	190mm long	SC-GI/049	180mm long
SC-GI/046	205mm long	SC-GI/050	205mm long
SC-GI/047	240mm long	SC-GI/051	240mm long

Roberts Artery Forceps

Screw joint, curved on flat, serrated



Dandy Artery Forceps

Box joint, curved in plane




Spencer Wells Artery Forceps

Box joint, serrated jaws



OTHER SIZES AVAILABLE

Carmalt Angiotribe Forceps



Clamps

Dieffenbach Bulldog Clamp



Straight, serrated

Code	Description	
SC-GI/112	1.2mm fine jaw, 40mm long	
SC-GI/113	1.2mm fine jaw, 45mm long	
SC-GI/108	1.8mm jaw, 30mm long	
SC-GI/109	1.8mm jaw, 35mm long	
SC-GI/110	2.5mm jaw, 40mm long	
SC-GI/111	2.5mm jaw, 45mm long	
SC-GI/107	6mm heavy jaw, 60mm long	
SC-GI/455	20mm jaw, 75mm long	
SC-GI/455a	30mm jaw, 85mm long	
SC-GI/455b	45mm jaw, 105mm long	
SC-GI/455c	65mm jaw, 120mm long	

Straight, plain			
Code	Description		
SC-GI/129	1.2mm fine jaw, 40mm long		
SC-GI/131	1.2mm fine jaw, 45mm long		
SC-GI/121	1.8mm jaw, 30mm long		
SC-GI/123	1.8mm jaw, 35mm long		
SC-GI/125	2.5mm jaw, 40mm long		
SC-GI/127	2.5mm jaw, 45mm long		



Curved, serrated

Code	Description
SC-GI/119	1.2mm fine jaw, 40mm long
SC-GI/120	1.2mm fine jaw, 45mm long
SC-GI/115	1.8mm jaw, 30mm long
SC-GI/116	1.8mm jaw, 35mm long
SC-GI/117	2.5mm jaw, 40mm long
SC-GI/118	2.5mm jaw, 45mm long
SC-GI/114	6mm heavy jaw, 60mm long
SC-GI/456	20mm jaw, 75mm long
SC-GI/456a	30mm jaw, 85mm long
SC-GI/456b	45mm jaw, 105mm long
SC-GI/456c	65mm jaw, 120mm long

Curved, plain

Code	Description
SC-GI/130	1.2mm fine jaw, 40mm long
SC-GI/132	1.2mm fine jaw, 45mm long
SC-GI/122	1.8mm jaw, 30mm long
SC-GI/124	1.8mm jaw, 35mm long
SC-GI/126	2.5mm jaw, 40mm long
SC-GI/128	2.5mm jaw, 45mm long

Blalock Bulldog Clamp



Serrated, with hole in handle

Straight

Code	Description	Code	Description
SC-GI/480	11mm jaw, 38mm long	SC-GI/478	11mm jaw, 38mm lor
SC-GI/482	14mm jaw, 48mm long	SC-GI/479	14mm jaw, 48mm lor
SC-GI/493	17mm jaw, 60mm long	SC-GI/491	17mm jaw, 60mm Ion
SC-GI/481	27mm jaw, 70mm long	SC-GI/492	27mm jaw, 70mm lor

Curved

Plain, with hole in handle

Straight C				Cur
	Code	Description		Coc
	SC-GI/486	11mm jaw, 38mm long		SC-
	SC-GI/487	14mm jaw, 48mm long		SC-
	SC-GI/488	17mm jaw, 60mm long		SC-
	SC-GI/489	27mm jaw, 70mm long		SC-

urved				
Code	Description			
SC-GI/483	11mm jaw, 38mm long			
SC-GI/484	14mm jaw, 48mm long			
SC-GI/485	17mm jaw, 60mm long			
SC-GI/490	27mm jaw, 70mm long			

An Overview On

Forcep Jaw Types

You will see that forceps in our catalogue will have either a traditional serrated jaw or a different design which is known as atraumatic (or non-traumatic).

Traditional Serrated Jaw

The traditional serrated jaw is generally more robust and used on forceps where the areas are not as delicate.

Atraumatic Jaw

Atraumatic jaws will be mainly found in delicate surgery areas such as vascular and cardiothoracic. The design is such, that the surface area gripped is reduced, ensuring less trauma without the loss of grip.

Spring Forceps


Adsons Dissecting Forceps **DeBakey Atraumatic Dissecting Forceps** Straight, atraumatic jaws 127mm long Description 1.5mm 2.5mm SC-GI/257 1-2 teeth Code Description Code Description SC-GI/258 Lightweight, micro SC-GI/259 Serrated SC-GI/238a 150mm long SC-GI/238b 150mm long SC-GI/240b 240mm long SC-GI/260 Serrated, lightweight, micro SC-GI/240a 240mm long 2mm 2.8mm **Block End Dissecting Forceps** Code Description Code Description SC-GI/238 150mm long SC-GI/238c 150mm long SC-GI/241 180mm long SC-GI/240c 240mm long SC-GI/241c 180mm long **Dissecting Forceps** Serrated Turnover end, serrated jaws Description SC-GI/255 127mm long SC-GI/256 152mm long SC-GI/256a 203mm long **Bonney Dissecting Forceps** Code Description Code Description SC-GI/242 127mm long SC-GI/244 178mm long SC-GI/243 152mm long SC-GI/245 203mm long **Dunhill Dissecting Forceps** Light, fine serrated jaws 178mm long Code Description SC-GI/267 Serrated, 1-2 teeth SC-GI/268 1-2 teeth Code Description SC-GI/457 145mm long **Charnley-McIndoe Dissecting Forceps** Fredrick Coole Dissecting Forceps 150mm long Extra fine serrated jaws 125mm long Code Description Code Description

SC-GI/066 Straight SC-GI/067 Curved

SC-OR/N045a 2-3 teeth

Gillies Dissecting Forceps

1-2 teeth



Jackson Burrows Dissecting Forceps


140mm long



Code	Description
SC-OR/N037	Tongue and groove, 1-2 teeth
SC-OR/N037a	Tongue and groove only

Jean Dissecting Forceps

Serrated jaws, 140mm long



SC-GI/269	Fine
SC-GI/270	1-2 teeth
SC-GI/271	2-3 teeth, 178mm long

Description

127mm long

152mm long

178mm long

Code

SC-GI/246

SC-GI/247

SC-GI/248


Lanes Dissecting Forceps

1-2 teeth on serrated jaws, 140mm long



McIndoe Dissecting Forceps

152mm long




McIndoe Modified Dissecting Forceps

152mm long, 1mm tip, modified for fine powerful dissection with strong action tip




SH/5961 Serrated

Treves Dissecting Forceps



SC-GI/251 127mm long SC-GI/252 127mm long SC-GI/252 152mm long SC-GI/254 152mm long				
SC-GI/252 152mm long SC-GI/254 152mm long	SC-GI/251	127mm long	SC-GI/253	127mm long
	SC-GI/252	152mm long	SC-GI/254	152mm long

Waugh Dissecting Forceps



Description


SC-GI/249 2-3 teeth 178mm long

SC-GI/250 2-3 teeth 203mm long

Tissue Forceps


Allis Tissue Forceps

Box joint



SC-GI/414 5-6 teeth


Duval Tissue Forceps



Littlewood Tissue Forceps

Box joint, 2-3 teeth

SC-GI/273 4-5 teeth



Lane Tissue Forceps

Box joint, 1-2 teeth, fenestrated




Babcock Tissue Forceps

Box joint, serrated jaws



Rutherford Morrison Tissue Forceps


Box joint



Tissue Forceps

Shaw Tissue Forceps

Box joint



Guys Gallabin Suture-Holding Forceps

Screw joint, angled flat blades



Russian Dissecting Forceps

With fenestrated end



Code	Description
SC-GI/SH-5503b	150mm long
SC-GI/SH-5503	200mm long
SC-GI/403	230mm long
SC-GI/SH-5503a	250mm long



Special Offer



Email us to place your order

info@surgicalholdings.co.uk

General Instruments

Sponge Holding and Sterilising Forceps

Rampley Sponge Holding Forceps

Box joint, straight, serrated jaws




Box joint, serrated tip



Code	Description
SC-GI/353	178mm long
SC-GI/354	241mm long

Foerster	Sponge	Holding	Forceps
----------	--------	---------	---------


Box joint, straight, serrated jaws





Cheatle Sterilizer Forceps

Screw joint






Harrison Bowl Steriliser Forceps

Screw joint, angled, double jaw

Code	Description
SC-GI/365	356mm long
SC-GI/366	457mm long



Hunter Splinter Forceps

Serrated grooved pointed jaws

Michel Forceps

Box joint, double-ended, for inserting and removing clips stainless steel



Code	Description
SC-GI/357	130mm long

Michel Skin Holding Forceps

With magazine for Michel clips


Code Description SC-GI/360 115mm long



Code	Description
SC-GI/358	1-2 teeth, 127mm long

Childe Approximating Forceps


With magazine for Michel clips



Code	Description
SC-GI/469	180mm long

Michel Clip Extracting Forceps

Box joint



Robins Anchoring Forceps

140mm long



Code	Description
SC-GI/453	5mm with holder for diathermy leads
SC-GI/453a	13mm with holder for suction tubing


Premium Quality Stainless Steel Scissors

Dressing Scissors Metzenbaum Scissors C **Blunt/Blunt** Straight Curved Straight Curved Code Descr Code Code Description Description SC-GI/133 130mm long SC-GI/142 130 SC-GI/405 SC-GI/407 180mm Iona 185mm Iona SC-GI/136 146mm long SC-GI/143 146mm long SC-GI/406 147mm long SC-GI/408 147mm long SC-GI/139 165mm long SC-GI/144 165mm long SC-GI/406a 200mm long SC-GI/408a 200mm long SC-GI/145 185mm long SC-GI/406b 230mm long SC-GI/408b 230mm long SC-GI/404 147mm long SC-GI/406c 250mm long SC+GI/408c 250mm long With one serrated blade SC-GI/406d 280mm long SC-GI/408d 280mm long Blunt/Sharp Sharp/Sharp Straight Straight Code Code **Mayo Scissors** SC-GI/134 130mm long SC-GI/135 130mm long SC-GI/137 SC-GI/138 146mm long 146mm long SC-GI/140 165mm long SC-GI/141 165mm long Curved Code Description SC-GI/146 130mm long Straight Curved **Dressing Scissors, Left-Handed** Description Code With one black bow, 145mm long SC-GI/409 SC-GI/410 140mm long 140mm long SC-GI/409a 160mm long SC-GI/410a 160mm long **Blunt/Blunt** SC-GI/411 190mm long SC-GI/412 190mm long SC-GI/411a 230mm long SC-GI/412a 230mm long Code Descriptio SC-GI/192 Straight SC-GI/193 Curved **Fergusson Scissors Blunt/Sharp** Sharp/Sharp Angled on flat Descrip Code Desc SC-GI/194 Straight SC-GI/196 Straight SC-GI/195 Curved SC-GI/197 Curved Code Descrip SC-GI/148 180mm long Littler Dissecting Scissors Fine pattern, with eye for suture **Pavnes Suture Scissors** TA Code Description SC-PL/001f 115mm long SC-GI/152 100mm long **Kilner Dissecting Scissors** Ligature Stitch Scissors Fine pattern Fine pattern Ð Code Description Description SC-PL/001a 115mm long SC-GI/153 115mm long

Lister Bandage Scissors




McIndoe Scissors



Code	Description	
SC-GI/149	Straight	
SC-GI/150	Curved on flat	
SC-GI/150sc	Supercut, curved	

Heath Stitch Scissors



Tubing Scissors




Lister Bandage Scissors

Angled, with probe end



Thomas Strabismus Scissors

114mm long



 Code
 Description

 SC-OP/030
 Straight

 SC-OP/029
 Curved

"Careless" Scissors





Tungsten Carbide Hard Insert Scissors

Tungsten carbide hard insert, blunt/blunt



Straight		Curved		
Code	Description	Code	Description	
SC-GI/430	130mm long	SC-GI/431	130mm long	
SC-GI/432	145mm long	SC-GI/433	145mm long	
SC-GI/434	165mm long	SC-GI/435	165mm long	
SC-GI/436	185mm long			

Mayo Dressing Scissors

Tungsten carbide hard insert				
Straight	Straight Curved			
Code	Description	Code	Description	
SC-GI/419	140mm long	SC-GI/420	140mm long	
SC-GI/419a	160mm long	SC-GI/420a	160mm long	
SC-GI/421	190mm long	SC-GI/422	190mm long	
SC-GI/421a	230mm long	SC-GI/422a	230mm long	

Dressing Scissors

Tungsten carbide hard insert, blunt/sharp				
Straight Curved				
Code	Description		Code	Description
SC-GI/423	130mm long		SC-GI/424	130mm long
SC-GI/425	145mm long		SC-GI/426	145mm long

Mayo-Harrington Dressing Scissors

Tungsten carbide hard insert



Dressing Scissors

SC-GI/427 165mm long

SC-GI/429 185mm long

Tungsten carbide hard insert sharp/sharp




SC-GI/428 165mm long

Straight		Curved		
Code	Description	Code	Description	
SC-GI/437	130mm long	SC-GI/438	130mm long	
SC-GI/439	145mm long	SC-GI/440	145mm long	
SC-GI/441	165mm long	SC-GI/442	165mm long	
SC-GI/443	185mm long			


Baby Metzenbaum Scissors

Tungsten carbide, hard insert, 145mm long

Code



Metzenbaum Scissors



long

long

long

long long

-61/445	145mm long	SC=01/440	145111111
-GI/448	185mm long	SC-GI/449	185mm
-GI/471a	230mm long	SC-GI/449a	200mm
-GI/471b	250mm long	SC-GI/449b	230mm
-GI/471c	280mm long	SC-GI/449c	250mm
		SC-GI/449d	280mm



SC

SC

SC

SC

Tungsten Carbide Supercut Scissors

Mayo Harrington Scissors



Mayo Scissors

Supercut, tungsten carbide



Metzenbaum Dissecting Scissors

Supercut, tungsten carbide



Straight		Curved	
Code	Description	Code	Description
SC-GI/188tc	185mm long	SC-GI/191tc	147mm
SC-GI/188tca	200mm long	SC-GI/189tc	185mm
SC-GI/188tcb	230mm long	SC-GI/186tca	200mm
SC-GI/188tcc	250mm long	SC-GI/188tcb	230mm
SC-GI/188tcd	280mm long	SC-GI/188tcc	250mm
		SC-GI/188tcd	280mm

Dressing Scissors

Supercut, tungsten carbide



Blunt/Blunt

Straight		Curved	
Code	Description	Code	Description
SC-GI/174	130mm long	SC-GI/175	130mm long
SC-GI/176	145mm long	SC-GI/177	145mm long
SC-GI/178	165mm long	SC-GI/179	165mm long
SC-GI/180	185mm long	SC-GI/180a	185mm long

Blunt/Sharp

Straight				
Code	Description			
SC-GI/167	130mm long			
SC-GI/169	145mm long			
SC-GI/171	165mm long			
SC-GI/173	185mm long			

c	u	r	V	e	d	

Code	Description
SC-GI/168	130mm long
SC-GI/170	145mm long
SC-GI/172	165mm long
SC-GI/173a	185mm long

Sharp/Sharp

Straight		Curved	
Code	Description	Code	Description
SC-GI/181	130mm long	SC-GI/182	130mm long
SC-GI/183	145mm long	SC-GI/184	145mm long
SC-GI/185	165mm long	SC-GI/186	165mm long
SC-GI/187	185mm long	SC-GI/187a	185mm long

Guide for Gold and Black Bows

It's important to know the difference for the intended use of each scissor:



1 Black and 1 Gold Bow Supercut tungsten carbide 1 Black Bow (left ring) Left-handed instrument

Supercut Operating Scissors

Dressing Scissors



Dressing Scissors



Dressing Scissors







Metzenbaum Dissecting Scissors

Supercu	t			0		
Straight (I	ength)		Curved (le	ngth)		
Straight (l	ength) Description		Curved (le	ngth) Description		
Code	Description		Code	Description		
Code SC-GI/190	Description 147mm		Code SC-GI/191	Description 147mm		
Code SC-GI/190 SC-GI/188	Description 147mm 185mm		Code SC-GI/191 SC-GI/189	Description 147mm 185mm		
Code SC-GI/190 SC-GI/188 SC-GI/188a	Description 147mm 185mm 200mm		Code SC-GI/191 SC-GI/189 SC-GI/186a	Description 147mm 185mm 200mm		

Stevens Tenotomy Scissors




Curved/sharp

Curved/blunt
Code Description

Supercut

Code	Description	Code	Description
SC-OP/424	114mm long	SC-OP/425	114mm long
SC-OP/424a	125mm long	SC-OP/425a	140mm long
SC-OP/424gb	135mm long	SC-OP/425b	160mm long

Code	Description	With cutout, 110mm long	
SC-OP/426	Straight/blunt, 114mm long	Code	Description
SC-OP/427	Straight/sharp, 114mm long	SC-OP/427a	Straight/blunt
SC-OP/427c	Straight/sharp, 180mm long	SC-OP/427b	Curved/blunt





Needle Holders

Bruce Clarke Needle Holder

Screw joint



Lawrence Needle Holder

Fine jaw, 150mm long



Crile-Wood Needle Holder

Straight, baby, fine



Code Description SC-GI/199 150mm long

Kilner Needle Holder

Box joint, cross-serrated jaw



Code	Description
SC-GI/200	165mm long
SC-GI/201	133mm long

Needle Holder

Surgical Holdings pattern, box joint, lightweight, with curved plain jaws and one offset bow for thumb



MacPhail Needle Holder

Screw joint, with one copper jaw



Code	Description
SC-GI/206	159mm long
SC-GI/207	178mm long
SC-GI/208	209mm long

Needle Holders

Mathieu Needle Holder **Naughton Morgan Needle Holder** Screw joint Box joints, cross serrated jaws Code SC-GI/210 172mm long SC-GI/211 209mm long 243mm long SC-GI/212 Tungsten carbide **Mayo Needle Holder** Code Description Description Box joint, straight SC-GI/401 185mm long SC-GI/399 185mm long SC-GI/402 254mm long SC-GI/400 254mm long **Rochester-Mayo Needle Holder** Screw joint CONTINUE (1) Wide cross serrated jaws Narrow cross serrated jaws Code Description Description Code D 159mm long SC-GI/225 159mm long SC-GI/213 178mm long SC-GI/224 SC-GI/398 203mm long

Choosing the most suitable needle size

Our tungsten carbide inserts are available in a variety of styles. Choosing the most suitable needle size for each procedure helps give the surgeon maximum grip and control.

Coarse 0.65mm

For powerful needle holder patterns and wire twisters. Maximum holding capabilities for heavy needles or wire.

Standard 0.5mm Used, for example, in Mayo-Hegar Mathieu Needle Holders. For needles 3-0 and larger.

Fine 0.4mm For delicate needle holders such as Crile-Wood and DeBakey. For needles and sutures between 4-0 and 6-0.

Micro 0.2-0.3mm

For micro needle holders such as Castroviejo, Microvascular and Mini-Ryder. For delicate needles and sutures between 7-0 and 10-0.

Smooth


mooth tungsten carbide surface educes damage to delicate needles nd suture of 6-0 and finer.



Need advice? Our friendly team of experts will be happy to assist you


Berry Needle Holder/Wire Twister

Tungsten carbide jaws




Finochietto Needle Holder

Tungsten carbide jaws



Crile-Wood Needle Holder

Tungsten carbide insert, box joint, satin finish, straight




Halsey Needle Holder

Tungsten carbide insert, box joint, satin finish, straight



Derf Needle Holder




Kilner Needle Holder

Tungsten carbide insert

SC-GI/221 140mm long

Tungsten Carbide Needle Holders

Masing Needle Holder



Nievert Needle Holder




DeBakey Needle Holder

250mm long



Ryder Needle Holder

Tungsten carbide insert, box joint, satin finish, straight




0.2mm jaw pitch

Code	Description	Code	Description
SC-GI/SH-5197a	135mm long	SC-GI/SH-5198a	135mm long
SC-GI/SH-5197b	155mm long	SC-GI/SH-5198b	155mm long
SC-GI/SH-5197c	175mm long	SC-GI/SH-5198c	175mm long
SC-GI/SH-5197d	210mm long	SC-GI/SH-5198d	210mm long

0.4mm jaw pitch

Sarot Needle Holder

Diamond jaw, tungsten carbide, gold plated bows



Tungsten Carbide Needle Holders

Features and Benefits

Jaw Fin Tungsten carbide inserts tu bonded into the jaws for a long lasting, secure grip Bevelled Edges Va Prevents suture from snagging OC Laser Etching Date and catalogue number are etched for easy re-ordering. Laser etch will not harm

are etched for easy re-order Laser etch will not harm instrument like traditional methods

Satin Finish -Reduces glare Fine, hard teeth, diamond cut into an insert of hardened tungsten carbide for superior needle grasping and suture tying capabilities

Teeth

 Vacuum Brazing Combines bonding of the tungsten carbide insert with the simultaneous vacuum hardening of the instrument for greater durability and longevity


Shanks Precision balanced for comfort and control

Gold Handles

Ratchets Precisely cut for effortless locking/unlocking Gold Handles For easy identification

Mayo-Hegar Needle Holder

Tungsten carbide insert, box joint, satin finish, straight




Combined Scissor and Needle Holders

Olsen-Hegar Scissor/Needle Holder


With tungsten carbide insert

Code



Olsen-Hegar Scissor/Needle Holder

Screw joint



Gillies Scissor/Needle Holder

Right-hand, with tungsten carbide insert



Heavy pattern, 158mm long

Gillies Scissor/Needle Holder



Gillies Scissor/Needle Holder

158mm long



Foster Gillies Scissor/Needle Holder

Fine pattern



Bag and Towel Clips

Spring Clip

Serrated pattern, stainless steel, for fixing slings to leg splints, also for towel clipping



Code	Description
SC-GI/344	63mm wide

Spring Bag Clip

Serrated






Jones Towel Clip

Superior quality spring forcep type



Backhaus Towel Clip

Superior quality, box joint



Moynihan Towel Clip

Screw joint, 2-2 teeth, for sides of incision



 Code
 Description

 SC-GI/470
 191mm long

Ball and Socket Towel Forceps

Box joint

SC-GI/350 89r SC-GI/351 115r

SC-GI/352

140mm long



Code	Description
SC-GI/346	89mm long

Shardles Towel Clip

Spring curved, cross action




Code	Description
SC-GI/347	77mm long

Retractors

Goligher Retractor

Complete with 2 centre blades: 1 x 90mm wide x 50mm deep 1 x 90mm wide x 90mm deep

and 2 pairs of lateral blades: 1 pair 25mm wide x 45mm deep 1 pair 50mm wide x 90mm deep



SC+GI/305

Balfour Abdominal Retractor

Trivalve, frame 203mm

Child


102mm opening, 51mm deep lateral blades, 38mm deep centre blade

Adult

178mm opening, 64mm deep lateral blades, 44mm deep centre blade

Extra large

178mm opening, 102mm deep lateral blades. 108mm deep centre blade




Code	Description
SC-GI/373	Child
SC-GI/374	Adult
SC-GI/375	Extra large

Miles Retractor

Complete with 2 centre blades: 1 x 90mm wide x 50mm deep 1 x 90mm wide x 75mm deep

and 2 pairs of lateral blades: 1 pair 80mm wide x 35mm deep 1 pair 80mm wide x 75mm deep



Bivalve, baby

SC-GI/372 Fully open 89mm



SC-GI/461 Blades 35mm x 60mm deep, 215mm long

Turner Warwick Retractors

Turner Warwick Universal Self Retaining Retractor

Code	Description
SC-GI/520	With frames A/B/C/D and blades A/B/C/D/E/F/K/L
SC-GI/521-A	Circular frame A only 178mm
SC-GI/521-B	Circular frame B only 229m
SC-GI/521-C	Circular frame C only 305mm
SC-GI/521-D	Oval frame D only 475mm x 330mm
SC-GI/521-E	Oval frame E only 370mm x 300mm

Turner Warwick Needleholder with Tungsten Carbide Jaws

Curved shanks with one short bow


SH-5835 23	Omm long
SH-5835A 20	Omm long

Blade A




Code	Description
SC-GI/522-A	Blade only - 19mm x 76mm


Blade C



Code	Description
SC-GI/522-C	Blade only - 57mm x 89mm




Blade B



SC-GI/522-B Blade only - 38mm x 102mm

Blade D




Turner Warwick Retractors

(faal

Blade E

Code	Description
SC-GI/522-E	Blade only - curved stainless steel blade 51mm x 38mm

Blade G



Blade J	
---------	--


Blade F



SC-GI/522-F Blade only - curved stainless steel blade 76mm x 60mm

Code	Description
SC-GI/522-G	Blade only - flat stainless steel blade 89mm x 60mm


Blade K



Code	Description	
SC-GI/522-K	Blade only - flat with clip 55mm x 90mm	

Code	Description
SC-GI/522-J	Blade only - single hook, 13mm deep

Blade L




Code Description
SC-GI/522-L Blade only - flat with clip 55mm x 120mm

Retractors

Collin Retractor

110mm fully open, 235mm long



Joll Thyroid Retractor

Self retaining, with lock-in nut



SC-GI/363 152mm long

Code SC-GI/322 Complete with centre blade 55 x 40mm and one pair of blades 38mm x 60mm SC-GI/322a Without centre blade

Gelpi Retractor

Descriptio

Self retaining




()i=		
Code	Description	
SC-GI/303	Large, fully open 185mm long	
SC-GI/304	Small, fully open 95mm long	

Gossett Abdominal Retractor


Code Desc SC-GY/013a 140mm long SC-GY/013b 178mm long

Pozzi Abdominal Retractors




Bowlby Dilator

Curved




Travers Retractor

Self retaining with internal rack, 210mm long



Irwin Retractor

Self retaining, with cam rack, 4 x 4 teeth



Code	Description
SC-GI/324	178mm long

Norfolk and Norwich Retractor


4 x 5 teeth 50mm blades x 223mm long



Code	Description	
SC-GI/325	Standard rack	

West Retractor

Self retaining, 140mm long



 Code
 Description

 SC-GI/326
 3 x 4 teeth, 60mm fully open

 SC-GI/326cr
 3 x 4 teeth, 60mm fully open, with cam rack

Weitlander Retractor

3-4 teeth



135mm long				
	Code	Description		
	SC-GI/496	Blunt		
	SC-GI/496a	Sharp		
	SC-GI/496b	Sharp, 2-3 teeth		
	SC-GI/496c	Blunt, 2-3 teeth		

10511111101	9
Code	Description
SC-GI/497	Blunt
SC-GI/497a	Sharp

165mm long

Kocher Retractor



230mm long

 Code
 Description

 SC-GI/378
 36 x 20mm blade

250mm long			
Description			
67 x 25mm blade			
80 x 42mm blade			
80 x 50mm blade			
80 x 60mm blade			

Kocher Langenbeck Retractor

215mm long



Code	Description	Code	Description	
SC-GI/336sp	25 x 6mm blade	SC-THO/cutsp015	55 x 11mm blade	
SC-GI/379a	35 x 8mm blade	SC-GI/379d	70 x 15mm blade	
SC-GI/379b	35 x 11mm blade	SC-GI/379e	80 x 12mm blade	
SC-GI/379c	35 x 15mm blade	SC-GI/379f	80 x 16mm blade	
SC-GI/357sp	41 x 11mm blade			

Langenbeck Retractor


Single ended



Code	Description
SC-GI/336	Blade 22mm x 6mm
SC-GI/337	Blade 45mm x 13mm
SC-GI/337a	Blade 65mm x 13mm
SC-GI/338a	Blade 75mm x 25mm
SC-GI/338a1	Blade 65mm x 25mm

Canny Ryall Navy Retractor

Double-ended, with one single hook blade 5mm wide, and one fenestrated blade




\d

couc	beschption
SC-GI/343	Fenestrated blade 16mm wide, 190mm long
SC-GI/343a	Fenestrated blade 19mm wide, 190mm long
SC-GI/343b	Fenestrated blade 25mm wide, 220mm long
SC-GI/343c	Fenestrated blade 30mm wide 220mm long

Canny Ryall Army Retractor


Double prong, 13mm wide, 165mm long



Czerny Retractor



Landon Retractor



Lane Retractor

Double ended blade, 229mm long



Retractors

Morris Retractor

Single blade, 229mm long

Morris Retractor

Double ended, 216mm long

C)				
				Code	Description
Code	Description			SC-GI/327	32 x 51mm and 38 x 51m
SC-GI/330				SC-GI/328	44 x 51mm and 51 x 51m
SC-GI/331				SC-GI/329	64 x 51mm and 38 x 51n
SC-GI/332					
Moy	niham Retrac	tor	Deaver Abdom Retractors	inal	



Elliot (Modified Deaver) Retractor



Deaver Abdominal Retractors

With round handles



Code	Description
SC-GI/312	Set of 5, 305mm overall length
SC-GI/313	25mm wide
SC-GI/314	25mm wide
SC-GI/315	38mm wide
SC-GI/316	51mm wide
SC-GI/317	76mm wide

Dyball and Lock Retractor



Code

SC-GI/306

SC-GI/307

SC-GI/308

SC-GI/309

SC-GI/310

SC-GI/311

Description

25mm wide

25mm wide

38mm wide

51mm wide

76mm wide

Set of 5, 315mm overall length

Kelly Retractor Blade

267mm long



Retractors

Durham Retractor

SC-GI/302 20mm wide

With fenestrated blade, 216mm long

G		
Code	Description	
SC-GI/300	10mm wide	
SC-GI/301	14mm wide	

Ballance (Volkmann) Retractor

203mm long



Sharp		Blunt	
Code	Description	Code	Description
SC-GI/463	1 prong	SC-GI/463e	1 prong
SC-GI/463a	2 prong	SC-GI/463f	2 prong
SC-GI/463b	3 prong	SC-GI/463g	3 prong
SC-GI/463c	4 prong	SC-GI/463h	4 prong
SC-GI/463d	6 prong	SC-GI/463i	6 prong

Volkmann Retractor



Code	Description
SC-OR/400	Sharp, 220mm long
SC-OR/400a	Blunt, 220mm long
SC-OR/400b	Sharp, 245mm long

Kocher Retractor



Code	Description
SC-OR/N036c	Sharp, 205mm long
SC-OR/N036	Sharp, 220mm long
SC-OR/N036b	Blunt, 220mm long

Senn-Miller Retractor

170mm long



Single Hook Retractor

Sharp

Code	Description
SC-GI/SH5671	152mm long
SC-PL/065b	Half size, fine


Gillies Skin Hook

Solid handle, 190mm long

Code	Description
SC-PL/065a	Standard, fine
SC-PL/065b	Half size, fine

Kilner Retractor

Double ended, cat's paw



Dissectors and Probes

MacDonald Dissector Double Ended



Brodie Director Probe Ended, Silver Plated



Watson Cheyne Dissector

Fine, double-ended, with probe end

Probe, Sinus Silver

-		
Code	Description	
SC-GI/464	152mm long	

Probe, Fish Tail Silver Plated



Code	Description
SC-GI/293	127mm long
SC-GI/294	152mm long
SC-GI/295	203mm long
SC-GI/296	254mm long

Jobson Horne Probe Dressing (Aural)

With one flat tapered end, serrated edge and one small ring end, stainless steel

Code	Description
SC-GI/287	127mm long
SC-GI/288	191mm long
SC-GI/288a	230mm long
SC-GI/288b	255mm long

Probe, Malleable Silver Plated

Code	Description
SC-GI/289a	127mm long
SC-GI/289	152mm long
SC-GI/290	178mm long
SC-GI/291	203mm long
SC-GI/292	255mm long
SC-GI/292a	305mm long

Bookwalter[®] Retractor Kit III

A table-fixed retractor system designed by surgeons to efficiently and reliably hold surgical instruments in place during surgery while allowing greater visibility.




Kit includes Bookwalter Kit: Table Post and Coupling, Horizontal Flex Bar, Segmented Rings, Tilt Ratchet Mechanism, Kelly Retractors, Blades and Ratchets

Contact us for full details at info@surgicalholdings.co.uk

Miscellaneous

Aneurysm Needle

Small, for infusions



Syme Aneurysm Needle






Ring Cutter



Mayo Instrument Pin



Code Description SC-GI/355 102mm long SC-GI/356 152mm long

Baron Scalpel Handle

130mm long



000000

82

Code	Description
SC-PM/003	No 3
SC-PM/004	No 4

Scalpel Handle

127mm long

Description
No 3
No 4
No 4 L
No 7

Endoscope Re-fibre

Surgical Holdings manufacture their own fibre bundles in house. We cover small and larger diameter scopes and only use high quality German fibre optic material and stainless steel tubing.

Quality components, globally sourced: Vew solderable sapphire

- ✓ Medical grade stainless steel tubing ✔ German engineered new
- ✓ New fibre cone as standard
- distal window
- ✓ New light post colour ring
- ✓ New objective in 95% of cases
- Replacement 'O' ring to prevent leakage





Before repair

Contact us with your requirements info@surgicalholdings.co.uk

Hollow Ware



Silicone Matting & Sterilising Baskets 5 symmetry surgical

Silicone Matting


Medical grade for storage, sterilisation and transportation of instruments.

Cushion, protect, and organise delicate surgical instruments during storage, transportation, and sterilisation. Perforations allow for ventilation during the sterilisation process. The silicone mat can be cut to fit any size container.




502mm long x 242mm wide 7mm holes, 10mm between holes

Code BES/1004M



515mm long x 240mm wide 3mm holes, 6mm between holes





570mm long x 550mm wide 4mm holes, 13mm between holes



Sterilising Baskets

DIN Basket

Wire mesh base with perforated sides and two carrying handles

Code	Description	Code	Description
SC-GI/CO006	230 x 130 x 60mm	SC-GI/CO013	480 x 250 x 70m
SC-GI/CO007	240 x 130 x 50mm	SC-GI/CO016	480 x 400 x 60n
SC-GI/CO008	240 x 130 x 60mm	SC-GI/CO017	480 x 480 x 50m
SC-GI/CO009	240 x 250 x 50mm	SC-GI/CO018	480 x 480 x 60m
SC-GI/CO010	240 x 250 x 60mm	SC-GI/CO019	500 x 300 x 50m
SC-GI/CO014	300 x 300 x 50mm	SC-GI/CO020) 550 x 430 x 50m
SC-GI/CO015	350 x 350 x 60mm	SC-GI/CO021	600 x 300 x 60n
SC-GI/CO011	480 x 250 x 50mm		
SC-GI/CO012	480 x 250 x 60mm		



DIN Basket

Wire mesh basket with two carrying handles

Code	Description	Code	Description
SC-GI/CO022	240 x 250 x 50mm	SC-GI/CO026	480 x 250 x 60mm
SC-GI/CO023	240 x 250 x 60mm	SC-GI/CO027	480 x 250 x 70mm
SC-GI/CO024	240 x 250 x 70mm	SC-GI/CO028	480 x 250 x 100mm
SC-GI/CO031	300 x 300 x 60mm	SC-GI/CO029	480 x 480 x 50mm
SC-GI/CO025	480 x 250 x 50mm	SC-GI/CO030	480 x 480 x 70mm

DIN Basket Lids

For perforated and wire mesh baskets

Code	Description
SC-GI/CO032	240 x 250mm
SC-GI/CO033	300 x 300mm
SC-GI/CO034	480 x 250mm
SC-GI/CO035	480 x 480mm

Did you know...

Organising and fixing instruments in place can aid the cleaning process.

Help keep your instruments in good condition for longer with our baskets and trays. Available in custom sizes, to suit your requirements perfectly.

info@surgicalholdings.co.uk

ENT Basket

Wire mesh base with perforated sides and two carrying handles, supplied with lid




Sieve



Arthroscope Basket

Wire mesh base with two fixing clamps



SC-GI/CO001 290 x 80 x 52mm

Endoscope Basket

With small lens basket

	Code	Description	
	SC-GI/CO002	6400 x 150 x 77mm	
Wire mesh base with two fixing clamps			

	abe first the fixing claim
Code	Description

SC-GI/CO003	S - 460 x 80 x 52mm
SC-GI/CO004	D - 480 x 80 x 52mm
SC-GI/CO005	D - 680 x 80 x 52mm

Fine Mesh Basket

Symmetry Surgical Quadlock Containers



Surgical Holdings offer the unique Quadlock container. The Patented design offers 4 locks, as opposed to the traditional 2, to ensure extra security and peace of mind.

Contact us for full spec and to order info@surgicalholdings.co.uk



With hinged take off / lockable lid

Code	Description
SC-GI/CO041	120 x 155 x 40mm
SC-GI/CO044	120 x 155 x 50mm
SC-GI/CO042	235 x 155 x 40mm
SC-GI/CO045	235 x 155 x 50mm
SC-GI/CO043	375 x 260 x 40mm

With separate lid

Code	Description
SC-GI/CO038	120 x 155 x 40mm
SC-GI/CO039	235 x 155 x 40mm
SC-GI/CO040	375 x 260 x 40mm

() WE ALSO MAKE CUSTOM BASKETS AND TRAYS, CONTACT US TO DISCUSS YOUR REQUIREMENTS


Optical Fibre Cables

Our optical fibre cables are compatible with all major scopes and illuminators. Every light lead is high temperature xenon compatible, rated to the highest level of thermal resistance and tested to ensure the highest light output possible.



How to order

Select an option from each of the following steps to create your part number.



General Instruments

Optical Fibre Cables

Step 2: Distal All end tips are stainless steel Wolf Snap On 6 Storz ACMI Snap-On 8 Ø 0.355 Ø 0.302 9 10 Olympus Wolf Clip-On Pilling 11 Ø 0.355 12 13 U ACMI Clip-On Mini Olympus Universal Ø 0.302 2.8 mm & 4.25 mm only 3/8-24 or 5/16-24 thread Step 3: Bundle Size Place gauge on bundle to determine diameter 2.2mm 22 2.8mm 28 3.0mm 30 3.8mm 38 4.25mm 425 4.8mm 48 0 0 Mini Olympus only Mini Olympus only Step 4: Length Step 5: Connector Angle Step 6: Sheathing Type 10ft 6ft 8ft 180° 45° 6" Reinforced* Non-reinforced 180 15 0 6 8 10 *Reinforced options 90° 90 Fully Reinforced* FR are not available with 🔇 Call us for custom lengths clear sheathing Step 7: Sheathing Colour Black Blue Clear Green Gray Red Pink White Yellow GR GY 4.8mm BK BL RD ΡK wн YL. only

Α

Adson Artery Forceps	14
Adsons Dissecting Forceps	18
Allis Tissue Forceps	20
Aneurysm Needle	44
Arthroscope Basket	47

Childe Approximating Forceps	23
Collin Retractor	38
Crile Artery Forceps	14
Crile-Wood Needle Holder	29, 31
Cushing Artery Forceps	14
Czerny Retractor	40

D

Dandy Artery Forceps	16
Deaver Abdominal Retractors	41
DeBakey Atraumatic Dissecting Forceps	18
DeBakey Needle Holder	32
Derf Needle Holder	31
Dieffenbach Bulldog Clamp	17
DIN Basket Lids	46
DIN Basket	46
Dissecting Forceps	18
Dressing Scissors	24, 26, 27, 28
Dressing Scissors, Left-Handed	24
Dunhill Artery Forceps	14
Dunhill Dissecting Forceps	18
Durham Retractor	42
Duval Tissue Forceps	20
Dyball and Lock Retractor	41

В

Babcock Tissue Forceps	20
Baby Metzenbaum Scissors	26
Backhaus Towel Clip	34
Balfour Abdominal Retractor	35
Ball and Socket Towel Forceps	34
Ballance (Volkmann) Retractor	42
Baron Scalpel Handle	44
Berkeley Abdominal Retractor	35
Berry Needle Holder/Wire Twister	31
Blalock Bulldog Clamp	17
Block End Dissecting Forceps	18
Bonney Dissecting Forceps	18
Bowlby Dilator	39
Boyd Scissors	28
Brodie Director Probe Ended, Silver Plated	43
Bruce Clarke Needle Holder	29

C

"Careless" Scissors	25
Cairns Artery Forceps	14
Canny Ryall Army Retractor	40
Canny Ryall Navy Retractor	40
Carmalt Angiotribe Forceps	16
Charnley-McIndoe Dissecting Forceps	18
Cheatle Jar	45
Cheatle Sterilizer Forceps	22

E

Elliot (Modified Deaver) Retractor	41
ENT Basket	47

$\left| \mathbf{F} \right|$

Fergusson Scissors	24
Fine Mesh Basket	47
Finochietto Needle Holder	31
Foerster Sponge Holding Forceps	22

Index

Foster Gillies Scissor/Needle Holder	33
Fraser Kelly Artery Forceps Curved	15
Fredrick Coole Dissecting Forceps	18

G

Gallipot	45
Gelpi Retractor	38
Gillies Dissecting Forceps	19
Gillies Scissor/Needle Holder	33
Gillies Skin Hook	42
Goligher Retractor	35
Gossett Abdominal Retractor	38
Guys Gallabin Suture-Holding Forceps	21

H

Halsey Needle Holder	31
Halstead Mosquito Forceps	15
Harrison Bowl Steriliser Forceps	22
Heath Stitch Scissors	25
Heiss Artery Forceps	15
Hunter Splinter Forceps	23

I

Instrument Tray	45
Iris Scissors	28
Irwin Retractor	39

J

Jackson Burrows Dissecting Forceps	19
Jean Dissecting Forceps	19
Jobson Horne Probe Dressing (Aural)	43
Joll Thyroid Retractor	38
Jones Towel Clip	34

K

Kelly Artery Forceps	15
Kelly Retractor Blade	41
Kidney Dish	45
Kilner Dissecting Scissors	24
Kilner Needle Holder	29, 31
Kilner Retractor	42
Kocher Artery Forceps	15
Kocher Langenbeck Retractor	40
Kocher Retractor	40, 42

L.

Landon Retractor	40
Lane Retractor	40
Lane Tissue Forceps	20
Lanes Dissecting Forceps	19
Langenbeck Retractor	40
Lawrence Needle Holder	29
Lawson Tait Artery Forcep	14
Ligature Scissors	26
Ligature Stitch Scissors	24
Lister Bandage Scissors	25
Lister Sinus Forceps	22
Littler Dissecting Scissors	24
Littlewood Tissue Forceps	20
Lotion Bowl	45

М

MacDonald Dissector Double Ended	43
MacPhail Needle Holder	29
Masing Needle Holder	32
Mathieu Needle Holder	30
Mayo Dressing Scissors	26
Mayo Harrington Scissors	27

Index

Mayo Instrument Pin	44
Mayo Needle Holder	30
Mayo Scissors	24, 27, 28
Mayo-Harrington Dressing Scissors	26
Mayo-Hegar Needle Holder	32
McIndoe Dissecting Forceps	19
McIndoe Modified Dissecting Forceps	19
McIndoe Scissors	25
Measuring Jug	45
Metzenbaum Dissecting Scissors	27, 28
Metzenbaum Scissors	24, 26
Michel Clip Extracting Forceps	23
Michel Forceps	23
Michel Skin Holding Forceps	23
Miles Retractor	35
Morris Retractor	41
Moyniham Retractor	41
Moynihan Artery Forceps	16
Moynihan Towel Clip	34

Ν

Naughton Morgan Needle Holder	30
Needle Holder	29
Nievert Needle Holder	32
Norfolk and Norwich Retractor	39

0

Olsen-Hegar Scissor/Needle Holder	33
Optical Fibre Cables	48-49

P

Paynes Suture Scissors	24
Petri Dish	45

Pozzi Abdominal Retractors	38
Probe, Fish Tail Silver Plated	43
Probe, Malleable Silver Plated	43
Probe, Sinus Silver	43

R

Rampley Sponge Holding Forceps	22
Ring Cutter	44
Roberts Artery Forceps	16
Robins Anchoring Forceps	23
Rochester-Mayo Needle Holder	30
Rochester-Pean Artery Forceps	16
Russian Dissecting Forceps	21
Rutherford Morrison Tissue Forceps	20
Ryder Needle Holder	32

S

Sarot Needle Holder	32
Scalpel Handle	44
Senn-Miller Retractor	42
Shardles Towel Clip	34
Shaw Tissue Forceps	21
Sieve	47
Silicone Matting	46
Single Hook Retractor	42
Spencer Wells Artery Forceps	16
Spring Bag Clip	34
Spring Clip	34
Sterilising Baskets	46-47
Stevens Tenotomy Scissors	28
Strabismus Scissors	28
Syme Aneurysm Needle	44



Т

Thomas Strabismus Scissors	25
Travers Retractor	39
Treves Dissecting Forceps	19
Tubing Scissors	25
Turner Warwick Blades	36-37
Turner Warwick Needleholder with Tungsten Carbide Jaws	36
Turner Warwick Universal Self Retaining Retractor	36

U

Utility Scissors	25

V

Volkmann Retractor	42
--------------------	----

w

Watson Cheyne Dissector	43
Waugh Dissecting Forceps	19
Weitlander Retractor	39
West Retractor	39

Our Manufacturing Partners

We're proud to be approved UK suppliers of the top global brands in surgical instruments, giving us access to thousands of both big-name and some of the more obscure, but nevertheless important, instruments.



Our strong relationships with these businesses go back many years. This means we can provide you with quality instruments from the brands you know, at very competitive prices.

BOOKWALTER.



Instrument Sales

- Surgical Holdings Instruments
 Symmetry Surgical
 Hermann Laparoscopic
- b Bissinger Bipolar
- ∠ Light Leads
- **DIN Baskets & Containers**

Refurbish & Repair

Rigid Endoscopes

- Surgical Instruments
 Laparoscopic Instrumentation
 Electrosurgery Devices
- <complex-block>

Surgical Holdings Parkside Centre, Temple Farm Ind. Est., Southend-on-Sea, Essex, SS2 5SJ E. info@surgicalholdings.co.uk T. +44 (0)1702 602050 F. +44 (0)1702 460006 surgicalholdings.co.uk



Based in the City of **Southend**