

Sterile Arthrodesis & Kirschner Wires



Why Choose Sterile?



Safer for Patients

The Scottish Executive Health Department¹ found evidence that reprocessed products “retain or acquire organic and/or chemical residues during use and reprocessing. There is evidence to suggest that reprocessed steel implants can cause inflammatory reactions to a much greater degree than ‘pristine’ devices as delivered by the manufacturer.”



Improved Performance

Repeated reprocessing creates concerns about corrosion, weakening and surface oxidation, which could affect implant failures¹.

Accredited to
ISO 13485 and
MDD 93/42/EEC

GS1
barcode

Clear
labelling

Manufactured in
Southend-on-Sea,
United Kingdom

Simple to use
and open

Peel-off labels
for patient record
traceability

Blister pouch insert
with clear product
visibility

Machine packed
protective tubing
pouch protects
kirschner wire

Available in Various Sizes and Styles








These are our most popular sizes but others are available. Please contact us with your requirements.

Size	Trocar/Trocar Arthrodesis Wire	Trocar/Plain K-Wire
0.8mm x 70mm	KWHS325A-S	KWHS325-S
0.9mm x 70mm	KWHS037-S	KWHS043-S
0.9mm x 75mm	KWHS055-S	KWHS061-S
0.9mm x 100mm	KWHS091-S	KWHS097-S
0.9mm x 125mm	KWHS127-S	KWHS133-S
0.9mm x 140mm	KWHS145-S	KWHS-151-S
0.9mm x 150mm	KWHS0009-S	KWHS0014A-S
0.9mm x 180mm	KWHS163-S	KWHS169-S
1.0mm x 70mm	KWHS038-S	KWHS044-S
1.0mm x 75mm	KWHS056-S	KWHS062-S
1.0mm x 100mm	KWHS092-S	KWHS098-S
1.0mm x 125mm	KWHS128-S	KWHS134-S
1.1mm x 90mm	KWHS269-S	KWHS270-S
1.1mm x 100mm	KWHS272-S	KWHS273-S
1.1mm x 115mm	KWHS275-S	KWHS276B-S
1.1mm x 125mm	KWHS278-S	KWHS279-S
1.1mm x 140mm	KWHS281-S	KWHS282-S
1.1mm x 150mm	KWHS0010C-S	KWHS285-S
1.1mm x 200mm	KWHS290-S	KWHS291-S
1.25mm x 150mm	KWHS0019B-S	KWHS0015A-S
1.2mm x 70mm	KWHS039-S	KWHS045-S
1.2mm x 100mm	KWHS093-S	KWHS099-S
1.2mm x 125mm	KWHS129-S	KWHS135-S
1.2mm x 140mm	KWHS147-S	KWHS153-S
1.2mm x 150mm	KWHS0010B-S	KWHS0019-S

Size	Trocar/Trocar Arthrodesis Wire	Trocar/Plain K-Wire
1.2mm x 180mm	KWHS165-S	KWHS171-S
1.2mm x 200mm	KWHS183-S	KWHS189-S
1.4mm x 70mm	KWHS040-S	KWHS046-S
1.4mm x 100mm	KWHS094-S	KWHS100-S
1.4mm x 115mm	KWHS112-S	KWHS118-S
1.4mm x 125mm	KWHS130-S	KWHS136-S
1.4mm x 140mm	KWHS148-S	KWHS154-S
1.4mm x 150mm	KWHS0011-S	KWHS0021-S
1.4mm x 180mm	KWHS166-S	KWHS172-S
1.4mm x 200mm	KWHS184-S	KWHS190-S
1.6mm x 50mm	KWHS005-S	KWHS011-S
1.6mm x 70mm	KWHS041-S	KWHS047-S
1.6mm x 100mm	KWHS095-S	KWHS101-S
1.6mm x 115mm	KWHS113-S	KWHS119-S
1.6mm x 125mm	KWHS131-S	KWHS137-S
1.6mm x 150mm	KWHS0012-S	KWHS0017-S
1.6mm x 180mm	KWHS167-S	KWHS173-S
1.6mm x 200mm	KWHS185-S	KWHS191-S
1.8mm x 200mm	KWHS283-S	KWHS475-S
2.0mm x 100mm	KWHS096-S	KWHS102-S
2.0mm x 115mm	KWHS114-S	KWHS120-S
2.0mm x 125mm	KWHS132-S	KWHS138-S
2.0mm x 150mm	KWHS0013-S	KWHS0018-S
2.0mm x 180mm	KWHS168-S	KWHS174-S
2.0mm x 200mm	KWHS186-S	KWHS192-S

- ✓ Highest grade implant quality stainless steel to ISO 5832-1
- ✓ Reduces risk of sharps injuries to Sterile Services personnel
- ✓ Surgical Holdings holds all popular sizes in stock and can manufacture any custom size, including threaded wires

Instrument Sales

-  Symmetry Surgical
-  Hermann Laparoscopic
-  Bissinger Bipolar
-  Light Leads
-  Din Baskets & Containers
-  Sterile Packed K-Wires
-  Custom Made Instruments

Refurbish & Repair

-  Surgical Instruments
-  Laparoscopic Instrumentation
-  Electrosurgery Devices
-  Rigid Endoscopes

 **Laser and Bar Coding Service
Instrument Care Training**



Surgical Holdings Parkside Centre,
Temple Farm Ind. Est., Southend-on-Sea, Essex, SS2 5SJ
E. info@surgicalholdings.co.uk
T. +44 (0)1702 602050 **F.** +44 (0)1702 460006
surgicalholdings.co.uk

