

forceps page 1

haemostatic forceps page 17

vascular clamps page 19

needle holders page 27

scissors page 41

dissectors page 55

miscellaneous page 59

retractors page 61

Dressing Forceps

Without teeth

1813102 13cm

1813103 15cm

1813104 16cm

1813105 18cm

1813106 20cm

1813107 23cm

1813108 25cm



1813102 1/1



1813104 1/1



1813106 1/1

Forceps


DeBakey Vascular and Cardiac Dissecting Forceps

Flat handle, 2.0mm atraumatic tips

1801270 Straight, 15cm

1801271 Straight, 20cm

1801272 Straight, 24cm



1801271


DeBakey Vascular and Cardiac Dissecting Forceps

Flat handle, 1.5mm atraumatic tips

1801274 Straight, 15cm

1801275 Straight, 20cm

1801276 Straight, 24cm



1801274



1801275


DeBakey Vascular and Cardiac Dissecting Forceps

Flat handle, 1.0mm atraumatic tips

1801280 Straight, 15cm

1801281 Straight, 20cm

1801282 Straight, 24cm



1801281


DeBakey Vascular and Cardiac Dissecting Forceps

Round handle, 2.0mm atraumatic tips

1801333 Straight, 15cm

1801335 Straight, 20cm

1801337 Straight, 24cm




1801335

DeBakey Vascular and Cardiac Dissecting Forceps

Round handle

- 1801342** Straight, 1.0mm atraumatic tips, 15cm
- 1801343** Straight, 1.0mm atraumatic tips, 20cm
- 1801344** Straight, 1.0mm atraumatic tips, 24cm
- 1801345** Angled, 1.2mm atraumatic tips, 21cm




1801343

1801345

DeBakey Vascular and Cardiac Dissecting Forceps

Round handle, Straight, 1.0mm atraumatic tips

- 1801350** Straight, 18cm
- 1801351** Straight, 21cm
- 1801352** Straight, 23cm




1801351

DeBakey Needle Pulling Tissue Forceps

2.0mm atraumatic tips, Tungsten carbide coated behind the 6mm length atraumatic jaw for secure needle pulling action

- 1801287** Straight, Flat handle, 20cm
- 1801346** Straight, Round handle, 20cm




1801287

1801346

Gerald DeBakey Tissue Forceps

- 1801283** Straight, 1.5mm atraumatic tips, 20cm
- 1801284** Straight, 2.0mm atraumatic tips, 20cm




1801284

Forceps

Gerald Tissue Forceps

1801285 Straight, 1.0mm atraumatic tips, 18cm

1801286 Straight, 1.0mm atraumatic tips, 20cm




1801286

Gerald Tissue Forceps

1801002 Straight, 1.0mm Platforms, 18cm

1801004 Straight, 1.0mm Platforms, 20cm




1801004

Gerald Tissue Forceps

Tungsten carbide coated tips, 1.0mm ring tips

1801241 Straight, 18cm

1801242 Straight, 20cm



1801242

Dual Vascular and Cardiac Tissue Forceps

1.2mm beak nose tip, with 10mm tungsten carbide coated platforms

1801347 20cm

1801348 23cm




1801348

Rhoton Dural Forceps

Tungsten carbide coated tips, 17.8cm


1801349 1.0mm tip width, 1×2teeth



Tissue Forceps

1801354 0.4mm, 1×2teeth, 6mm platforms, 18cm


1801355 0.4mm, 1×2teeth, 7mm platforms, 21cm



1801355

Adson Forceps

1801435 0.4mm, 1×2teeth, 15cm



1/1

Cushing Forceps

1801288 1mm serrated tips, Straight, Scraper end, 18cm

1801289 1.4mm serrated tips, Angled, 18cm

1801290 1mm serrated tips, Straight, 21cm



1801288



1801289

Barraya Tissue Forceps

1801291 2x3teeth, 20cm



Forceps

Micro Dissecting Forceps

- 1801293** 0.8mm tips, 1 x 2teeth,
Tungsten carbide coated tips, 11.5cm
- 1801294** 0.8mm serrated tips, 11.5cm

1801293

1/1



1801294

1/1



Graefe Forceps

Tungsten carbide coated tips, Straight, 10cm

- 1801295** 0.5mm delicate tips, 1 x 2 teeth
- 1801296** 0.8mm tips, 1 x 2 teeth
- 1801297** 0.5mm delicate serrated tips
- 1801298** 0.8mm serrated tips

1801295

1/1



Graefe Forceps

Tungsten carbide coated tips, Curved, 10cm

- 1801299** 0.5mm delicate tips, 1 x 2 teeth
- 1801300** 0.8mm tips, 1 x 2 teeth
- 1801301** 0.5mm delicate serrated tips
- 1801302** 0.8mm serrated tips

1801302

1/1



Graefe Forceps

Tungsten carbide coated tips, Strong curved, 10cm

- 1801303** 0.8mm tips, 1 x 2 teeth
- 1801304** 0.5mm delicate serrated tips
- 1801305** 0.8mm serrated tips

1801303

1/1



Graefe Forceps

Tungsten carbide coated tips, Angled, 10cm

- 1801306** 0.8mm tips, 1 x 2 teeth
- 1801307** 0.5mm delicate serrated tips
- 1801308** 0.8mm serrated tips

1801306

1/1



Resano Vascular Dissecting Forceps

Narrow double row of fine serrations on each jaw

- 1801309 18cm
- 1801310 23cm
- 1801311 32cm



Resano Vascular Forceps

1801312 With 3 articulations, 28cm



Adson-Brown Tissue Forceps

1801313 1.8mm tips, 7×7 teeth, 12.1cm



Adson Tissue Forceps

Tungsten carbide coated tips, 12.1cm

- 1801315 1.0mm tip width, 1×2 teeth
- 1801316 1.0mm tip width, 1×2 teeth,
with tying platforms
- 1801317 1.3mm tip width,
2×3 teeth



Adson Tissue Forceps

Tungsten carbide coated tips, 12.1cm

- 1801003 1.0mm tips, without tying platforms
- 1801088 1.0mm tips, with tying platforms



Nugent Forceps

1801005 1.2mm tips, Tungsten carbide coated
tips, Angled, 10.8cm



Forceps

Thumb Forceps

Tungsten carbide coated tips

1801006 1mm tips, 16cm

1801007 1mm tips, 20cm


1801008 1mm tips, 24cm

1801009 2mm tips, 16cm

1801010 2mm tips, 20cm

1801011 2mm tips, 24cm


1801010



Adam Forceps

Very delicate jaws, with 10mm tungsten carbide coated platforms, Round handle, Straight

1801089 1.2mm tips, 20cm



Dressing Forceps

1801090 2mm tungsten carbide coated tips, 22.5cm




Micro Tying Forceps

Tungsten carbide coated platforms, Round handle, 15cm

1801091 Straight, 0.2mm tips

1801092 Straight, 0.1mm tips




1801091

Micro Tying Forceps

Tungsten carbide coated platforms, Round handle, 15cm

1801093 Straight, 0.3mm tips

1801094 Straight, 0.2mm tips




1801093

Micro Tying Forceps

Tungsten carbide coated platforms, Round handle, 15cm

1801095 Curved, 0.3mm tips

1801096 Curved, 0.2mm tips



1801095 1/1

Micro Tying Forceps

Tungsten carbide coated platforms, Straight, Round handle, 18cm

1801097 0.7mm tips

1801098 0.5mm tips

1801099 0.3mm tips

1801100 0.2mm tips



1801099

Micro Tying Forceps


Tungsten carbide coated platforms, Curved, Round handle, 18cm

1801101 0.7mm tips

1801102 0.5mm tips

1801103 0.3mm tips

1801104 0.2mm tips



1801103

Micro Tying Forceps

Tungsten carbide coated platforms, Round handle

1801105 Straight 0.5mm tips, 20cm

1801106 Straight 0.3mm tips, 20cm

1801109 Straight 0.7mm tips, 23cm

1801110 Straight 0.3mm tips, 23cm

1801113 Straight 1.0mm tips, 25cm

1801114 Straight 0.5mm tips, 25cm

1801117 Straight 1.0mm tips, 30cm

1801119 Straight 1.0mm tips, 33cm

1801107 Curved 0.5mm tips, 20cm

1801108 Curved 0.3mm tips, 20cm

1801111 Curved 0.7mm tips, 23cm

1801112 Curved 0.3mm tips, 23cm

1801115 Curved 1.0mm tips, 25cm

1801116 Curved 0.5mm tips, 25cm

1801118 Curved 1.0mm tips, 30cm

1801120 Curved 1.0mm tips, 33cm



1801105



1801107

Forceps

Micro Tying Forceps

With 10mm tungsten carbide coated platforms, Straight

1801121 0.2mm tips, 20cm

1801122 0.3mm tips, 23.5cm

1801123 0.3mm tips, 27.5cm

1801123



Dennis Micro Forceps

Tungsten carbide coated platforms

1801133 Straight, 0.3mm tips, 15cm

1801138 Straight, 0.5mm tips, 18cm

1801139 Straight, 0.3mm tips, 18cm

1801145 Straight, 0.5mm tips, 21cm

1801146 Straight, 0.3mm tips, 21cm

1801135 Curved, 0.3mm tips, 15cm

1801142 Curved, 0.5mm tips, 18cm

1801143 Curved, 0.3mm tips, 18cm

1801147 Curved, 0.5mm tips, 21cm

1801148 Curved, 0.3mm tips, 21cm

1801139



1801143



Micro Ring Tip Forceps

Tungsten carbide coated tips, 1.0mm ring tips, Round handle

1801210 15cm

1801214 18.1cm

1801215 20cm


1801216 21cm

1801217 23cm

1801218 25cm

1801219 27.5cm

1801210



1801219



Micro Ring Tip Forceps


Tungsten carbide coated tips, Round handle

1801220 1.3mm ring tips, 15cm

1801221 1.3mm ring tips, 18cm

1801228 1.5mm ring tips, 20cm









1801230 1.5mm ring tips, 23cm



1801228



Jewelers Forceps

<p>1801012 1# Straight, 0.2×0.12mm tips, 12cm</p>		<p>1/1</p>
<p>1801013 1# Straight, 0.1×0.06mm tips, 12cm</p>		<p>1/1</p>
<p>1801014 2# Straight, 0.34×0.14mm tips, 12cm</p>		<p>1/1</p>
<p>1801015 2# Straight, 0.17×0.07mm tips, 12cm</p>		<p>1/1</p>
<p>1801016 3# Straight, 0.17×0.1mm tips, 12cm</p>		<p>1/1</p>
<p>1801017 3# Straight, 0.08×0.04mm tips, 12cm</p>		<p>1/1</p>
<p>1801018 3c# Straight, 0.08×0.04mm tips, 11cm</p>		<p>1/1</p>
<p>1801019 4# Straight, 0.13×0.08mm tips, 11cm</p>		<p>1/1</p>

Forceps

Jewelers Forceps

1801020 4# Straight, 0.06×0.02mm tips,
11cm

1/1



1801022 5# Straight, 0.05×0.01mm tips,
11cm

1/1



1801023 5/45°# Bend, 0.05×0.01mm
tips, 11cm

1/1



1801025 7# Curved, 0.17×0.1mm tips,
11.7cm

1/1



1801026 0.15×0.12mm tips, 15.5cm

1/1



1801027 5a#, 0.1×0.06mm tips, 11.7cm

1/1



1801028 5a#, 0.13×0.08mm tips, 11.7cm

1/1



Yasargil Micro Forceps

Tungsten carbide coated tips, Bayonet style, Straight

- 1801600** 0.6mm tips, 21.6cm
- 1801604** 0.9mm tips, 24cm
- 1801606** 0.9mm tips, 25cm



Yasargil Bayonet Forceps

- 1801602** Straight, 1.0mm tips, 23cm
- 1801602A** Straight, 0.5mm tips, 23cm
- 1801650** Curved, 1.0mm tips, 23cm
- 1801650A** Curved, 0.5mm tips, 23cm
- 1801651** Angle up, 1.0mm tips, 23cm
- 1801651A** Angle up, 0.5mm tips, 23cm



Micro Bayonet Forceps

Bayonet style, Straight, Tungsten carbide coated tips

- 1801608** 0.4mm tips, 26.3cm




Forceps

Yasargil Micro Forceps

Tungsten carbide coated tips

1801610 1.3mm tips,
1×2teeth,
21.6cm

1801610A 0.7mm tips,
1×2teeth, 21.6cm



1801610

Yasargil Micro Forceps

Bayonet style, Straight, Tungsten carbide coated tips

1801620 1.0mm ring tips, 23cm



Bipolar Coagulation Forceps

1801630 Bayonet style, Straight 0.5mm tips, 22cm



Bipolar Coagulation Forceps

1801640 Bayonet style, Straight 0.7mm tips, 25cm




Carpentier Papillary Muscle Clamp

1800452 17cm

1800440 24cm

1800441 27cm



1800441


Clamps

O'shaughnessy Clamping Forceps

Ring handle



1800437B



18.5cm

20cm
22cm
24cm

28cm

Right Angle

1/2 Curved









Curved


1800435B	18.5cm	1800435A	18.5cm	1800435	18.5cm
1800436B	20cm	1800436A	20cm	1800436	20cm
1800437B	22cm	1800437A	22cm	1800437	22cm
1800438B	24cm	1800438A	24cm	1800438	24cm
1800439B	28cm	1800439A	28cm	1800439	28cm

Haemostatic Forceps

Halstead Mosquito Forceps

Very delicate serrated jaws

- | | | | |
|----------------|---|---|-----|
| 1814115 | 0.6mm tips,
Straight, 12.5cm |  | 1/1 |
| 1814116 | 0.6mm tips,
Slightly curved, 12.5cm |  | 1/1 |
| 1814117 | 0.6mm tips,
Curved, 12.5cm |  | 1/1 |
| 1814118 | 0.6mm tips,
Right angle, 12.5cm |  | 1/1 |
| 1814315 | Straight, 0.6mm tip width,
1 × 2teeth, 12.5cm |  | 1/1 |
| 1814316 | Slightly curved, 0.6mm tip width,
1 × 2teeth, 12.5cm |  | 1/1 |
| 1814317 | Curved, 0.6mm tip width,
1 × 2teeth, 12.5cm |  | 1/1 |
| 1814318 | Right angle, 0.6mm tip width
1 × 2teeth, 12.5cm |  | 1/1 |



Haemostatic Forceps

Halstead Mosquito Forceps

1800404 Straight, Serrated jaws, 12.5cm




1800405 Curved, Serrated jaws, 12.5cm



1800406 Straight, 1×2teeth, 12.5cm



1800407 Curved, 1×2teeth, 12.5cm



Haemostatic Forceps

Serrated jaws

1814113 Straight, 14cm




1814133 Curved, 14cm



1814313 Straight, 1×2teeth, 14cm



1814333 Curved, 1×2teeth, 14cm



Vascular Clamps

DeBakey-Diethrich Bulldog Clamp

Cross-action, Atraumatic tips, Tension 180gms

1820450 Straight, 2×11mm jaws, 55mm

1820455 Straight, 2×15mm jaws, 61mm

1820470 Straight, 2×20mm jaws, 69mm

1821450 Angled, 2×11mm jaws, 51mm

1821455 Angled, 2×15mm jaws, 57mm

1821470 Angled, 2×20mm jaws, 65mm

1821450



1/1

1821455



1/1

Diethrich Bulldog Clamp

Cross-action, Serrated tips, Tension 80gms

1822450 Angled, 2×11mm jaws, 51mm

1822455 Angled, 2×15mm jaws, 57mm

1822470 Angled, 2×20mm jaws, 65mm

1823450 Straight, 2×11mm jaws, 55mm

1823455 Straight, 2×15mm jaws, 61mm

1823470 Straight, 2×20mm jaws, 69mm

1822450



1/1

1822455



1/1

1823455




1/1

Dennis Multipurpose Clamp

1800280 6cm atraumatic jaw, 18.5cm

1800288 8cm atraumatic jaw, 20cm



1800280


1800280 1/1

1800288 1/1

Vascular Clamps

ck Infant Aorta Clamp


0366 Miniature aorta clamp, Atraumatic jaw, 29mm jaw length, 7mm jaw depth, 16.5cm



1/1

bramanian Infant Aorta Clamp

0368 Miniature aorta clamp, Atraumatic jaw, 30mm jaw length, 7mm jaw depth, 16cm



Top View



1/1

Vascular Clamps

DeBakey Aortic Exclusion Clamp


1800369 76mm jaw length, 14mm jaw depth, 18.8cm



DeBakey Aortic Clamp


1800370 Kay aorta exclusion clamp, Atraumatic jaw, 51mm jaw length, 32mm jaw depth, 20cm

1800372 Lambert-Kay aorta exclusion clamp, Atraumaic jaw, 51mm jaw length, 20mm jaw depth, 20cm




Cooley Aorta Clamp

1800374 Atraumatic jaw, 27cm




Vascular Clamps

1800376 DeBakey Aorta Clamp, Atraumatic jaw, 27cm



1800378 Glover Clamp, 23.2cm



1800388 DeBakey Aorta Clamp, 27cm

