



BB Retraktor-Set

| *BB Retractor set*






BB Retraktor-Set | BB Retractor set

F1

Edelstahl | Stainless steel


T1

Titan | Titanium




	max. Spreizweite max. opening	Armlänge Arm length
MB1-10-769-08 F1	65 mm	80 mm
MB1-10-769-07 T1	65 mm	80 mm

BB-Retraktor rechts
BB-Retractor right




	max. Spreizweite max. opening	Armlänge Arm length
MB1-10-769-13 F1	65 mm	80 mm
MB1-10-769-12 T1	65 mm	80 mm

BB-Retraktor links
BB-Retractor left




	max. Spreizweite max. opening	Armlänge Arm length
MB1-10-769-18 F1	80 mm	80 mm
MB1-10-769-17 T1	80 mm	80 mm

BB-Retraktor Offset
BB-Retractor offset



	max. Spreizweite max. opening	Armlänge Arm length
MB1-10-769-21 F1	60 mm	60 mm
MB1-10-769-22 T1	60 mm	60 mm

BB-Retraktor compression & distraction
BB-Retractor compression & distraction



	max. Spreizweite max. opening	Armlänge Arm length
MB1-10-769-23 F1	80 mm	80 mm
MB1-10-769-24 T1	80 mm	80 mm

BB-Retraktor mit Instrumentenkupplung
BB-Retractor frame with instrument coupling



Edelstahl | Stainless steel

Titan | Titanium

Kaltlicht | cold light

F1

TI

KL



F1	TI	CL	Tiefe Depth
MB1-10-767-50	MB1-10-769-50	MB1-10-759-50	34 mm
MB1-10-767-51	MB1-10-769-51	MB1-10-759-51	40 mm
MB1-10-767-52	MB1-10-769-52	MB1-10-759-52	45 mm
MB1-10-767-53	MB1-10-769-53	MB1-10-759-53	50 mm
MB1-10-767-54	MB1-10-769-54	MB1-10-759-54	55 mm
MB1-10-767-55	MB1-10-769-55	MB1-10-759-55	60 mm
MB1-10-767-57	MB1-10-769-57	MB1-10-759-57	70 mm
MB1-10-767-59	MB1-10-769-59	MB1-10-759-59	85 mm

Valve gezahnt, 22,5mm
blade with teeth, 22.5mm



F1	TI	CL	Tiefe Depth
MB1-10-767-40	MB1-10-769-40	MB1-10-759-40	40 mm
MB1-10-767-41	MB1-10-769-41	MB1-10-759-41	47 mm
MB1-10-767-42	MB1-10-769-42	MB1-10-759-42	55 mm

Tiefen Valve gezahnt, 26,5mm
Depth blade with teeth, 26.5mm



CL-Valven mit Kaltlichtanschluss zur zusätzlichen Ausleuchtung des OP-Feldes.
CL-Blades with cold light connector for additional illumination of the OR field.

F1	TI	CL	Tiefe Depth
MB1-10-767-60	MB1-10-769-60	MB1-10-759-60	34 mm
MB1-10-767-61	MB1-10-769-61	MB1-10-759-61	40 mm
MB1-10-767-62	MB1-10-769-62	MB1-10-759-62	45 mm
MB1-10-767-63	MB1-10-769-63	MB1-10-759-63	50 mm
MB1-10-767-64	MB1-10-769-64	MB1-10-759-64	55 mm
MB1-10-767-65	MB1-10-769-65	MB1-10-759-65	60 mm
MB1-10-767-67	MB1-10-769-67	MB1-10-759-67	70 mm

Valve stumpf, 22,5mm
blade blunt, 22.5mm



F1	TI	CL	Tiefe Depth
MB1-10-767-45	MB1-10-769-45	MB1-10-759-45	40 mm
MB1-10-767-46	MB1-10-769-46	MB1-10-769-46	47 mm
MB1-10-767-47	MB1-10-769-47	MB1-10-769-47	55 mm

Tiefen Valve stumpf, 26,5mm
Depth blade blunt, 26.5mm



Kaltlicht-Adapter, MB1-10-769-27
Cold-light-adaptor



Edelstahl | Stainless steel

Titan | Titanium

F1

T1



MB1-10-769-84 F1

MB1-10-769-83 T1

Lange Zahnstange mit Rändelschraube
Long rack



MB1-10-769-87 F1

BB Adapter für Haltearm
BB adapter for holding arm




MB1-10-769-30 F1

BB Valvenhalter (RS Design)
BB Blade handle (RS Design)



MB1-10-769-35 F1

Handgriff Standard
Handle standard



MB1-10-769-37 T1

Tiefenmessstab
Depth gauge



MB1-10-769-93 F1

MB1-10-769-92 T1

Spitze Valve
Spike blade



MB1-10-769-96 F1

MB1-10-769-95 T1

T-Valve
T-blade



MB1-10-769-98 F1

MB1-10-769-97 T1

Bandvalve
Ligament blade



MB1-10-769-90 F1

MB1-10-769-91 T1

Instrumentenkupplung
Instrument attachment coupling



F1



Edelstahl | Stainless steel

T1

Titan | Titanium

BB Retraktor-Set | BB Retractor set

medicalbees



- Zervikales Retraktorsystem für anteriore und posteriore Zervikalchirurgie
- Optimale Exponierung
- Produkte erhältlich in Titan und Edelstahl
- Edelstahlvalven schwarz beschichtet (Ti ALN)
- Made in Germany

- Cervical retractor system for anterior and posterior cervical surgery
- Optimal exposure
- Products available in titanium and stainless steel
- Blades Ti ALN coated with black finish
- Made in Germany

MB1-20-968-77	Siebkorb mit Inlays und Deckel Tray with inserts and lid	Länge Length	Breite Width	Tiefe Depth
		540 mm	253 mm	50 mm

MB1-10-769-00 F1	BB Set komplett (Edelstahlversion) mit Siebkorb 1-20-968-77 bestehend aus BB set complete (stainless steel version) with tray 1-20-968-77 consisting of
------------------	--

Set-Inhalt Set content	Menge Quantity	Set-Inhalt Set content	Menge Quantity
MB1-10-769-08 F1	1	MB1-10-767-55 F1	2
MB1-10-769-13 F1	1	MB1-10-767-57 F1	2
MB1-10-769-18 F1	1	MB1-10-767-59 F1	1
MB1-10-769-30 F1	1	MB1-10-767-60 F1	2
MB1-10-769-35 F1	1	MB1-10-767-61 F1	2
MB1-10-769-37 T1	1	MB1-10-767-62 F1	2
MB1-10-767-40 F1	4	MB1-10-767-63 F1	2
MB1-10-767-41 F1	4	MB1-10-767-64 F1	2
MB1-10-767-42 F1	4	MB1-10-767-65 F1	2
MB1-10-767-45 F1	2	MB1-10-767-67 F1	2
MB1-10-767-46 F1	2	MB1-10-769-90 F1	1
MB1-10-767-47 F1	2	MB1-10-769-93 F1	1
MB1-10-767-50 F1	2	MB1-10-769-96 F1	1
MB1-10-767-51 F1	2	MB1-10-769-98 F1	1
MB1-10-767-52 F1	2	MB1-10-769-84 F1	1
MB1-10-767-53 F1	2	MB1-20-968-77	1
MB1-10-767-54 F1	2		

MB1-10-769-01 T1	BB Set komplett (Titanversion) mit Siebkorb 1-20-968-77 bestehend aus BB set complete (titanium version) with tray 1-20-968-77 consisting of
------------------	---

Set-Inhalt Set content	Menge Quantity	Set-Inhalt Set content	Menge Quantity
MB1-10-769-07 T1	1	MB1-10-769-55 T1	2
MB1-10-769-12 T1	1	MB1-10-769-57 T1	2
MB1-10-769-17 T1	1	MB1-10-769-59 T1	1
MB1-10-769-37 T1	1	MB1-10-769-60 T1	2
MB1-10-769-30 F1	1	MB1-10-769-61 T1	2
MB1-10-769-35 F1	1	MB1-10-769-62 T1	2
MB1-10-769-40 T1	4	MB1-10-769-63 T1	2
MB1-10-769-41 T1	4	MB1-10-769-64 T1	2
MB1-10-769-42 T1	4	MB1-10-769-65 T1	2
MB1-10-769-45 T1	2	MB1-10-769-67 T1	2
MB1-10-769-46 T1	2	MB1-10-769-91 T1	1
MB1-10-769-47 T1	2	MB1-10-769-92 T1	1
MB1-10-769-50 T1	2	MB1-10-769-95 T1	1
MB1-10-769-51 T1	2	MB1-10-769-97 T1	1
MB1-10-769-52 T1	2	MB1-10-769-83 T1	1
MB1-10-769-53 T1	2	MB1-20-968-77	1
MB1-10-769-54 T1	2		



MB1-20-968-81	Siebkorb mit Silikonmatte und Deckel <i>Tray with silicone mat and lid</i>	Länge <i>Length</i>	Breite <i>Width</i>	Tiefe <i>Depth</i>
		243 mm	253 mm	50 mm

MB1-10-769-02	BB HWS Mikrosset komplett (Edelstahl) mit Siebkorb 1-20-968-81 bestehend aus BB Micro Dissectomy Set complete (stainless steel) <i>with screen tray 1-20-968-81 consisting of</i>
----------------------	--

MB1-10-769-03	BB HWS Mikrosset komplett (Titan) mit Siebkorb 1-20-968-81 bestehend aus BB Micro Dissectomy Set complete (titanium) <i>with screen tray 1-20-968-81 consisting of</i>
----------------------	---

Set-Inhalt <i>Set content</i>	Menge <i>Quantity</i>	Set-Inhalt <i>Set content</i>	Menge <i>Quantity</i>
MB1-10-769-23	1	MB1-10-767-52	1
MB1-10-769-35	1	MB1-10-767-53	1
MB1-10-769-93	1	MB1-10-767-55	1
MB1-10-769-96	1	MB1-10-767-57	1
MB1-10-769-98	1	MB1-10-767-59	1
MB1-10-767-51	1	MB1-20-968-81	1

Set-Inhalt <i>Set content</i>	Menge <i>Quantity</i>	Set-Inhalt <i>Set content</i>	Menge <i>Quantity</i>
MB1-10-769-24	1	MB1-10-769-52	1
MB1-10-769-35	1	MB1-10-769-53	1
MB1-10-769-92	1	MB1-10-769-55	1
MB1-10-769-95	1	MB1-10-769-57	1
MB1-10-769-97	1	MB1-10-769-59	1
MB1-10-769-51	1	MB1-20-968-81	1

