



# Battery & Pneumatic Ortho Power Tool Repairs

All repairs carried out in our  
Southend-on-Sea workshop



 Proud to be a family-run British company



Promoting education;  
sharing knowledge and  
offering training opportunities

# Contents

3M Maxi Driver Sets (L100 & G100)	4
3M Mini Driver Set (K200) & Minos Drill (A200)	5
3M Micro Driver Set (U100) & Attachments	6
Conmed/Linvitec Handpieces & Attachments	7
DeSoutter MultiDrive Set	8
DeSoutter Small Bone Series (D)	9
DeSoutter Sternudrive	10
Halls Drill/Reamer, Oscillator & Recipricator	11
Halls Sternum Saw & Trauma Drill, Zimmer Dermatome Set	12
Halls Surgairtome Two H/P & Micro 100 Range	13
MicroAire SmartDriver & Attachments	14
MicroAire Large Bone Series (7000) & Pulse Lavage System	15
MicroAire Air Motor Set (1000) & Pencil Grip H/P (1625)	16
Stryker System 6 (Battery) & Attachments	17
Stryker System 7 (Battery) & Attachments	18
Stryker System 8 (Battery) & Cordless Driver plus Attachments	19
Synthes AO Compact Air Drive plus Attachments	20
Synthes AO Air Pen Set & MicroDrive Set	21
Synthes AO Colibri Attachments	22



To order please call +44 (0)1702 602050  
or email [info@surgicalholdings.co.uk](mailto:info@surgicalholdings.co.uk)

# Help us support global sustainability by choosing to Repair over Replace

When a power tool fails, it is easy to conclude that it needs to be replaced. But often this isn't the case. With the right expertise, these tools can be repaired, returning it to full working condition.

We can all do our bit to go greener and help the environment. By opting to repair over replace, we take the best approach for global sustainability and make the smart choice for your budget too.

Our repairs are carried out by our highly skilled, experienced technicians, so you can be fully confident in the quality of the repair. Our specialists are well practiced in the repair of all the major brands of power tools here in the UK at our established factory in Southend-on Sea, Essex.

 **Accredited to ISO 13485**



  
**90 DAY  
WARRANTY  
ON OUR  
REPAIRS**

  
**FAST  
TURNAROUND  
ON REPAIR  
AND  
DELIVERY**

  
**LOAN SETS  
AVAILABLE  
TO YOU**

## Full Repair Service for Orthopaedic Power Tools

### Specialists in repair of Stryker System 6 & 7 Power Tools

- ✓ Repairs cover: performance testing, interior cleaning, replace bearings, o-rings, controller, motor, yoke, lubricate, function and safety test
- ✓ We also repair other brands of battery hand pieces and all common pneumatic handpieces and hoses including Hall, Microaire and 3M

## 3M Maxi Driver Set



L100



G100

## 3M Mini Driver Set (K200)



K200



3M Minos (A200) with Burr Guard

## 3M Micro-Driver Set



U100



### 3M Micro-Driver Attachments



## Conmed/Linvetec Powerpro



## Desoutter MultiDrive Set



Desoutter Pneumatic Power Tools



## Desoutter Small Bone Series (D)





With Air Hose

## Halls Drill/Reamer, Oscillator & Recipicator



Halls Drill/Reamer



Halls Oscillator



Halls Recipicator

## Halls Sternum Saw & Trauma Drill



## Zimmer Dermatome Set




## Halls Surgairtome Two Set



Halls Micro 100 Range



## MicroAire Pneumatic Power Tools



MicroAire SmartDriver (6640) & Various Attachments



Zimmer type Air Hose



MicroAire Air Hose

## MicroAire Pneumatic Power Tools



Large Bone 7000 series



MicroAire PulseLavage



## MicroAire Air Motor Set



1000 series



MicroAire Pencil Grip Wiredriver (1625)

## Stryker System 6 (Battery)



With attachments



## Stryker System 7 (Battery)



Handpieces with attachments





## Stryker System 8 (Battery)



Handpieces with attachments



Stryker System 8 Cordless Driver (Battery) & Attachments



Synthes Compact Air Drive with some Attachments



## Synthes Air Pen Drive, Air Hose & Attachments



# Synthes Colibri Attachments



# Manufacturer and Repairer of Surgical Instruments



SM<sup>✂</sup>RT REPAIR

SURGIFIX  
ONLINE PORTAL

Smart Repair is our unique repair service for surgical instruments, offering the promise of long-term value for money in choosing the sustainable option, and a product replacement if we are unable to return the instrument to an 'as new' condition\*.

In choosing our innovative Smart Repair service, customers benefit from dashboard access to our Surgifix Online Portal, enabling you to book, manage and monitor repairs online with your own assigned log in. Instantly communicate directly with the repair team to stay up-to-date with the live progress of repairs and the option to view historical scope repair reports.

\*For the price of the repair. Applies to the 50 most used and repaired products, correct at the time of going to press.



Promoting  
education; sharing  
knowledge and offering  
training opportunities





## Instrument Sales



Surgical Holdings Instruments



Symmetry Surgical



Hermann Laparoscopic



Bissinger Bipolar



Light Leads



DIN Baskets & Containers



Sterile Packed K-Wires



Custom-made Instruments

## Refurbish & Repair



Surgical Instruments



Laparoscopic Instrumentation



Electrosurgery Devices



Rigid Endoscopes



**Surgical Holdings** Parkside Centre, Temple Farm Ind. Est., Southend-on-Sea, Essex, SS2 5SJ  
E. [info@surgicalholdings.co.uk](mailto:info@surgicalholdings.co.uk) T. +44 (0)1702 602050 F. +44 (0)1702 460006  
[surgicalholdings.co.uk](http://surgicalholdings.co.uk)



Based in  
the City of  
**Southend**