

Validated steriliation baskets for CMR instruments and scopes

Ensure efficient cleaning and reprocessing with our CMR Surgical approved baskets

Baskets have been designed and validated in collaboration with the CMR Decontamination team for optimum cleaning efficacy.

Baskets are manufactured from 304 Stainless Steel and electro-polished, with silicone fittings.

Hand fitted; our highly skilled technicians bench fit the baskets at our factory in Southend-on-Sea.

Our focus is to produce high quality baskets to keep your valuable equipment protected, organised and clean for use.

SC-SH001590 - Single instrument basket

510 x 125 x 75mm high with lid and appropriate fittings. 20mm wire mesh base, no handles.


SC-SH005831 - Scope basket

480 x 250 x 75mm high with lid and appropriate fittings. 20mm wire mesh base, no handles.



SC-SH001590 - Single instrument basket

SC-SH005831 - Scope basket



Manufacturer and Repairer of Surgical Instruments

SMART REPAIR

SURGIFIX
ONLINE PORTAL











Smart Repair is our unique repair service for surgical instruments, offering the promise of long-term value for money in choosing the sustainable option, and a product replacement if we are unable to return the instrument to an 'as new' condition*

In choosing our innovative Smart Repair service, customers benefit from dashboard access to our Surgifix Online Portal, enabling you to book, manage and monitor repairs online with your own assigned log in. Instantly communicate directly with the repair team to stay up-to-date with the live progress of repairs and the option to view historical scope repair reports.

*For the price of the repair. Applies to the 50 most used and repaired products, correct at the time of going to press.

Suppliers of

-  Surgical Holdings Instruments
-  Symmetry Surgical
-  Hermann Laparoscopic
-  Bissinger Bipolar
-  Light Leads
-  Din Baskets & Containers
-  Sterile Packed K-Wires
-  Custom Made Instruments

Refurbishment & Repairs of

-  Surgical Instruments
-  Laparoscopic Instrumentation
-  Orthopaedic Power Tools
-  Electrosurgery Devices
-  Rigid Endoscopes
-  Full Tray Refurbishment

