


# Hip Instruments Catalogue



 Proud to be a family-run British company

  
Promoting education;  
sharing knowledge and  
offering training opportunities

## Contents

About Surgical Holdings	2
Bone Curette	3
Elevators	3-4
Hip Set	5
Hooks	5
Mallets	6-7
Osteotomes, Chisels and Gouges	8-9
Pushers, Impactors and Skids	10
Retractor	11

## How To Order

- 1 Choose your desired product from our extensive range or send your detailed tray list to us and we'll do the work.
- 2 Email [info@surgicalholdings.co.uk](mailto:info@surgicalholdings.co.uk) or call +44 (0)1702 602050 with the product codes and we will provide you with a quote.
- 3 Upon approval of the quote we will process your order.



### NEED A BESPOKE INSTRUMENT?

CONTACT OUR FRIENDLY TEAM TO DISCUSS YOUR REQUIREMENTS

## How To Pay

30 days credit is available (subject to status checks) or we accept Visa, Mastercard, Amex cards and cheque payments.



**We're proud to present a selection of our reusable instrumentation which may be suitable for the compilation of your hip surgery sets. The sets can be tailored to your requirements.**



**We have over 50 years of knowledge and experience in our workshop. Our team of highly skilled technicians undertake a minimum apprenticeship of 5 years and are trained to manufacture instruments, rather than just repair them. This gives us 360 degree understanding of the challenge and the importance of correctly functioning medical devices.**

We care immensely about the work we do, never forgetting that instruments we create and repair are responsible for saving lives.

We insist on controlling the quality of every item we manufacture, from the materials used to the rigorous procedures we apply to testing. Our range of supporting products are serviced and repaired in-house, ensuring we are always able to exceed our customers' expectations.

Customers can receive training, participate in factory visits or simply call their local representative for guidance with instrumentation or management.


Skill-sharing has been at the heart of our business for four generations and we frequently arrange workshops and attend events at which we can educate others.



## Bone Curette

### Charnley Bone Curette

With T handle, 380mm long




Code	Description
SC-OR/N040	Large, 20 x 14mm
SC-OR/N041	Medium, 17 x 12mm
SC-OR/N042	Small, 11 x 8mm

## Elevators

### Hohmann-Muller Spike Elevator


With narrow spade



Code	Description
SC-OR/307f	220mm long, 10mm wide
SC-OR/307	230mm long, 12mm wide
SC-OR/308	240mm long, 18mm wide

### Hohmann-Aldinger Bone Elevator


Curved



Code	Description
SC-OR/304	270mm long x 24mm wide

### Hohmann Bone Elevator


Curved



Code	Description
SC-OR/305	285mm long x 34mm wide

### Hohmann Bone Lever

Curved, Leicester Royal Infirmary pattern




Code	Description
SC-OR/LGH-HOH	340mm long x 23mm wide

## Elevators

### Schumacher Bone Elevator


290mm long



Code	Description
SC-OR/307k	35mm wide, 18mm at edge

### Lange-Hohmann Bone Elevator


250mm long



Code	Description
SC-OR/307i	24mm wide

### Muller Soft Tissue Retractor

40mm wide



Code	Description
SC-OR/307aNOR	410mm long
SC-OR/307bNOR	300mm long

### Schumacher Bone Elevator


345mm long



Code	Description
SC-OR/307j	37mm wide, 18mm at edge

### Wagner Bone Elevator

220mm long



Code	Description
SC-OR/307h	17mm wide

## Sterile Arthrodesis & Kirschner Wires

- ✓ Highest grade implant quality stainless steel to ISO 5832-1
- ✓ Reduces risk of sharps injuries to Sterile Services personnel
- ✓ All popular sizes in stock and can manufacture any custom size, including threaded wires



## Hip Set

### Femoral Dislocation Lever

285mm



Code	Description
SC-OR/N059A	Spoon elevator 25mm with blue acetal suregrip handle - large

### Small Gouge

240mm



Code	Description
SC-OR/N059E	8mm cranked blue acetal suregrip handle

### Muller Curette


240mm



Code	Description
SC-OR/N059B	Curette 27mm with blue acetal suregrip handle - large

## Hooks

### Dislocation Hook



Code	Description
SC-OR/SH-5703	253mm long

### Muller Curette


270mm



Code	Description
SC-OR/N059C	Curette 9mm with blue acetal suregrip handle - small

### Dislocation Hip Hook

Blue acetal "Suregrip" handle



Code	Description
SC-OR/420	185mm long

### Curved Osteotome

260mm




Code	Description
SC-OR/N059D	With blue acetal suregrip handle - 2

## Mallets

### Surgical Holdings Orthopaedic Mallet


Stainless steel with blue acetal "Suregrip" handle  
200mm long



Code	Description
SC-OR/414BH	660 grams
SC-OR/415BH	790 grams
SC-OR/416BH	1020 grams
SC-OR/417BH	1400 grams

### Surgical Holdings Orthopaedic Mallet

Stainless steel with blue acetal "Suregrip" handle  
270mm long



Code	Description
SC-OR/410	720 grams
SC-OR/411	800 grams
SC-OR/412	900 grams
SC-OR/413	970 grams

Did you know...

**We repair power tools,  
pneumatic tools and hoses**




- ✓ Stryker, Synthes, Hall, Microaire, De Soutter, and all other common brands repaired
- ✓ All dental tools can be repaired, including bearings and rotor replacement
- ✓ We also repair Phaco Handpieces and Ultrasound Handpieces

[info@surgicalholdings.co.uk](mailto:info@surgicalholdings.co.uk)

### Orthopaedic Mallet


Surgical Holdings pattern, stainless steel,  
170mm long



Code	Description
SC-OR/N055b	800 grams
SC-OR/N055	880 grams
SC-OR/N055a	1020 grams

### Surgical Holdings Mallet

Stainless steel with blue acetal "Suregrip"  
handle domed faced head, 200mm long




Code	Description
SC-OR/418BH	650 grams

## Mallets

### Mallet

Blue acetal "Suregrip" handle, 275mm long




#### Code

SC-OR/270A-BH

### Mallet

Stainless steel, 270mm long




#### Code

SC-OR/270 640 grams

SC-OR/271 1020 grams

### Mallet

Blue acetal, 255mm long




#### Code

SC-OR/22a 520 grams, vertical end

SC-OR/22 650 grams, horizontal end

### Heath Mallet

Stainless steel, 135mm long



#### Code

SC-OR/N044c 650 grams

SC-OR/N044b 920 grams

SC-OR/N044 1150 grams

SC-OR/N044a 1340 grams

## Refurbishment of orthopaedic sets

Refurbishing  
instruments  
can save up to  
**80%**

of the cost of buying new,  
with compelling evidence to  
demonstrate effectiveness

Independent research<sup>1</sup> showed that instruments  
refurbished by Surgical Holdings had a 4x reduction in  
the amount of residual protein.

The smooth and polished surface has an increased cleaning  
efficacy, and all refurbishments come with an 18 month  
warranty as standard. Standard turnaround times are  
1 week and an express 2 day service is also available.

Reference: <sup>1</sup> Data available on file.

Book your orthopaedic set refurbishment  
on [info@surgicalholdings.co.uk](mailto:info@surgicalholdings.co.uk)



# Osteotomes, Chisels and Gouges

## Lexer Osteotome

Blue acetal "Suregrip" handle, 180mm long



Code	Description
SC-OR/SH5347a	4mm

## Capener Osteotome

Blue acetal "Suregrip" handle, 250mm long



### Shallow curved

Code	Description
SC-OR/384	Width 6mm
SC-OR/385	Width 10mm
SC-OR/386	Width 12mm
SC-OR/387	Width 15mm
SC-OR/388	Width 20mm
SC-OR/389	Width 25mm
SC-OR/390	Width 30mm

### Acute curved

Code	Description
SC-OR/390a	Width 6mm
SC-OR/390b	Width 10mm
SC-OR/390c	Width 15mm
SC-OR/390d	Width 20mm
SC-OR/390e	Width 25mm

## Capener Chisel


Blue acetal "Suregrip" handle straight, 250mm long



Code	Description	Code	Description
SC-OR/368	Width 6mm	SC-OR/371	Width 15mm
SC-OR/369	Width 10mm	SC-OR/372	Width 20mm
SC-OR/370	Width 12mm	SC-OR/373	Width 25mm

## Capener Osteotome

Blue acetal "Suregrip" handle, straight, 255mm long



Code	Description	Code	Description
SC-OR/374	Width 6mm	SC-OR/379	Width 15mm
SC-OR/375	Width 8mm	SC-OR/380	Width 18mm
SC-OR/376	Width 10mm	SC-OR/381	Width 20mm
SC-OR/377	Width 12mm	SC-OR/382	Width 25mm
SC-OR/378	Width 13mm	SC-OR/383	Width 30mm

## Moore Hollow Osteotome

Rectangular, for curved standard stems



Code	Description
SC-OR/N008	135mm long
SC-OR/N008a	255mm long

## Box Chisel

Blue acetal "Suregrip" handle, 285mm long




Code
SC-OR/N008b



## Osteotomes, Chisels and Gouges

### Capener Lamina Gouge

Cranked, blue acetal "Suregrip" handle, 250mm long



With reversed edge

Code	Description	Code	Description
SC-OR/373c	Width 6mm	SC-OR/373cre	Width 6mm
SC-OR/373d	Width 8mm	SC-OR/373dre	Width 8mm
SC-OR/373j	Width 10mm	SC-OR/373jre	Width 10mm
SC-OR/373e	Width 11mm	SC-OR/373ere	Width 11mm
SC-OR/373f	Width 16mm	SC-OR/373fre	Width 16mm
SC-OR/373k	Width 17mm	SC-OR/373kre	Width 17mm

OTHER LENGTHS AND SIZES AVAILABLE

### Capener Lamina Gouge

Swan neck, blue acetal "Suregrip" handle, 250mm long



Code	Description
SC-OR/373g	Width 8mm

### Capener Lamina Gouge


Blue acetal "Suregrip" handle, straight, 265mm long



Code	Description	Code	Description
SC-OR/362	Width 6mm	SC-OR/365	Width 15mm
SC-OR/363	Width 10mm	SC-OR/366	Width 20mm
SC-OR/364	Width 12mm	SC-OR/367	Width 25mm

### Capener Lamina Gouge

Reverse cranked, blue acetal "Suregrip" handle, 245mm long



Code	Description	Code	Description
SC-OR/373ar	Width 11mm	SC-OR/3730r	Width 6mm
SC-OR/373nr	Width 16mm	SC-OR/3731r	Width 8mm
SC-OR/373jr	Width 17mm	SC-OR/373mr	Width 10mm

## Arthroscope Repair

Surgical Holdings' highly skilled and experienced technicians will return scopes to the manufacturer's original specification. All parts are brand new and have full traceability. Surgical Holdings also manufacture their own fibre optic light bundles.



18 month guarantee on the repair



OEM optical specification followed where possible



Personal collection and delivery



Laser welding for sheaths



ISO 13485 accredited as recommended by MHRA\* advice around endoscope storage and care



Free educational visits to our scope workshop and see all our repair facilities

\*MHRA Managing Medical Devices - April 2016

## Pushers, Impactors and Skids

### Charnley Cup Pusher

With nylon head



Code	Description
SC-OR/N020	295mm long

### Hardinge Cement Restrictor

With T handle



Code	Description
SC-OR/HR001	335mm long

### Charnley Prosthesis Pusher

Stainless steel handle with nylon insert



Code	Description
SC-OR/N021	180mm long

### Moore Femoral Head Extractor


With T handle



Code	Description
SC-OR/N005	150mm long
SC-OR/N006	205mm long

### Judet Impactor

Stainless steel handle with nylon head



Code	Description
SC-OR/N009	140mm long

### Murphy Skids

Double ended, 355mm long



Code	Description
SC-OR/N007	Uncoated
SC-DIA/111	Coated

### Femoral Pusher

Head impactor, white Pomalux



Code	Description
SC-OR/N053	31mm ø x 180mm long

### Murphy Skids


Single ended, 330mm long



Code	Description
SC-OR/N007a	Coated

# Retractor

## Charnley Initial Incision Retractor



Code	Description
SC-OR/N038a	Frame only
SC-OR/N038b	Weight only, 75mm x 50mm ø
SC-OR/N038s	Goal post pins and chain









### Short blades

Code	Description
SC-OR/N038r	Short blade, 25mm deep
SH006084	Short blade, 35mm deep
SC-OR/N038c	Short blade, 50mm deep
SC-OR/N038l	Short blade, 65mm
SC-OR/N038e	Short blade, 75mm deep
SC-OR/N038q	Short blade, 80mm deep
SC-OR/N038g	Short blade, 100mm deep
SH003879	Short blade, 125mm
SH003880	Short blade, 150mm





### Long blades

Code	Description
SH006085	Long blade, 35mm deep
SC-OR/N038d	Long blade, 50mm deep
SC-OR/N038k	Long blade, 65mm
SC-OR/N038p	Long blade, 70mm deep
SC-OR/N038f	Long blade, 75mm deep
SC-OR/N038h	Long blade, 100mm deep
SH003877	Long blade, 125mm
SH003878	Long blade, 150mm
SC-OR/N038m	Long round blade, 50mm deep

## Instrument Sales

-  Surgical Holdings Instruments
-  Symmetry Surgical
-  Hermann Laparoscopic
-  Bissinger Bipolar
-  Light Leads
-  DIN Baskets & Containers
-  Sterile Packed K-Wires
-  Custom-made Instruments

## Refurbish & Repair

-  Surgical Instruments
-  Laparoscopic Instrumentation
-  Electrosurgery Devices
-  Rigid Endoscopes



**Surgical Holdings** Parkside Centre, Temple Farm Ind. Est., Southend-on-Sea, Essex, SS2 5SJ  
E. [info@surgicalholdings.co.uk](mailto:info@surgicalholdings.co.uk) T. +44 (0)1702 602050 F. +44 (0)1702 460006  
[surgicalholdings.co.uk](http://surgicalholdings.co.uk)

