

# Full Repair Service for Orthopaedic Power Tools



**90 DAY  
WARRANTY  
ON OUR  
REPAIRS**

**FAST  
TURNAROUND  
ON REPAIR  
AND  
DELIVERY**

**+ LOAN SETS AVAILABLE**

## Help us support global sustainability by choosing to Repair over Replace

When a power tool fails, it is easy to conclude that it needs to be replaced. But often this isn't the case. With the right expertise, these tools can be repaired, returning it to full working condition.

We can all do our bit to go greener and help the environment. By opting to repair over replace, we take the best approach for global sustainability and make the smart choice for your budget too.

Our repairs are carried out by our highly skilled, experienced technicians, so you can be fully confident in the quality of the repair. Our specialists are well practiced in the repair of all the major brands of power tools here in the UK at our established factory in Southend-on Sea, Essex.

### Specialists in repair of Stryker System 6 & 7 Power Tools

- ✓ Repairs cover: performance testing, interior cleaning, replace bearings, o-rings, controller, motor, yoke, lubricate, function and safety test
- ✓ We also repair other brands of battery hand pieces and all common pneumatic handpieces and hoses including Hall, Microaire and 3M

 **Accredited to ISO 13485**

Contact the experts in orthopaedic power tools to book your repair  
**+ 44 (0) 1702 965245 | [info@surgicalholdings.co.uk](mailto:info@surgicalholdings.co.uk)**

 **surgicalholdings**

# Manufacturer and Repairer of Surgical Instruments

**SMART REPAIR**

**SURGIFIX**  
ONLINE PORTAL










Smart Repair is our unique repair service for surgical instruments, offering the promise of long-term value for money in choosing the sustainable option, and a product replacement if we are unable to return the instrument to an 'as new' condition\*

In choosing our innovative Smart Repair service, customers benefit from dashboard access to our Surgifix Online Portal, enabling you to book, manage and monitor repairs online with your own assigned log in. Instantly communicate directly with the repair team to stay up-to-date with the live progress of repairs and the option to view historical scope repair reports.


\*For the price of the repair. Applies to the 50 most used and repaired products, correct at the time of going to press.

## Suppliers of

-  Surgical Holdings Instruments
-  Symmetry Surgical
-  Hermann Laparoscopic
-  Bissinger Bipolar
-  Light Leads
-  Din Baskets & Containers
- | Sterile Packed K-Wires
-  Custom Made Instruments

## Refurbishment & Repairs of

-  Surgical Instruments
-  Laparoscopic Instrumentation
-  Orthopaedic Power Tools
-  Electrosurgery Devices
-  Rigid Endoscopes
-  Full Tray Refurbishment

 Proud to be a family-run British company  
**+44 (0)1702 602050 | [info@surgicalholdings.co.uk](mailto:info@surgicalholdings.co.uk)**  
**[surgicalholdings.co.uk](http://surgicalholdings.co.uk)**



Promoting  
education; sharing  
knowledge and offering  
training opportunities