

Trays, Baskets, Fixation and Matting



Proud to be a family-run
British company



Contents

Introduction_3

Rigid Scope Baskets_4

Optimising Surgical Instrument Cleaning_5

Stainless Steel Din Baskets_6-7

Silicone Matting_8-9

Quad Lock_10

Our Services_11



If you have any questions about the issues raised in this booklet, please contact your Surgical Holdings representative for more information on +44 (0)1702 602050

This booklet has been compiled to help theatres and sterile services manage the storage and care of surgical instruments in an organised, safe and efficient way, with optimal cleaning efficacy.

Surgical Holdings' is a family business based in Southend-on-Sea, close to London and 5 minutes from Southend Airport. We are accredited to ISO 13485.




We understand that problems with instrumentation including breakages, disorganisation and issues with washing and sterilisation can be expensive and cause frustrating delays.

This guide will provide some useful solutions to help best manage these valuable assets, and help prevent unnecessary expenditure.



Sterilising Baskets and Trays

Rigid Scope Baskets



Endoscope/Arthroscope basket

With small lens basket



SC-GI/CO002	6400 x 150 x 77mm
-------------	-------------------

Wire mesh base with two fixing clamps

SC-GI/CO001	290 x 80 x 52mm
SC-GI/CO003 S	460 x 80 x 52mm
SC-GI/CO004 D	480 x 80 x 52mm
SC-GI/CO005 D	680 x 80 x 52mm

Optimising Surgical Instrument Cleaning

- ✓ Protect patients and your valuable assets, by reducing damage to instruments during reprocessing
- ✓ Fixation means instruments can be held in an open position for improved cleaning efficacy
- ✓ Use of silicone fixings rather than matting, improves drying performance of the basket
- ✓ Tray can be designed around your set and requirements
- ✓ Cleaning of complex surgical instruments made easier
- ✓ Ideal for ENT sets, ophthalmic, plastics, micro and orthopaedic sets where sharpness of the instrument is essential
- ✓ Laparoscopic basket sets and Da Vinci can also be manufactured



Our skilled technicians determine the best way to fixate instruments into the basket by looking at what is included on the set, depending on the weight, size and type of the surgical instruments we would choose either medical grade acetal or silicone.

Baskets are manufactured from high grade 304 stainless steel. The real importance should of course be patient safety, therefore ensuring the best cleaning efficacy possible in the tray by appropriately positioning the surgical instruments and using the correct materials.

Stainless Steel DIN Baskets

Available in various shapes and sizes

DIN basket

Wire mesh base with perforated sides and two carrying handles

SC-GI/CO006	230 x 130 x 60mm	SC-GI/CO012	480 x 250 x 60mm
SC-GI/CO007	240 x 130 x 50mm	SC-GI/CO013	480 x 250 x 70mm
SC-GI/CO008	240 x 130 x 60mm	SC-GI/CO016	480 x 400 x 60mm
SC-GI/CO009	240 x 250 x 50mm	SC-GI/CO017	480 x 480 x 50mm
SC-GI/CO010	240 x 250 x 60mm	SC-GI/CO018	480 x 480 x 60mm
SC-GI/CO014	300 x 300 x 50mm	SC-GI/CO019	500 x 300 x 50mm
SC-GI/CO015	350 x 350 x 60mm	SC-GI/CO020	550 x 430 x 50mm
SC-GI/CO011	480 x 250 x 50mm	SC-GI/CO021	600 x 300 x 60mm

Wire mesh basket with two carrying handles

SC-GI/CO022	240 x 250 x 50mm	SC-GI/CO026	480 x 250 x 60mm
SC-GI/CO023	240 x 250 x 60mm	SC-GI/CO027	480 x 250 x 70mm
SC-GI/CO024	240 x 250 x 70mm	SC-GI/CO028	480 x 250 x 100mm
SC-GI/CO031	300 x 300 x 60mm	SC-GI/CO029	480 x 480 x 50mm
SC-GI/CO025	480 x 250 x 50mm	SC-GI/CO030	480 x 480 x 70mm

DIN basket lids

For perforated and wire mesh baskets

SC-GI/CO032	240 x 250mm	SC-GI/CO034	480 x 250mm
SC-GI/CO033	300 x 300mm	SC-GI/CO035	480 x 480mm





Fine mesh basket

With hinged take off / lockable lid

SC-GI/CO041	120 x 155 x 40mm
SC-GI/CO044	120 x 155 x 50mm
SC-GI/CO042	235 x 155 x 40mm
SC-GI/CO045	235 x 155 x 50mm
SC-GI/CO043	375 x 260 x 40mm

With separate lid

SC-GI/CO038	120 x 155 x 40mm
SC-GI/CO039	235 x 155 x 40mm
SC-GI/CO040	375 x 260 x 40mm



Remember: when choosing baskets and trays, ensure the steel and material is compatible with your other basket stock and instrumentation




Silicone Matting

Surgical Holdings medical grade matting can be cut to size to fit any container. Available in a number of designs, this helps cushion and protect instruments during storage, transportation and sterilisation.



**All matting
is latex
free as
standard**



Perforations allow for ventilation during the sterilisation process...




Creating air flow



Preventing bacterial growth



Keeping patients safe



Symmetry Surgical Quadlock Containers

Advanced engineered containerisation for heavy equipment and instrumentation

Filter Lids

Available in six colours. Deep-seated grooves for secure stacking. Unique 4 lock system for improved security over traditional 2 locks.

Basket

Sturdy lightweight holders for your instrument sets.

Wire Base

Prevents instruments, protective mats and baskets from touching the bottom of the container.

Security Lid

Provides an additional layer of protection for the filter.

Protective Mat

Helps protect the most delicate instruments.

Base

Available in a selection of sizes and styles. Ergonomic handles with 90° stops to protect fingers during transportation.



Manufacturer and Repairer of Surgical Instruments

SMART REPAIR

SURGIFIX
ONLINE PORTAL

Smart Repair is our unique repair service for surgical instruments, offering the promise of long-term value for money in choosing the sustainable option, and a product replacement if we are unable to return the instrument to an 'as new' condition*.

In choosing our innovative Smart Repair service, customers benefit from dashboard access to our Surgifix Online Portal, enabling you to book, manage and monitor repairs online with your own assigned log in. Instantly communicate directly with the repair team to stay up-to-date with the live progress of repairs and the option to view historical scope repair reports.











Promoting education; sharing knowledge and offering training opportunities



*For the price of the repair. Applies to the 50 most used and repaired products, correct at the time of going to press.

Suppliers of

-  Surgical Holdings Instruments
-  Symmetry Surgical
-  Hermann Laparoscopic
-  Bissinger Bipolar
-  Light Leads
-  Din Baskets & Containers
-  Sterile Packed K-Wires
-  Custom Made Instruments

Refurbishment & Repairs of

-  Surgical Instruments
-  Laparoscopic Instrumentation
-  Orthopaedic Power Tools
-  Electrosurgery Devices
-  Rigid Endoscopes
-  Full Tray Refurbishment



Questions about the issues raised in this booklet?
Contact your Surgical Holdings representative



Surgical Holdings Parkside Centre, Temple Farm Ind. Est., Southend-on-Sea, Essex, SS2 5SJ
E. info@surgicalholdings.co.uk T. +44 (0)1702 602050 F. +44 (0)1702 460006
surgicalholdings.co.uk



Version 2



Based in
the City of
Southend