



Hangzhou Kangji Medical Instrument Co., Ltd.



KANGJI MEDICAL

Hangzhou Kangji Medical Instrument Co., Ltd.

Add: No. 1668 Chunjiang East Road, Economic Development Zone, Tonglu,
Hangzhou, Zhejiang 311501, China

Marketing Department

Tel 86-571-69900059
Fax 86-571-69900069

International Department

Tel 86-571-69900010
Fax 86-571-69900020

E-mail kangji@hzkangji.com

Web www.hzkangji.com / www.kangjimed.com

Our reserves all the right to fixed product mode and specification and words



Disposable Trocar

- Durable universal seal, no reducer needed
- Threaded design maximizes abdominal wall retention
- Recessed stopcock valve
- Autosheild / Optical / Bladeless




Disposable Trocar



Disposable Trocar (Hasson)

- Durable universal seal, no reducer needed
- Threaded design maximizes abdominal wall retention
- Recessed stopcock valve
- Optical



Obturator For Hasson Trocar



Optical 101Y902 5×100mm
101Y905 10×100mm
101Y907 12×100mm

Disposable Thoracoscopic Trocar



101Y510 Φ6×65mm



101Y512 Φ10×65mm



101Y513 Φ12×65mm

Disposable Orthopedics Trocar



101Y550 Φ6.5×85mm



101Y551 Φ6.5×100mm



101Y552 Φ8.5×85mm



101Y553 Φ8.5×100mm

Disposable Laparoscopic Instrument

- Light, well balanced
- Smooth rotation of jaws
- Comfortable easy open and close handle
- Double hinge design insure maximum jaw opening and accurate jaw motion
- Colors of reel could be customized.



106Y.101 Curved Scissor, High-insulation(17mm)

106Y.103 Curved Scissor, High-insulation Short Blade(12mm)



107Y.101 Maryland(19mm)



108Y.101 Clinch(22mm)



108Y.102 Fenestrated Grasper(17mm)



108Y.103 Babcock(17mm)



108Y.104 Fenestrated Babcock(21mm)



108Y.105 Blunt(19mm)



108Y.106 Fenestrated Grasper, Long(22mm)



108Y.107 Fenestrated Bowel Grasper(35mm)



108Y.108 Fenestrated Grasping Forceps (22mm)

Disposable Monopolar HF Electrode



Hook 110Y.101

Disposable Bipolar

- Unified handle for same hand feeling
- Exceptional safety and security
- With cable



109Y.201 Straight



109Y.202 Curved

Disposable Bipolar with Cutting


- Unified handle for same hand feeling
- Exceptional safety and security
- Cutting blade, unique design
- With ratchet system
- With cable



109Y.601 Straight


Reusable Instrument (Disposable Shaft + Reusable Handle)

- Economical
- Perfect compatible with reusable handle



Disposable Retrieval Bag

- Nitinol spring for extra-strength
- TPU bag for superior durability
- Exceptional safety and security



Disposable Suction / Irrigation

- Quick probe exchange
- Phthalate free of PVC tubing set
- Transparent hand-piece for visualization of fluid
- Superior ergonomically designed for both right and left hand

