




# Ophthalmic/ Micro

Forceps	132
Scissors	134
Needle Holders	136
Retractors	137
Speculum	137
Cannula	137
Miscellaneous	138

# Forceps

## Iris Forceps


1-2 teeth, 114mm long



Code	Description
SC-OP/01	Straight
SC-SC-OP/02	Curved

## Iris Forceps


Serrated, 114mm long



Code	Description
SC-OP/03	Straight
SC-OP/04	Curved

## Iris Forceps


Extra fine, serrated, 94mm long



Code	Description
SC-OP/06	Curved

## Iris Fixation Forceps


Extra fine, 1-2 teeth, 102mm long



Code	Description
SC-OP/05	Curved

## Barraquer Corneal Forceps


76mm long



Code	Description
SC-OP/023	1-2 teeth

## Epilation Forceps


89mm long



Code	Description
SC-OP/07	Oblique jaw

## Epilation Forceps


89mm long



Code	Description
SC-OP/08	Straight ends, 3mm wide
SC-OP/08a	Round ends, 3mm wide

## Suture/Dissecting Forceps


102mm long



Code	Description
SC-OP/09	Extra fine

## Barraquer Fixation Forceps


102mm long



Code	Description
SC-OP/09a	9-10 teeth

## Jameson Evans Cilia Forceps


120mm long



Code
SC-OP/012

## Castroviejo Corneal Suture Forceps


94mm long



Code	Description
SC-OP/013	1-2 teeth

## Tarsal Cyst Forceps


89mm long



Code	Description
SC-OP/015	Round ends

## Arruga Intra Capsule Forceps


For removal of residual lens material during extra-capsula extraction, 120mm long



Code	Description
SC-OP/024	Curved tips

## St Martin Corneal Suturing Forceps


1-2 teeth, 89mm long



Code	Description
SC-OP/022	For fine work

## Desmarre Entropian Forceps


89mm long



Code	Description
SC-OP/014	Right eye
SC-OP/014a	Left eye

## St Martin Micro Corneal Suturing Forceps


1-2 teeth, 89mm long



Code	Description
SC-OP/022a	Fine pattern

## Wilde Entropian Forceps

3 sizes, 89mm long




Code	Description
SC-OP/015a	With screw action

# Scissors



FOR INFORMATION ON  
BLACK AND GOLD BOWS  
SEE PAGE 23


## de Wecker Iris Scissors



Code	Description
SC-OP/031	114mm long


## Castroviejo Scissors

102mm long



Code	Description
SC-OP/025	Curved


## Barraquer Iris Scissors



Code	Description
SC-OP/031a	62mm long

## Iris Scissors


102mm long



Code	Description
SC-OP/032	Straight
SC-OP/033	Curved

## Pooley Conjunctival Scissors


89mm long



Code	Description
SC-OP/017	For fine dissection

## Thomas Strabismus Scissors


114mm long



Code	Description
SC-OP/029	Straight
SC-OP/030	Curved

## Strabismus Scissors


Supercut, 114mm long



Code	Description
SC-GI/475	Straight
SC-GI/475a	Curved

## Iris Scissors

Supercut, 114mm long



Code	Description
SC-OP/0476	Straight
SC-OP/0476a	Curved



## Westcott Scissors

120mm long




Code	Description
SC-OP/251	Curved, blunt
SC-OP/249	Curved, sharp
SC-OP/250	Straight, blunt
SC-OP/248	Straight, sharp

## Vannas Scissors

### Scissors

Code	Description
SC-OP/241	Blunt, curved 45°, 110mm
SC-OP/245	Curved 85mm
SC-OP/244	Straight 85mm
SC-OP/246	Strongly curved 85mm



### Micro Scissors

Code	Description
SC-OP/278	Delicate, straight, 85mm
SC-OP/280	Delicate, curved, 85mm
SC-OP/279	Delicate, strongly curved in plane, 85mm

## Did you know...

## The history of our family business

We've had a long tradition of surgical instrument makers in the Coole family, going back way before our Managing Director, Daniel, was even born.

Fred Coole, started his ophthalmic instrument making business back in 1940. Instrument workshops used to be in the West End of London, close to key hospital sites.


Instrument makers would go to work dressed in 3 piece suits, as they would have to walk through showrooms where surgeons would be viewing instruments.



# Needle Holders

## Castroviejo Needle Holder

Code	Description
SC-OP/021	140mm long, with catch
SC-OP/020	140mm long, without catch



### Diamond dust jaws


Code	Description
SC-THO/060	180mm long, without catch
SC-THO/061	200mm long, without catch
SC-THO/062	215mm long, without catch



ALSO AVAILABLE WITH CATCH - ADD SUFFIX CTH TO CODE

### Tungsten carbide jaws

Code	Description
SC-THO/058a	180mm long
SC-THO/058	200mm long
SC-THO/059	215mm long




ALSO AVAILABLE WITH CATCH - ADD SUFFIX CTH TO CODE

## Silcock Needle Holder

140mm long


Code	Description
SC-OP/21a	With catch
SC-OP/20	Without catch



## Micro Needle Holder

Round knurled handle, 120mm long


Code	Description
SC-OP/YA4020	Curved, without catch
SC-OP/YA4050	Curved, with catch
SC-OP/WA2010	Straight, without catch



## Micro Needle Holder

Curved, without catch


Code	Description
SC-OP/YA4010	110mm long



# Retractors, Speculums and Cannulas

## Lang Eyelid Retractor

76mm long



Code	Description
SC-OP/011b	Left eye
SC-OP/011c	Right eye

## Desmarre Eyelid Retractor


Single, 140mm long



Code	Description
SC-OP/013	12.7mm wide blade
SC-OP/014	15.9mm wide blade
SC-OP/015	19mm wide blade

## Barraquer Eye Speculum


Stainless steel wire, 44mm long



Code	Description
SC-OP/006	Standard
SC-OP/007	Small

## Speculum


Wire, 30mm wide x 45mm long



Code	Description
SC-OP/8004012	Titanium
SC-OP/8004012ss	Stainless

## Speculum


Wire, 25mm wide x 40mm long



Code	Description
SC-OP/8004013	Titanium
SC-OP/8004013ss	Stainless

## Speculum


Wire, 20mm wide x 30mm long



code	Description
SC-OP/8004013a	Titanium
SC-OP/8004013a-ss	Stainless

## Clark Speculum

Wire, 35mm wide x 84mm long



Code	Description
SC-SH/5746	
SC-SH/5746L	Infant - left
SC-SH/5746R	Infant - right

## Lang Lachrymal Cannula

38mm long



Code	Description
SC-OP/026	Cannula = 1mm ø

## Cannula


Curved, 44mm long



Code	Description
SC-OP/026	Flat ended

## Harrison Butler Lachrymal Cannula

38mm long




Code	Description
SC-OP/027	Cannula = 0.6mm ø

# Miscellaneous

## Spring Bag Clip


Serrated



Code	Description
SC-GI/345	26mm wide

## Castroviejo Calliper


89mm long



Code	Description
SC-OP/028	Stainless
SC-OP/8006007	Titanium

## Graefe Knife

127mm long



Code
SC-OP/010

## Tooke Knife

122mm long



Code	Description
SC-OP/035	6 sizes

## Chalazian Sharp Spoons

122mm long



Code	Description
SC-OP/034	6 sizes

## Snellen Vectis

Blunt, long loop, 140mm long



Code	Description
SC-OP/018	Straight
SC-OP/019	Angled