



Plastic/Oral

Scissors	181
Needle Holders	182
Forceps	183
Raspatories and Elevators	184
Raspatories	185
Spring Forceps, Hooks and Picks	187
Retractors	188
Osteotomes, Chisels and Gouges	190
Miscellaneous	191




FOR INFORMATION ON
BLACK AND GOLD BOWS
SEE PAGE 23

Scissors

Plastic/Oral

Kilner Scissors


Rounded blades, 115mm long



Code	Description
SC-PL/001a	Straight, sharp points
SC-PL/001b	Curved, blunt points

McIndoe Scissors

Curved on flat, round tapering ends



Code	Description
SC-PL/007	190mm long

Stevens Tenotomy Scissors


Supercut, 114mm long



Code	Description
SC-OP/424	Curved/blunt
SC-OP/425	Curved/sharp
SC-OP/426	Straight/blunt
SC-OP/427	Straight/sharp

Aufricht Dissecting Scissors


Heavy pattern, triangular blades, 140mm long



Code	Description
SC-PL/002a	Straight
SC-PL/002d	Curved

Cartilage Cutting Scissors


Fine points, serrated edges



Code	Description
SC-PL/009	190mm long

Aufricht Dissecting Scissors


Light pattern, 150mm long



Code	Description
SC-PL/002b	Straight
SC-PL/002c	Curved

Gorney Supercut Scissors


Rhytidectomy scissor, tungsten carbide, one micro serrated and one razor sharp edge



Code	Description
SC-PL/085a	125mm long
SC-PL/085	200mm long, straight
SC-PL/085b	200mm long, curved

Scissors


Curved, sharp points



Code	Description
SC-PL/004	130mm long

Scissors


Angular, for cutting cartilage, wire cutting notch



Code	Description
SC-PL/006	127mm long

Fickling (Cathcart) Scissors

Angled



Code	Description
SC-PL/049	180mm long

Needle Holders



FOR INFORMATION ON
BLACK AND GOLD BOWS
SEE PAGE 23

Foster Scissor/Needle Holder


Fenestrated jaw, 125mm long



Code	Description
SC-PL/010a	Left

Lawrence Needle Holder

Fine jaw, 150mm long



Code	Description
SC-PL/012a	
SC-PL/012b	Tungsten carbide, with two gold bows

Derf Needle Holder


Tungsten carbide insert,
with two gold bows



Code	Description
SC-PL/013	115mm long
SC-PL/013L	125mm long


Reverdin Suture Needle

With automatic eye, 203mm long



Masing Needle Holder

Tungsten carbide insert,
one angled bow




Code	Description
SC-GI/009R	130mm long

Code	Description
SC-PL/013a	Fully curved
SC-PL/013b	Curved
SC-PL/013c	Slightly curved
SC-PL/013d	Straight

Castroviejo Needle Holder

Code	Description
SC-OP/021	140mm long, with catch
SC-OP/020	140mm long, without catch



Diamond dust jaws


Code	Description
SC-THO/060	180mm long, without catch
SC-THO/061	200mm long, without catch
SC-THO/062	215mm long, without catch



ALSO AVAILABLE WITH CATCH - ADD SUFFIX CTH TO CODE

Tungsten carbide jaws


Code	Description
SC-THO/058a	180mm long
SC-THO/058	200mm long
SC-THO/059	215mm long



ALSO AVAILABLE WITH CATCH - ADD SUFFIX CTH TO CODE

Mustarde Suture Holding Forceps


Convex and concave jaw



Code	Description
SC-PL/011	133mm long

Walsham Septum Straightening Forceps


220mm long



Code	Description
SC-PL/030a	Straight
SC-PL/030b	Right
SC-PL/030c	Left

Fickling Artery Forceps


190mm long



Code	Description
SC-PL/044a	Serrated
SC-PL/044b	1-2 teeth


Bone Cutting Forceps

Straight, compound action



Code	Description
SC-PL/033	146mm long


Asch Septum Forceps



Code	Description
SC-PL/031	235mm long

Ward Bone Cutting Forceps


Slightly angled on flat



Code	Description
SC-PL/038	178mm long

Jansen-Middleton Septum Forceps


Compound action, 178mm long



Code	Description
SC-PL/032	5mm bite

Symmetry Vein Strippers

- Proven, cost effective removal of saphenous varicose veins
- Eliminate the need to clean and process them after use
- 30-year heritage of quality



info@surgicalholdings.co.uk

Raspatories and Elevators


Cleft Palette Raspatory



Code	Description
SC-PL/091	Curved downwards
SC-PL/092	Curved upwards

Cleft Palette Raspatory

Curved upwards



Code	Description
SC-PL/091	Left
SC-PL/092	Right

McIndoe Raspatory

For cleft palette, spoon end, 170mm long



Code	Description
SC-/PL088	Small
SC-/PL089	Medium
SC-/PL090	Large

Farabeuf Elevator

Light, 160mm long



Code	Description
SC-PL/024	12mm curved
SC-PL/024a	8mm curved

Howarth Elevator



Code	Description
SC-PL/025	220mm long

Syme Periosteal Elevator


Double ended



Code	Description
SC-PL/046	170mm long

Cryer Elevator


150mm long



Code	Description
SC-PL/069	Right
SC-PL/069a	Left

Coupland Elevator


Curved, 150mm long



Code	Description
SC-PL/068	3mm wide
SC-PL/068a	3.5mm wide
SC-PL/068b	4mm wide

Aufricht Raspatory

9mm wide, 220mm long




Code	Description
SC-PL/095	Up stroke
SC-PL/095a	Down stroke



Durogrip Raspatory

7mm wide, 205mm long




12/14 teeth

Code	Description
SC-PL/096	12/14 teeth
SC-PL/096a	16/20 teeth
SC-PL/096b	24/30 teeth
SC-PL/096c	36/44 teeth
SC-PL/096d	52/62 teeth

Gallagher Raspatory


190mm long



Code
SC-PL/098

Joseph Raspatory


Chequered, 165mm long



Code	Description
SC-PL/099	Fine, 7.5mm
SC-PL/099a	Coarse, 9mm

Lewis Raspatory


8mm wide, 175mm long



Code	Description
SC-PL/100	Fine teeth
SC-PL/100a	Coarse teeth

Fomon Raspatory

6.5 and 7.5mm wide, 215mm long




Code	Description
SC-PL/097	Chequered
SC-PL/097a	Knotched

Raspatories

Maltz Raspatory

180mm long


Code	Description
SC-PL/101	Straight, 9mm
SC-PL/101b	Straight, 5mm
SC-PL/101c	Curved, 5mm




Miller Raspatory

180mm long

Code	Description
SC-PL/102	4.3 & 6mm
SC-PL/102a	5 & 6mm
SC-PL/102b	5 & 8mm



Passivation Process



Did you know...

Stainless steel is self passivating which provides natural corrosion resistance.

Passivation is the natural process that occurs to a surgical instrument containing chromium. The chromium in the stainless steel naturally reacts with oxygen in the air, creating the very thin passive layer which protects the surface of the steel.

Spring Forceps, Hooks and Picks

Gillies Dissecting Forceps



Code	Description
SC-PL/052	1 x 2 teeth, 152mm long

McIndoe Dissecting Forceps

152mm long



Code	Description
SC-PL/053a	1-2 teeth
SC-PL/053b	Serrated

Waugh Dissecting Forceps

Fine, serrated



Code	Description
SC-GI/261	152mm long
SC-GI/262	178mm long
SC-GI/263	203mm long

College Pattern Stopping and Dressing Forceps

152mm long



Code	Description
SC-PL/048	Serrated

Mustarde Dissecting and Suture Holding Forceps


Combined double hooked jaws



Code	Description
SC-PL/050	156mm long

Childe Approximating Forceps

With magazine for Michel clips



Code	Description
SC-PL/073	180mm long

Mitchell Trimmer

165mm long



Code
SC-ENT/013

Dental Pick



Code	Description
SC-PL/067	160mm long

Gillies Skin Hook

190mm long



Code	Description
SC-PL/065a	Standard
SC-PL/065b	Half Size



Joseph Double Skin Hook

165mm long



Code
SC-PL/056a

Kleinert Kutz Skin Hooks




Code	Description
SC-PL/109	16cm large
SC-PL/108	16cm medium
SC-PL/107	16cm small

Retractors

Alms Skin Retractor

Self retaining



Code	Description
SC-PL/060	70mm long
SC-PL/060a	100mm long

Kilner Retractor

Double ended, cat's paw



Code	Description
SC-PL/056	152mm long

Kilner Alae Retractor

85mm long



Code	Description
SC-PL/063a	10mm
SC-PL/063b	13mm

Austin Kay Retractor


85mm long



Code	Description
SC-PL/074	16mm wide

Bowdler Henry Retractor


170mm long



Code	Description
SC-PL/076	12mm wide

Tongue Retractor


190mm long



Code	Description
SC-PL/071	12mm wide

Cheek Retractor

140mm long



Code	Description
SC-PL/078	40mm wide

Ragnell Davis Retractor

140mm long




Code
SC-PL/074a

symmetry surgical
authorised distributor

Bovie

Bovie offer Electrosurgery units from 10 watts, suitable for Desiccation, all the way up to 300 units.


We recommend the A1250S as the all round unit, for plastic surgery as it has both cut and coagulate power, up to 120 watts.



For information on the full Bovie range email info@surgicalholdings.co.uk

Flap Retractor


180mm long



Code	Description
SC-PL/079	32mm wide

Kilner Cheek Retractor

Double ended, 167mm long



Code	Description
SC-PL/058a	Small, 19mm and 25mm wide
SC-PL/058b	Large, 29mm and 35mm wide

Kilner Nasal Retractor


152mm long



Code	Description
SC-PL/061a	7mm wide

Laster Cheek Retractor


110mm long



Code	Description
SC-PL/081	35mm wide

Hajak Retractor


Upper lip, 165mm long



Code	Description
SC-PL/059	20mm wide

Freece Oral Retractor

150mm long



Code	Description
SC-PL/080	33mm wide

Symmetry Skin Markers

Designed to provide a permanent skin marking solution for intact skin that won't wash off while preparing a patient for a procedure

- Protects risk of cross contamination
- Sterile Skin Markers require no preoperative preparation
- Packaged with a flexible, plastic ruler that follows body contours for easy marking and reference photos
- Available in standard tip for broad clear lines as well as in micro tips for extra precision
- Packaged 36 pre-sterile Skin Markers per box



Choice of
standard tips
or micro tips



info@surgicalholdings.co.uk

Osteotomes, Chisels and Gouges

Ward (East Grinstead) Osteotome

Square handle, 178mm long



Code	Description
SC-PL/040aeg	3mm wide
SC-PL/040beg	5mm wide

Silver Nasal Chisel

182mm long



Code	Description
SC-PL/029a	Straight
SC-PL/029b	Curved to right
SC-PL/029c	Curved to left


Gillies (French Pattern) Osteotome

140mm long



Code	Description
SC-PL/026	3mm wide
SC-PL/026a	5mm wide
SC-PL/026b	7mm wide
SC-PL/026c	9mm wide
SC-PL/026g	11mm wide
SC-PL/026d	16mm wide

Pterigoideus (Maxilla) Chisel



Code	Description
SC-PL/087	255mm long

Jenkins Chisel

140mm long



Code	Description
SC-PL/027a	2mm wide
SC-PL/027b	4mm wide
SC-PL/027c	6mm wide
SC-PL/027d	8mm wide
SC-PL/027e	10mm wide

Tilley Septum Gouge


V shaped, straight



Code	Description
SC-PL/043	127mm long

McIndoe Nasal Chisel


152mm long



Code	Description
SC-PL/028a	5mm wide
SC-PL/028b	7mm wide
SC-PL/028c	9mm wide
SC-PL/028d	11mm wide


Fickling Dental Mallet

191mm long



Code	Description
SC-PL/044	25mm ø, tapering head

Silver Skin Graft Knife



Code	Description
SC-PL/015a	190mm long

Eckhoff Mapping Pen



Code	Description
SC-PL/019	145mm long
SC-PL/019a	Nibs packs of 10

Sommerland Surgical Marking Pen


Adjustable four-pointed tip



Code	Description
SC-PL/020	152mm long


Gillies Suction Tube

Control valve, 180mm long




Code	Description
SC-PL/075a	2.5mm
SC-PL/075b	3mm
SC-PL/075c	3.5mm

Doyen Mouth Gag




Code	Description
SC-ENT/88d	Adult
SC-ENT/88e	Child

Fergusson Ackland Mouth Gag




Code	Description
SC-ENT/209	Adult, 210mm
SC-ENT/209a	Child, 155mm

Dental Cartridge Syringe



Code	Description
SC-PL/072	115mm long

Watson Skin Graft Knife




Code	Description
SC-PL/017	Right-handed
SC-PL/017L	Left-handed

Miscellaneous

Freeman Areola Nipple Sizer


1.5mm

Code	Description
SC-PL/NIP25	25mm ø
SC-PL/NIP40	40mm ø
SC-PL/NIP45	45mm ø
SC-PL/NIP50	50mm ø
SC-PL/NIP55	55mm ø



Kelsey Fry Bone Awl


148mm long



Code	Description
SC-PL/039a	Straight
SC-PL/039b	Curved

Joseph Nasal Saw

Blade 35mm long



Horizontal shaft

Code	Description
SC-PL/034a	Cranked right
SC-PL/034b	Cranked left

Vertical shaft

Code	Description
SC-PL/034c	Cranked right
SC-PL/034d	Cranked left

Unna Comedo Expressor Double Ended Curette

2mm / 1mm holes



Code	Description
SC-PL/035	135mm long

Lang Steel Scoop


Double ended



Code	Description
SC-PL/064	7mm x 4mm, 140mm long

McHardy Curette

136mm long



Code	Description
SC-PL/036	1mm, 1.5mm, and 2mm ø

MacEwen Curette


Combined with seeker, 165mm long



Code	Description
SC-PL/045a	Oval
SC-PL/045b	Round

Nero Smooth Cut Scissors

A ceramic-coated range of scissors with unbeatable durability, life cycle cost, anti-glare, cleaning efficacy, and the smoothest ever performance.



Contact us for more details on the Nero range

info@surgicalholdings.co.uk