



STILLE

Surgical Instruments



Surgical Holdings Selected Patterns

The world's best instruments, for the world's best surgeons

For more than 180 years, we have developed and manufactured the best surgical instruments for the world's most demanding surgeons. We would like to extend a sincere thank you to all our current customers and welcome our future customers to experience the STILLE difference. In this catalog we present a selection of STILLE original surgical instruments.

Precision, durability and feel are characteristic qualities of all STILLE instruments. The majority are handcrafted by our highly skilled artisans in Eskilstuna, Sweden. Our instrument manufacturing process from melting steel rods to finished product is a long one, involving multiple stages of development. Our method of crafting STILLE instruments involves only the highest-grade steel, which provides a unique feel and reliability to the surgeon in any surgical specialty.

In this catalogue you will find a wide selection of instruments, including the original SuperCut scissor developed in 1984 and our Diamond SuperCut design. Quality goes into our STILLE

forceps, needle holders, clamps, rongeurs, curettes and retractors, as well as, a fine selection of stainless and titanium micro instruments. Common to all of our instruments is that they are smooth and consistently perform to your expectations. Our outstanding balance, solid design and ergonomics provides the surgeon with a reliable instrument throughout long and complex procedures.

We at STILLE are proud to have developed the world's first double-action bone rongeur. We also developed the world's first osteotome with the classic 8-sided handle. Many STILLE instruments have been copied, but a copy cannot measure up to the original.

“Dr. Mark L. Jewell

My patients deserve nothing but the best. STILLE has been the natural extension of my hands for over 30 years.”

- Dr. Mark L. Jewell, M.D.
Plastic Surgeon
www.markjewellmd.com

STILLE Surgical Instruments Surgical Holdings Selected Patterns

Art. no: 950-W100-0001-EN-01

Photos: Åke Gunnarsson, Vue and Jonas Frid, Studio Roland, Fanny Eriksson Palm

Printing: Arkitektkopia AB

© STILLE AB 2026

SYMBOLS

 SuperCut  TC-inserts



 Diamond Jaws

Table of Contents

6	Scissors	28	Bone Rongeurs
12	Forceps	29	Skin Hooks and Retractors
15	Needle Holders	31	Elevators
17	Micro Scissors	32	Osteotomes
19	Micro Forceps	33	Mallets
22	Micro Needle Holders	34	Various Instruments
25	Wire Cutters		





STILLE SUPERCUT EYE SCISSORS

SuperCut eye scissors with sharp or blunt tips and extremely sharp edges. Beveled blades.


ART. NUMBER	BLADES	TIPS	LENGTH	
101-8010-1	Straight	Sharp	10 cm 4"	•
101-8015-1	Curved	Sharp	10 cm 4"	•



KAYE STILLE DISSECTING SCISSORS

Scissors with serrated blades for "holding" the tissue, making it easier to cut exactly as far as intended.


ART. NUMBER	BLADES	TIPS	LENGTH
101-8022	Curved	Sharp serrated	12.5 cm 5"



JABALEY STILLE SUPERCUT DISSECTING SCISSORS

SuperCut scissors for dissection of tough tissue in hand surgery. Beveled blades.

ART. NUMBER	BLADES	LENGTH	
101-8060-1	Straight	12.5 cm 5"	•
101-8062-1	Curved	12.5 cm 5"	•



STILLE-MAYO SUPERCUT DISSECTING SCISSORS

Sturdy SuperCut scissors for dissection of tough tissue. Rounded blades.

ART. NUMBER	BLADES	LENGTH	
101-8085-17	Straight	17 cm 6.75"	•
101-8087-17	Curved	17 cm 6.75"	•







METZENBAUM STILLE SUPERCUT DISSECTING SCISSORS

The popular Metzenbaum dissecting scissors with SuperCut edges. Rounded/beveled blades.

ART. NUMBER	BLADES	LENGTH	✓
101-8175-15	Curved	15 cm 6"	•
101-8175-18	Curved	18 cm 7"	•
101-8175-20	Curved	20 cm 8"	•
101-8175-23	Curved	23 cm 9"	•



STEVENS-STILLE SUPERCUT TENOTOMY SCISSORS

Very popular SuperCut scissors for precise cutting and dissection of fine tissues. Rounded blades.

ART. NUMBER	BLADES	LENGTH	✓
101-8166-1	Straight	12.5 cm 5"	•
101-8167-1	Curved	12.5 cm 5"	•
101-8167-151	Curved	15 cm 6"	•
101-8167-181	Curved	18 cm 7"	•



MAYO-HARRINGTON STILLE SUPERCUT OPERATING SCISSORS

Operating SuperCut scissors useful for several surgical disciplines. Beveled blades.

ART. NUMBER	BLADES	LENGTH	✓
101-8099-20	Curved	20 cm 8"	•
101-8099-23	Curved	23 cm 9"	•





101-8115-1


JOSEPH STILLE SUPERCUT NASAL SCISSORS

SuperCut nasal scissors for the lower lateral cartilage. Useful in dissecting and cutting tough tissue. Beveled blades.

ART. NUMBER	BLADES	LENGTH	
101-8115-1	Curved	15 cm 6"	•



101-8433-20



101-8435-20

GORNEY STILLE SUPERCUT FACE-LIFT SCISSORS

SuperCut double-beveled face-lift scissors with flat tips and semi-sharp outer blades. Two models with serrated blades to prevent tissue slippage.

ART. NUMBER	BLADES	SERRATED	LENGTH	
101-8433-20	Curved	No	20 cm 8"	•
101-8435-20	Curved	Yes	20 cm 8"	•



101-8458

MATARASSO STILLE SUPERCUT FACE-LIFT SCISSORS

SuperCut face-lift scissors with double-beveled semi-sharp outer blades.

ART. NUMBER	BLADES	LENGTH	
101-8458-17	Curved	17 cm 6.75"	•
101-8458-20	Curved	20 cm 8"	•




101-8164

LITTLER STILLE SUTURE CARRYING SCISSORS

Suture carrying scissors. Rounded blades.

ART. NUMBER	BLADES	LENGTH
101-8164	Curved	12.5 cm 5"






ADSON STILLE TISSUE FORCEPS

Tissue forceps with very delicate points and 1x2 fine teeth. Used in most specialties for delicate work.


ART. NUMBER	TEETH	TIPS	LENGTH
102-9019	1x2	Very delicate	12 cm 4.75"
102-9022	1x2	Cross-serrated	12 cm 4.75"
102-9022-15	1x2	Cross-serrated	15 cm 6"



ADSON-BROWN STILLE TISSUE FORCEPS

Tissue forceps with side-grasping teeth for secure, atraumatic grasping of soft tissue or cartilage.


ART. NUMBER	TEETH	LENGTH
102-9023-12	Side-grasping	12 cm 4.75"



GILLIES STILLE TISSUE FORCEPS

Tissue forceps with 1x2 fine teeth and stabilizing pin. An excellent forceps for the surgeon who wants an Adson-size tip, but with greater length.

ART. NUMBER	TEETH	TIPS	LENGTH
102-9027	1x2 fine	Serrated	15.5 cm 6.25"



STILLE DELICATE TISSUE FORCEPS

Tissue forceps with 1x2 fine teeth and stabilizing pin. Delicate model.

ART. NUMBER	TEETH	TIPS	LENGTH
102-9028	1x2 fine	Serrated	15.5 cm 6.25"






WAUGH-BROPHY TISSUE FORCEPS

Tissue forceps with 1x2 fine teeth. Very fine tips, like Adson. Useful when long and slender instruments are needed.

ART. NUMBER	TEETH	TIPS	LENGTH
102-9252-20	1x2	Cross-serrated	20 cm 8"



MCINDOE STILLE DRESSING FORCEPS

Dressing forceps with stabilizing pin and delicate jaws for grasping fine tissue.

ART. NUMBER	JAWS	LENGTH
102-9025	Serrated	15.5 cm 6.25"



MAYO-HEGAR STILLE NEEDLE HOLDER

Commonly used general purpose needle holder with TC jaws. Intended for 0 suture and finer.

ART. NUMBER	LENGTH	TC
104-15561-161	16 cm 6.25"	•
104-15561-181	18 cm 7"	•



CRILE-WOOD STILLE NEEDLE HOLDER

General purpose needle holder with TC jaws. Intended for 2/0 suture and finer.

ART. NUMBER	LENGTH	TC
104-15565-151	15 cm 6"	•






OLSEN-HEGAR STILLE NEEDLE HOLDER

Popular needle holder with scissors for suturing and cutting with the same instrument. Intended for 0 suture and finer. 3+3 ratchets, 104-15590-171 has 2+2 ratchets.


ART. NUMBER	JAWS	LENGTH	TC
104-15590-121	Scissors	12 cm 4.75"	•
104-15590-131	Scissors	13.5 cm 5.25"	•



WEBSTER STILLE NEEDLE HOLDER

Very delicate needle holder for very fine sutures. Smooth jaws or Diamond jaws. Intended for 5/0 suture and finer.


ART. NUMBER	JAWS	LENGTH
104-15546	Very delicate Smooth	12 cm 4.75"



STILLE SUPERCUT MICRO SCISSORS

Fine micro scissors with straight or curved blades with sharp or blunt tips. All models with SuperCut blades.

ART. NUMBER	TIPS	BLADES	LENGTH	TC
117-A110-15	Sharp	Straight	15 cm 6"	•
117-A110-18	Sharp	Straight	18 cm 7"	•
117-A111-15	Blunt	Curved	15 cm 6"	•
117-A111-18	Blunt	Curved	18 cm 7"	•
117-A113-18	Sharp	Curved	18 cm 7"	•




STILLE SUPERCUT MICRO VASCULAR SCISSORS


SuperCut micro scissors for fine vascular surgery. Blade length 11 mm.

ART. NUMBER	BLADES	LENGTH	TC
117-A120-45	Angled 45°	18 cm 7"	•
117-A120-125	Angled 125°	18 cm 7"	•





117-A125-45



117-A125-90



117-A125-125

STILLE SUPERCUT MICRO FINE VASCULAR SCISSORS

Extra fine blade version of STILLE SuperCut vascular micro scissors. Blade length 7.5 mm.

ART. NUMBER	BLADES	LENGTH	
117-A125-45	Angled 45°	18 cm 7"	•
117-A125-90	Angled 90°	18 cm 7"	•
117-A125-125	Angled 125°	18 cm 7"	•



117-B110

STILLE MICRO FORCEPS

Well balanced and ergonomic micro forceps of ultra high quality. With tying platform and 0.4 mm tips.

ART. NUMBER	JAWS	LENGTH
117-B110-15	Straight	15 cm 6"
117-B110-18	Straight	18 cm 7"


117-B130

STILLE DIAMOND JAWS COUNTER-BALANCED MICRO FORCEPS

Micro forceps with Diamond Jaws and Counter-balance. 0.4 mm tips.

ART. NUMBER	JAWS	LENGTH	
117-B130-18	Straight	18 cm 7"	•
117-B130-21	Straight	21 cm 8.25"	•





117-B140

STILLE MICRO RING FORCEPS

Counter-balanced micro forceps with fine Diamond Jaws ring tips.

ART. NUMBER	RING TIPS	JAWS	LENGTH	
117-B140-18	1.0x0.5 mm	Straight	18 cm 7"	•



STILLE MICRO DEBAKEY FORCEPS

Round handle, with 1.2 mm fine atraumatic DeBakey tips.

ART. NUMBER	TIPS	JAWS	LENGTH
117-B160-18	DeBakey	Straight	18 cm 7"
117-B160-21	DeBakey	Straight	21 cm 8.25"



STILLE MICRO DEBAKEY COUNTER-BALANCED FORCEPS

Round handle, with 1.2 mm fine atraumatic DeBakey tips, counter-balanced.

ART. NUMBER	TIPS	JAWS	LENGTH
117-B165-18	DeBakey	Straight	18 cm 7"
117-B165-21	DeBakey	Straight	21 cm 8.25"



STILLE TITANIUM MICRO FORCEPS

Round handle, with tying platform, 0.4 mm tips.

ART. NUMBER	JAWS	LENGTH	
117-B120-18-TI	Straight	18 cm 7"	•



STILLE TITANIUM MICRO RING FORCEPS

Round handle, with fine ring tips, counter-balanced.

ART. NUMBER	RING TIPS	JAWS	LENGTH	
117-B140-18-TI	1.0x0.5 mm	Straight	18 cm 7"	•
117-B140-21-TI	1.0x0.5 mm	Straight	21 cm 8.25"	•





117-C120

**STILLE MICRO NEEDLE HOLDER
REGULAR BOX LOCK WITH RATCHET**

Round handle, fine 0.8 mm Diamond Jaws and Regular Box Lock for sutures 7-0 and smaller. With ratchet.

ART. NUMBER	JAWS	LENGTH	
117-C120-21	Straight	21 cm 8.25"	•

**STILLE MICRO NEEDLE HOLDER
STREAMLINE BOX LOCK WITH RATCHET**

Round handle, 1.2 mm Diamond Jaws and Streamline Box Lock for sutures 5-0 and smaller. With ratchet.


ART. NUMBER	JAWS	LENGTH	
117-C125-18	Straight	18 cm 7"	•
117-C125-23	Straight	23 cm 9"	•
117-C126-18	Curved	18 cm 7"	•
117-C126-23	Curved	23 cm 9"	•




**STILLE TITANIUM MICRO NEEDLE HOLDER
REGULAR BOX LOCK**

Round handle, fine 0.8 mm Diamond Jaws and Regular Box Lock for sutures 7-0 and smaller. With or without ratchet. 117-C106-series with ultra fine 0.4 mm smooth jaws.

ART. NUMBER	RATCHET	JAWS	LENGTH	
117-C106-15-TI	Without	Curved	15 cm 6"	
117-C106-18-TI	Without	Curved	18 cm 7"	
117-C111-15-TI	Without	Curved	15 cm 6"	•
117-C111-18-TI	Without	Curved	18 cm 7"	•
117-C120-21-TI	With	Straight	21 cm 8.25"	•
117-C120-23-TI	With	Straight	23 cm 9"	•






117-C140

STILLE MICRO NEEDLE HOLDER HEAVY MODEL

Extra heavy model for handling of larger needles with TC-inserts. Intended for 2/0 suture and finer.

ART. NUMBER	RATCHET	JAWS	LENGTH	TC
117-C140-23	With	Straight	23 cm 9"	•
117-C140-25	With	Straight	25 cm 10"	•



STILLE SINGLE-ACTION JOINT WIRE CUTTER | Trådavbitare STILLE enkelledad

Single joint wire cutter with TC-inserts. For hard wire up to 0.4 mm and soft wire up to 0.6 mm.

ART. NUMBER	LENGTH	TC
105-16215	12.5 cm 5"	•



STILLE DOUBLE-ACTION JOINT WIRE CUTTER | Trådavbitare STILLE dubbelledad

Double joint wire cutter with TC-inserts. For hard wire up to 1.7 mm and soft wire up to 2.0 mm.

ART. NUMBER	LENGTH	TC
105-16223	17.5 cm 7"	•



STILLE EXTRA STRONG DOUBLE-ACTION JOINT WIRE CUTTER | Trådavbitare STILLE dubbelledad extra kraftig

Strong double joint wire cutter with TC-inserts. For hard wire up to 2.5 mm and soft wire up to 2.8 mm.

ART. NUMBER	LENGTH	TC
105-16225	22 cm 8.75"	•




105-17261

LEKSELL-STILLE BONE RONGEUR

Bone rongeur with original STILLE double-action joint and sideways slightly curved jaws.

ART. NUMBER	JAWS	LENGTH
105-17261	Curved narrow	24 cm 9.5"



108-25011-1

JOSEPH STILLE SKIN HOOK

Skin hook with sharp prongs. Available in five different models. Popular in plastic surgery.

ART. NUMBER	HOOK	SIZE	LENGTH
108-25011-1	Single		16 cm 6.25"
108-25012-7	Double	7 mm	16 cm 6.25"




STILLE BREAST RETRACTOR

Useful double-sided retractor for breast augmentation. Serrated blades in two different sizes and shapes. Ergonomic handle for maximum comfort.

ART. NUMBER	TIPS	HANDLE	LENGTH
108-26038	Double ended	Black coated	23 cm 9"





108-24320

STILLE NERVE HOOK

Blunt nerve hook on long flat shaft.

ART. NUMBER	HOOK	TIP	LENGTH
108-24320	Slightly curved	Blunt	22 cm 8.75"

“The steel is hard enough to deal with hard bone tissue and to remain sharp, yet flexible enough to avoid damaging implants.”



FREER STILLE ELEVATOR

Double-ended elevator sharp/blunt with rounded, delicate handle.

ART. NUMBER	SIZE	LENGTH
111-22140-4	4.5 mm	19.5 cm 7.75"
111-22140-6	6.5 mm	19.5 cm 7.75"



STILLE HAND OSTEOTOME

Delicate osteotome with phenolic handle. Intended for hand surgery, but also useful in other fine orthopaedic work.

ART. NUMBER	SIZE	LENGTH
111-21456-2	2 mm	16 cm 6.25"
111-21456-4	4 mm	16 cm 6.25"
111-21456-6	6 mm	16 cm 6.25"
111-21456-8	8 mm	16 cm 6.25"
111-21456-10	10 mm	16 cm 6.25"
111-21456-12	12 mm	16 cm 6.25"

“Stille’s products are famous for their high and consistent quality, with a 30-year warranty on all surgical instruments.”

STILLE TEFLON INSERT Mallet | Hammare med tefloninlägg STILLE

Mallet with stainless steel head, phenolic handle and interchangeable teflon insert.

ART. NUMBER	TYPE	WEIGHT	LENGTH
111-21825	Mallet	145 g 5 oz	17 cm 6.75"





STILLE BARRON SCALPEL HANDLE

8-sided handle for precise handling and rotation.

ART. NUMBER	BLADE	LENGTH
119-50012	10-15	14 cm 5.5"

“An investment in a STILLE instrument is an investment that will last throughout your career.”

- Dr. Jordan Deschamps-Braly, M.D.
Plastic Surgeon



**For a complete presentation
of all STILLE instruments,**

Scan the QR code to visit our website and
see our product catalogue



950-W100-0001-EN-01 © 2026-04-15 STILLE. All rights reserved

Surgical perfection. For life.



Surgical perfection. For life.

DISTRIBUTED EXCLUSIVELY IN THE UK BY:

Surgical Holdings
Parkside Centre, Temple Farm Ind. Est.,
Southend-on-Sea, Essex, SS2 5SJ
☎ +44 (0)1702 602050
✉ info@surgicalholdings.co.uk
🌐 <https://surgicalholdings.co.uk>

STILLE AB

Ekbacken 11
SE-644 30 Torshälla, Sweden
☎ +46 8 588 580 00
✉ info@stille.se
🌐 www.stille.se

STILLE SURGICAL INC.

840 Parkview Blvd
Lombard, IL 60148
United States
☎ 1.800.6651614