

# RapidClean® Convertible Kerrisons

## Ceramic Coated

- Smooth, wear-resistant surface
- Reduces glare
- Optimal gliding
- Stays sharp longer



59-3200J-BD

| 40° Up | Without Ejector | Bite (mm)          | Opening (mm) | 7" (178 mm) | 8" (203 mm) | 9" (229 mm) |
|--------|-----------------|--------------------|--------------|-------------|-------------|-------------|
|        |                 | Thin Footplate     |              |             |             |             |
|        |                 | 1                  | 9            | 58-1180-BD  | 58-1200-BD  | 58-1230-BD  |
|        |                 | 2                  | 9            | 58-2180-BD  | 58-2200-BD  | 58-2230-BD  |
|        |                 | 3                  | 9            | 58-3180-BD  | 58-3200-BD  | 58-3230-BD  |
|        |                 | 4                  | 9            | 58-4180-BD  | 58-4200-BD  | 58-4230-BD  |
|        |                 | 5                  | 9            | 58-5180-BD  | 58-5200-BD  | 58-5230-BD  |
|        |                 | 6                  | 9            |             |             | 58-6230-BD  |
|        |                 | Standard Footplate |              |             |             |             |
|        |                 | 1                  | 9            | 59-1180-BD  | 59-1200-BD  | 59-1230-BD  |
|        |                 | 2                  | 9            | 59-2180-BD  | 59-2200-BD  | 59-2230-BD  |
|        |                 | 3                  | 9            | 59-3180-BD  | 59-3200-BD  | 59-3230-BD  |
|        |                 | 4                  | 9            | 59-4180-BD  | 59-4200-BD  | 59-4230-BD  |
|        |                 | 5                  | 9            | 59-5180-BD  | 59-5200-BD  | 59-5230-BD  |
|        |                 | 2                  | 15           |             | 58-2262-BD  |             |
|        |                 | 3                  | 15           |             | 58-2263-BD  |             |
|        |                 | 4                  | 15           |             | 58-2264-BD  |             |
|        |                 | 5                  | 15           |             | 58-2265-BD  |             |
|        | With Ejector    | Thin Footplate     |              |             |             |             |
|        |                 | 3                  | 9            | 58-3180J-BD | 58-3200J-BD | 58-3230J-BD |
|        |                 | 4                  | 9            | 58-4180J-BD | 58-4200J-BD | 58-4230J-BD |
|        |                 | 5                  | 9            | 58-5180J-BD | 58-5200J-BD | 58-5230J-BD |
|        |                 | Standard Footplate |              |             |             |             |
|        |                 | 3                  | 9            |             | 59-3200J-BD |             |
|        |                 | 4                  | 9            |             | 59-4200J-BD |             |
|        |                 | 5                  | 9            |             | 59-5200J-BD |             |